




# Battery Chargers




Lester 4 Amp. Battery Charger  
**R32 / R51**



Lester 8 Amp. Battery Charger  
**R32LX / R50LX / R51LX / R51LXP**



CTE Battery Charger  
**R32 / R51**  
and  
**R32LX / R50LX / R51LX / R51LXP**



# Battery Chargers

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	C	1061411-NLA	NLA - Charger, Battery, 4 Amp. 115vac, Off Board, Dual Mode, Round Plug, Lester	1	1
1	D,G	1053161-NLA	NLA - Charger, Battery, 8 Amp, 115vac, Off Board Dual Mode Round Plug, Lester	1	1
1	E	1123249	Charger, Battery, 8 Amp. 115vac (Off Board) Single Mode Round Plug w/ AC Power Line Cord	1	1
1	E	1123249E-NLA	NLA - WARR RPLMT - Charger, Battery, 8 Amp, 115VAC (Off Board) Single Mode Round Plug w/ AC Power Line Cord	1	1
2	F	1089627-NLA	NLA - Cord, AC Power, 115VAC, 2 Prong Off-Board (11' 6" L)	1	1
2		1121086	Cord, AC Power, 115VAC, 3 Prong On-Board (6' 6-1/2" L)	1	1
3		1083929-NLA	NLA - Cable, Charger Extension, Round Plug (24" L)	1	1
4	A	1098892-NLA	NLA - Fuse - Buss ATC-15 (15 AMP)	1	1
4	B	1090117-NLA	NLA - Fuse - Buss Mini ATM-15 (15 AMP)	1	2
<p>NOTE: A - Before 4/1/03            B - After 3/31/03            C - CTE 8 Amp. charger replaced the obsoleted Lester 4 Amp. charger as of 4/1/04 .            D - Lester 8 Amp. charger was standard on R32LX / R50LX / R51LX /R51LXP models before 4/1/04 and optional after 3/31/04, when the CTE charger became standard.            E - CTE charger introduced and standard on-chair 4/1/04 for all models. Used on Gel Cell or AGM (Absorbed Glass Mat) batteries. Serial number starts with 04A00000001. For wet cell batteries refer to the Lester 8 Amp. charger.            F - No longer available as a service part. Please use 1121086 as its replacement.            G - Use Battery Charger part number 1123249 as a replacement.</p>					