


Front Casters and Forks



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1114572-NLA	NLA - Fork and Stem Assembly	1	1
	B,C	1006049	Bearing (7/16")	1	2
	B,C	1020592	Spacer, Bearing (7/16 x 5/8 x 1-3/4")	1	1
2	A,D,E	1118250-NLA	NLA - Kit, Front Wheel Mounting Hardware	1	2
3			NLA - Spacer (7/16 x 5/8 x 11/16")	1	2
4			Locknut (7/16-20)	1	1
5			Bolt, Hex Head Cap (7/16-20 x 4-3/4")	1	1
6	D	1070811	Wheel Assembly, Semi-Pneumatic, Urethane, Light Grey (8" x 1-3/4")	1	1
6	E	1101108	Wheel Assembly w/ Bearings, Semi-Pneumatic (8" x 2-1/4")	1	1

NOTE: A - Includes: Items 2-5.
 B - Not Shown.
 C - These items are included with the caster (part number 1070811 & 1101108).
 D - Used for Zoom 300.
 E - Used for Zoom 400.