


Seat Post and Hardware



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1120441-NLA	NLA - Seat Post Assembly with Hardware	1	1
			Locknut, Low Profile (3/8-16)	1	1
			Washer, Flat (3/8 x 11/16 x 1/16")	1	2
			NLA - Label, Patent Pending	1	1
			NLA - Screw, Hex Head (3/8-16 x 1-3/4")	1	1
2			NLA - Seat Post, Pivot	1	1
NOTE: A - Used on both the Zoom 300 and Zoom 400.					