



# Rear Wiring Harness



# Rear Wiring Harness

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1119608-NLA	NLA - Harness, Main Wiring (Invacare 300 Scooter)	1	1
1	A,C	1125363-NLA	NLA - Harness, Main Wiring (Invacare 400 Scooter)	1	1
2		1117809-NLA	NLA - Circuit Breaker (40 AMP)	1	1
3		1118873-NLA	NLA - Fuse, 5 Amp.	1	1
4		1117817-NLA	NLA - Fuse (10 AMP)	1	1
	B	1095503-NLA	NLA - Package, Tie Wrap (5-5/8" L)	10	1
NOTE: A - Includes: Items 2 - 4. B - Not Shown C - This part plugs into the Controller, Motor, Battery and Tiller Harness.					