


Base Frame, Seat Frame, and Hardware

Transport Ready Chairs requiring Frames and/or Frame Hardware parts must be returned to Invacare for service.



Base Frame, Seat Frame, and Hardware

Transport Ready Chairs requiring Frames and/or Frame Hardware parts must be returned to Invacare for service.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		See Chart	#1 Base Frame	1	1
1		See Chart	#1 Base Frame Heavy Duty 250 to 350 Lbs.	1	1
2-6		See Chart	#2-6 Insert Tube and Tube Spacer	1	1
7		See Chart	#7 Seat Side Frame	1	1
8		1025765	Locknut (1/4-28)	1	10
9		1055309	Screw, Socket Button Head (1/4-28 x 1-1/4")	1	8
10		1071249	Plug Button, Dome, Black (1-1/4")	1	2
11		15680X013	Anti-Rattle (1")	1	2
12		1083099	Screw, Socket Button Head (1/4-28 x 1-1/2")	1	2
13		1085111	Nut, Hex (1/4-28)	1	2
14		1074640	Screw, Allen Head (1/4-28 x 1-1/4")	1	2
15		1052957	Screw, Hex Head w/ Patch (1/4-20 x 2")	1	2
16		1026907	Washer (1/4 x 1/2 x 1/16")	1	2
17		1086389	Tube, Support, Front Rigging Hanger, Black - Right	1	2
17		1086390	Tube, Support, Front Rigging Hanger, Black - Left	1	2
18	A	1054571	Kit, Sector Block - Right	1	2
18	A	1054572	Kit, Sector Block - Left	1	2

NOTE: A - Includes Items 15 and 16.

#1 Base Frame

	Short 15-17" Deep Right	Short 15-17" Deep Left	Medium 18-19" Deep Right	Medium 18-19" Deep Left	Long 20-22" Deep Right	Long 20-22" Deep Left
Color	8881079560	8881079561-NL A	8881079564-NL A	8881079565-NL A	8881079568-NL A	8881079569-NL A
	See Note(s) A	See Note(s) A	See Note(s) A	See Note(s) A	See Note(s) A	See Note(s) A
Chrome	1079562-NLA	1079563-NLA	1079566-NLA	1079567-NLA	1079570-NLA	1079571-NLA

NOTE: A - Specify frame color when ordering.

Base Frame, Seat Frame, and Hardware

Transport Ready Chairs requiring Frames and/or Frame Hardware parts must be returned to Invacare for service.

#1 Base Frame Heavy Duty 250 to 350 Lbs.

	Short 15-17" Deep Right	Short 15-17" Deep Left	Medium 18-19" Deep Right	Medium 18-19" Deep Left	Long 20-22" Deep Right	Long 20-22" Deep Left
Color	8881092896 See Note(s) A.	8881092897 See Note(s) A.	8881092898 See Note(s) A.	8881092899 See Note(s) A.	8881092900 See Note(s) A.	8881092901 See Note(s) A.
Chrome	1092902-NLA	1092903-NLA	1092879-NLA	1092880-NLA	1092881-NLA	1092882-NLA

NOTE: A. - Specify frame color when ordering.

#2-6 Insert Tube and Tube Spacer

Chair Width	#2 Insert Tube Front (Qty. 1)	#3 Insert Tube Rear (Qty. 4)	#4 Tube Spacer Front (Color) (Qty. 1)	#4 Tube Spacer Front (Chrome) (Qty. 1)	#5 Tube Spacer Rear (Color) (Qty. 4)	#5 Tube Spacer Rear (Chrome) (Qty. 4)	#6 Tube Spacer Contracture Platform
12"	1120000-NLA	1120002-NLA	8881120006-NLA See Note(s) A.	N/A	8881120008-NLA See Note(s) A.	N/A	N/A
13"	1120001-NLA	1120003-NLA	8881120007-NLA See Note(s) A.	N/A	8881120009-NLA See Note(s) A.	N/A	N/A
14"	1077051-NLA	1077060-NLA	8881085273-NLA See Note(s) A.	1085282-NLA	8881079610-NLA See Note(s) A.	1079611-NLA	1110441-NLA
15"	1077052-NLA	1077061-NLA	8881085274-NLA See Note(s) A.	1085283-NLA	8881079612-NLA See Note(s) A.	1079613-NLA	1110444-NLA
16"	1077053-NLA	1077062-NLA	8881085275-NLA See Note(s) A.	1085284-NLA	8881079614-NLA See Note(s) A.	1079615-NLA	1110447-NLA
17"	1077054-NLA	1077063-NLA	8881085276-NLA See Note(s) A.	1085285-NLA	8881079616-NLA See Note(s) A.	1079617-NLA	1110450-NLA
18"	1077055-NLA	1077064-NLA	8881085277-NLA See Note(s) A.	1085286-NLA	8881079618-NLA See Note(s) A.	1079619-NLA	1110453-NLA
19"	1077056-NLA	1077065-NLA	8881085278-NLA See Note(s) A.	1085287-NLA	8881079620-NLA See Note(s) A.	1079621-NLA	1110455-NLA
20"	1077057-NLA	1077066-NLA	8881085279-NLA See Note(s) A.	1085288-NLA	8881079622-NLA See Note(s) A.	1079623-NLA	1110459-NLA
21"	1077058-NLA	1077067-NLA	8881085280-NLA See Note(s) A.	1085289-NLA	8881079624-NLA See Note(s) A.	1079625-NLA	1110462-NLA

Transport Ready Chairs requiring Frames and/or Frame Hardware parts must be returned to Invacare for service.

Chair Width	#2 Insert Tube Front (Qty. 1)	#3 Insert Tube Rear (Qty. 4)	#4 Tube Spacer Front (Color) (Qty. 1)	#4 Tube Spacer Front (Chrome) (Qty. 1)	#5 Tube Spacer Rear (Color) (Qty. 4)	#5 Tube Spacer Rear (Chrome) (Qty. 4)	#6 Tube Spacer Contracture Platform
22"	1077059-NLA	1077068-NLA	8881085281-NLA See Note(s) A.	1085290-NLA	8881079626-NLA See Note(s) A.	1079627-NLA	1110465-NLA

NOTE: A. - Specify frame color when ordering.

Base Frame, Seat Frame, and Hardware

Transport Ready Chairs requiring Frames and/or Frame Hardware parts must be returned to Invacare for service.

#7 Seat Side Frame

	Short 15-17" Deep Right	Short 15-17" Deep Left	Medium 18-19" Deep Right	Medium 18-19" Deep Left	Long 20-22" Deep Right	Long 20-22" Deep Left
Color	8881079598-NL A	8881079599-NL A	8881079602-NL A	8881079603-NL A	8881079606-NL A	8881079607-NL A
Chrome	1079600-NLA	1079601-NLA	1079604-NLA	1079605-NLA	1079608-NLA	1079609-NLA