


Attendant Operated Wheel Locks

12" Wheels Only



Attendant Operated Wheel Locks 12" Wheels Only

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1025765	Locknut (1/4-28)	1	4
2		1080551-NLA	NLA - Strap, Wheel Lock	1	1
3		1028411-NLA	NLA - Spacer (1/4 x 1/2 x 5/32")	1	1
4		1055306-NLA	NLA - Screw, Button Head Cap (1/4-28 x 7/8")	1	1
5		1026103	Washer, Nylon (1/4 x 1/2 x 1/32")	1	8
6		1080615-NLA	NLA - Link, Wheel Lock, Black - Right	1	1
6		1080616-NLA	NLA - Link, Wheel Lock, Black - Left	1	1
7		1083099	Screw, Socket Button Head (1/4-28 x 1-1/2")	1	1
8		1055304	Screw, Socket Button Head (1/4-28 x 5/8")	1	2
9		1020363-NLA	NLA - Spacer (1/4 x 1/2 x 1/2")	1	1
10		1018733-NLA	NLA - Link, Stop, Black	1	1
11		1021405-NLA	NLA - Shoe, Wheel Lock	1	1
12		1085267-NLA	NLA - Wheel Lock Assembly - Right	1	1
12		1085268-NLA	NLA - Wheel Lock Assembly - Left	1	1
13	F	1112935-NLA	NLA - Foot Release Assembly	1	1
13	A	1112937-NLA	NLA - Foot Release Assembly	1	1
13	B	1112941-NLA	NLA - Foot Release Assembly	1	1
13	C	1112945-NLA	NLA - Foot Release Assembly	1	1
14		1055304	Screw, Socket Button Head (1/4-28 x 5/8")	1	4
15	D	47376X024	Tip, Rubber Wheel Lock Handle, Black	1	1
15	E	1128817	Package, Tip, Rubber, Black (Pkg. of 10)	1	1
16		90036X000-NLA	NLA - Screw, Allen Head (1/4-28 x 1")	1	2
17		1079550-NLA	NLA - Clamp, Non-Threaded	1	1
18		1079551-NLA	NLA - Clamp, Threaded	1	1
19		1101207-NLA	NLA - Screw, Set Socket Cone Point (1/4-20 x 1/4")	1	2
20		1085266-NLA	NLA - Wheel Lock Clamp Hardware	1	2
NOTE: A - 14-17" Wide chairs B - 18-21" Wide chairs C - 22" Wide chairs D - Before 9/02 E - After 8/02 F - 12-13" Wide chairs					