


Front Riggings - Pivot Slide Tube w/ Footplate

AHS4, AHL4



Front Riggings - Pivot Slide Tube w/ Footplate

AHS4, AHL4

AHL4A, AHL4C, AHS4A

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	B	1023905-NLA	NLA - Pivot/Slide Tube w/ Aluminum Footplate Assembly, Chrome - Right	1	1
1	B	1023906-NLA	NLA - Pivot/Slide Tube w/ Aluminum Footplate Assembly, Chrome - Left	1	1
1	C	1023907	Footplate w/Pivot Slide Tube - RH (18")	1	1
1	C	1023908	Footplate w/Pivot Slide Tube - LH (18")	1	1
1	D	1023909-NLA	NLA - Pivot/Slide Tube w/ Composite Footplate Assembly, Chrome - Right	1	1
1	D	1023910-NLA	NLA - Pivot/Slide Tube w/ Composite Footplate Assembly, Chrome - Left	1	1
2	A	1041526	Pivot/Slide Tube w/ Plug Button Assembly, Chrome	1	1
3	B	44200X059	Footplate, Aluminum, Small, Black (6-3/4" W x 6" L)	1	1
3	C	44200X027	Footplate, Aluminum, Medium, Black (7-3/4" W x 6" L)	1	1
3	D	1015178	Footplate, Composite, Black (6" W x 7-53/64" L)	1	1
4		1018133	Plug Button, Black (5/8")	1	1
5		1027097	Plug Button, Black (3/4")	1	1
6		1006981	Spring, Tension	1	1
7		40120X001	Screw, Phillips Pan Head (#8-32 x 3/8")	1	1
8		1016346	Cap, End, Tube, Black (7/8")	1	1

NOTE: A - Includes plug buttons and pivot slide tube.
 B - For use on AHS4A models.
 C - For use on AHL4A models.
 D - For use on AHL4C models.