

Front Riggings Reference Chart

PHAL4A

NOTE: Calf pads with attaching hardware not included on legrest.

NOTE: Pivot/ Slide tube attaching hardware not included.

NOTE: Pivot/ Slide tube and footplates do not include heel loops, see Footrest and Legrest Accessories page.

1650 = Aluminum footplate with "L" style pivot/slide tube.
1350 = Extra large footplate with "L" style pivot/slide tube.
1450 = Deluxe aluminum footplate with "L" style pivot/slide tube

Support Assembly For PHAL4A

Rigging Type	PHAL4A (Chrome Only)
Right	1031437-NLA
Left	1031438-NLA

Pivot/Slide Tube and Footplates For PHAL4A

Footplate Type	1650 16-22" Wide	1350 17-22" Wide	1450 16-22" Wide
Right	1023907	1035456	1023924-NLA See Note(s) A
Left	1023908	1035457	1023925-NLA See Note(s) A

NOTE: A - Option for the Pronto Series only