

Invacare® Colibri™ Scooter - Canada



Yes, you can.

Parts Catalogue Usage Guide

The information contained in this document is subject to change without notice. This catalog contains service parts for Form No. 16-044C. If you have any questions regarding this catalog, please contact one of our Customer Service Representatives.

For Customer Service Support Call:

1-800-668-5324 or 1-905-890-8300

For Faxed In Orders Call:

1-800-668-5478 or

1-905-890-5244

For Technical Service Call:

1-800-668-5324 x2655

For Orders:

canadaorders@invacare.com

Literature: www.invacare.ca or 1-800-668-5324

SYMBOL DEFINITION

Throughout this catalog you will notice letters to draw your attention to specific footnotes which are located at the bottom of the page or below specified charts.

Examples:

A,B,C etc. = Special Comments

NOTE: = General information pertaining to that specific page

NLA = No Longer Available

HARDWARE MEASUREMENTS

The measurements for the hardware will be listed as follows:

Screws, Bolts Diameter X Thread Count/Pitch X Length

Washers Inside Diameter X Outside Diameter X Thickness

Spacers Inside Diameter X Outside Diameter X Length

For your convenience, we have supplied an area for you to list your Invacare Representatives.


Name	Phone Numbers

© 2021 Invacare Corporation. All rights reserved. Republication, duplication or modification in whole or in part is prohibited without prior written permission from Invacare. Trademarks are identified by ™ and ®. All trademarks are owned by or licensed to Invacare Corporation or its subsidiaries unless otherwise noted.

TABLE of CONTENTS

Description	Page Number
01. Front Chassis and Front Axis	3
02. Front Shroud	5
03. Colour Shrouds	7
04. Rear Chassis	9
05. Rear Shrouds	11
06. Steering and Tiller Shrouds	13
07. Electronics and Controller	15
08. Seat and Armpads	17
09. Front Wheels	19
10. Rear Wheels	21
11. Motor and Transaxle	23
12. Batteries, Battery Box, and Chargers	25
13. Cables	27
14. Accessories	29

Front Chassis and Front Axis




SCP_031A

Front Chassis and Front Axis

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1577405	Front Frame	1	1
2		1577424	Front Axis, Complete	1	1
3		1577517	Kit, Bearing, Front Axis, Complete	1	1
4		1577539	Swing Assembly	1	1
		1577512	Steering Fork (3W)	1	1
		1577514	Front Frame (3W)	1	1

Front Shroud




SCP_032A

Front Shroud

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1577419	Front Shroud	1	1
1		1577515	Front Shroud (3W)	1	1
2		1577422	Kit, Front Shroud Mounting	1	1

Colour Shrouds

Colour Shrouds




SCP_034A

Colour Shrouds

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1577251	Shroud Kit, Jasper Red	1	1
2		1577253	Shroud Kit, Sapphire Blue	1	1
3		1577252	Shroud Kit, Jade Green	1	1
4		1577250	Shroud Kit, Marble White	1	1
5		1577254	Shroud Kit, Sunstone Orange	1	1

Rear Chassis

Rear Chassis




SCP_035A

Rear Chassis


Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1577434	Rear Frame	1	1
2		1577450	Anti-Tipper	1	1

Rear Shrouds




Rear Shrouds



See “Colour Shrouds”



See “Battery Boxes”




See “Colour Shrouds”

SCP_036A

Rear Shrouds

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1577436	Rear Shroud	1	1
2		1577435	Motor Protection Shroud	1	1


Steering and Tiller Shrouds



Steering and Tiller Shrouds

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1577432	Steering Bar	1	1
2		1577433	Steering Pivot	1	1
3		1577418	Steering Tube	1	1
4		1577417	Steering Bearing Kit	1	1
5		1577427	Tiller Shrouds Kit	1	1
6		1505401	Cover, Steering	1	1
7		1577429	Handle Bar Grip	1	1

Electronics and Controller




Electronics and Controller

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1577423	Controller, DR50 w/ Mounting Hardware	1	1
2		1577428	Wig Wag w/ Pot	1	1
3		1577426	Tiller, Electronic (Complete)	1	1
		1530212	Tiller Keys (Set of 2)	1	1

Seat and Armpads

Seats and Armpads





SCP_038B

Seat and Armpads

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1577444	Seat	1	1
2		1577445	Seat Support	1	1
3		1577446	Seat Post	1	1
4		1577443	Arm Pad (Left/Right)	1	1
5		1507005	Knob, Tightening	1	1

Front Wheels

Front Wheels




SCP_039A

Front Wheels

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1577425	Front Wheel (200 x 50mm), Complete	1	2
2		1577452	Front Wheel (210 x 65mm), Complete	1	2
		1577513	Front Wheel, 3W (200 x 55mm), Complete	1	1

Rear Wheels

Rear Wheels




SCP_040A

Rear Wheels

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1577448	Rear Wheel (200 x 50mm), Complete	1	2
2		1577451	Rear Wheel (210 x 65mm), Complete	1	2

Motor and Transaxle

Motor And Transaxle



SCP_041A


Motor and Transaxle

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1577447	Transaxle, Complete	1	1
2		1577449	Motor w/ Mounting Hardware	1	1
3		1577454	Key Kit (10 pcs.)	1	1


Batteries, Battery Box, and Chargers

Replace Batteries in pairs and with same type only. Complete Battery Boxes also require cables listed in the Category "Cables".


Batteries, Battery Boxes, and Chargers



12Ah



18Ah



SCP_042A


Batteries, Battery Box, and Chargers

Replace Batteries in pairs and with same type only. Complete Battery Boxes also require cables listed in the Category "Cables".

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1577440	Battery Box	1	1
2	A	1577439	Battery Lid, 18Ah, Complete	1	1
3	B	1577438	Battery Lid, 12Ah, Complete	1	1
4	A	1577516	Mounting Kit, 18Ah	1	1
5	B	1577518	Mounting Kit, 12Ah	1	1
6	A	1577442	Battery, 12v/18Ah, AGM	1	2
7	B	12AH-2	Batteries, 12 volt -12 Amp Hour (Pair)	1	1
8	A	1148830	Charger, Off-Board, 24 DC (5 AMP)	1	1
9	B	1577437	Charger, 2A	1	1
NOTE: A - 18AH Battery Box B - 12AH Battery Box Replace Batteries in pairs and with same type only. Complete Battery Boxes also require cables listed in the Category "Cables".					

Cables

Cables




SCP_043A

Cables

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1577421	Battery Connector	1	1
2		1577420	Motor/Battery Connector	1	1
3		1577431	Wire Harness	1	1
4		1577521	Circuit Breaker Cable (12Ah)	1	1
5		1577522	Circuit Breaker Cable (18Ah)	1	1
6		1577519	Battery Cable (12Ah)	1	2
7		1577520	Battery Cable (18Ah)	1	2

Accessories

Accessories



Accessories

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1577455	Rear Accesory Bracket	1	1
2	A	ACC220	Rear Oxygen Bottle Holder	1	1
3		1577430	Front Composite Basket	1	1
4	A	ACC230	Rear Basket	1	1
5	A	ACC200	Rear Forearm Crutch and Cane Holder	1	1
6	A	ACC140	Safety Flag, Rear (6')	1	1
7		1530212	Tiller Keys (Set of 2)	1	1
8		PRC571840	Storage Cover	1	1
9		1577212	Seat Suspension, Low (425mm)	1	1
10		1577213	Seat Suspension, Low (450mm)	1	1
11		1577214	Seat Suspension, Low (475mm)	1	1
NOTE: A - Requires 1577455, Rear Accessory Bracket.					

Invacare Canada LP

www.invacare.ca



570 Matheson Blvd E;8
Mississauga, Ontario,
L4Z 4G4 Canada
800-668-5324

Technical Support
8006685324 x2655
Customer Service
800-668-5324

Literature
800-668-5324
or www.invacare.ca

© 2021 Invacare Corporation. All rights reserved. Republication, duplication or modification in whole or in part is prohibited without prior written permission from Invacare. Trademarks are identified by ™ and ®. All trademarks are owned by or licensed to Invacare Corporation or its subsidiaries unless otherwise noted. (Revised 06/15/21)

Form No. 16-044C