

**Invacare® Platinum™  
Mobile Oxygen Concentrator  
POC1-100B, POC1-100C**



*Yes, you can.*

## Parts Catalog Usage Guide

The information contained in this document is subject to change without notice. This catalog contains service parts for Form No. 16-201. If you have any questions regarding this catalog, please contact one of our Customer Service Representatives.

### For Customer Service Support Call:

**1-800-333-6900**

### For Faxed In Orders Call:

**1-800-678-4682 or**

**1-440-366-0709**

### For Technical Service Call:

**1-800-832-4707**

### Literature Can Be Found at:

**www.invacare.com or**

**1-800-828-6272**

## SYMBOL DEFINITION

Throughout this catalog you will notice letters to draw your attention to specific footnotes which are located at the bottom of the page or below specified charts.

Examples:

A,B,C etc. = Special Comments

NOTE: = General information pertaining to that specific page

NLA = No Longer Available

## HARDWARE MEASUREMENTS

The measurements for the hardware will be listed as follows:

Screws, Bolts Diameter X Thread Count/Pitch X Length

Washers Inside Diameter X Outside Diameter X Thickness

Spacers Inside Diameter X Outside Diameter X Length

For your convenience, we have supplied an area for you to list your Invacare Representatives.



Name	Phone Numbers

© 2019 Invacare Corporation. All rights reserved. Republication, duplication or modification in whole or in part is prohibited without prior written permission from Invacare. Trademarks are identified by <sup>TM</sup> and ®. All trademarks are owned by or licensed to Invacare Corporation or its subsidiaries unless otherwise noted.

# TABLE of CONTENTS

Description	Page Number
POC1-100B, POC1-100C Concentrator .....	3

# POC1-100B, POC1-100C Concentrator



SRE\_018D

## POC1-100B, POC1-100C Concentrator

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1195090	Kit, Pad, Non-Skid	10	2
2		1195092	Kit, Battery Ring Pull (Pkg. of 2)	2	1
3		1195095	Kit, Filter, Patient Outlet (Pkg. of 5)	5	1
4	A	1195091	Kit, Shoulder Strap	2	2
5		1187454	Power Cord, Line, 2P, North America	1	1
6		1152805	Power Cord, Line, 2P, Europe	1	1
7		1160791	EU - Power Cord, Line, 2P, United Kingdom	1	1
8		1160792	Power Cord, Line, 2P, Australia	1	1
9		1195098	Fuse (12 Amp)	1	1
10	B	POC1-130	Adapter, AC w/ Power Cord (North America)	1	1
10	C	POC1-130-EU	Adapter, AC w/ Power Cord (Europe)	1	1
10	D	POC1-130-AZ	Adapter, AC w/ Power Cord (Australia/New Zealand)	1	1
11		POC1-140	DC Power Cord	1	1
12	E	POC1-110	Battery Pack	1	1
13	F	POC1-150	Carry Bag Assembly	1	1
14	G	60104617	Kit, USB Dongle, POC1-CONNECT	1	1
		See Chart	Literature, Labels, and Packaging	1	1

NOTE: A - Two straps are used when you carry the Concentrator in Backpack Mode.  
 B - Kit includes the Adapter and USA Power Cord.  
 C - Kit includes the Adapter, EU, and UK Power Cords.  
 D - Kit includes the Adapter and AZ Power Cord.  
 E - Kit includes Battery Pull Ring and Battery.  
 F - Kit includes two Straps for Messenger or Backpack Modes.  
 G - USB Dongle Kit is only available for Platinum Mobile Oxygen Concentrators beginning with Serial #17JF029716.

### Literature, Labels, and Packaging

**NOTE: Manuals are not offered through Invacare service parts. Refer to the Invacare Website under "Product Literature Search" and download a free PDF version copy.**

Region	User Manual Concentrator	Warning Label	User Manual Battery Charger	Carton w/ Inserts	Over Laminate Label
English (North American)	1194969	1187525 See Note(s) A	1187396	1195096	1168369
Canada	1187376	1187564	1193862	1195096	1168369

NOTE: A - Warning Label for POC1-100C is currently not available.

Invacare Corporation

[www.invacare.com](http://www.invacare.com)



***Yes, you can.***

One Invacare Way  
Elyria, Oh USA  
44035  
800-333-6900

Technical Support  
800-832-4707  
Customer Service  
800-333-6900

Literature Orders  
800-333-6900

© 2019 Invacare Corporation. All rights reserved. Republication, duplication or modification in whole or in part is prohibited without prior written permission from Invacare. Trademarks are identified by ™ and ®. All trademarks are owned by or licensed to Invacare Corporation or its subsidiaries unless otherwise noted.  
(Revised 04/23/18)

Form No. 16-201