

Küschall® Champion Wheelchair



Yes, you can.

Parts Catalog Usage Guide

The information contained in this document is subject to change without notice. This catalog contains service parts for Form No. 16-249. If you have any questions regarding this catalog, please contact one of our Customer Service Representatives.

**For Customer Service Support Call:
844-MY KUSCH (695-8724)**

**For Faxed In Orders Call:
866-577-3507**

**For Technical Service Call:
800-832-4707**

**Literature Can Be Found at:
www.invacare.com**

SYMBOL DEFINITION

Throughout this catalog you will notice letters to draw your attention to specific footnotes which are located at the bottom of the page or below specified charts.

Examples:

A,B,C etc. = Special Comments
NOTE: = General information pertaining to that specific page
NLA = No Longer Available

HARDWARE MEASUREMENTS

The measurements for the hardware will be listed as follows:

Screws, Bolts Diameter X Thread Count/Pitch X Length
Washers Inside Diameter X Outside Diameter X Thickness
Spacers Inside Diameter X Outside Diameter X Length

For your convenience, we have supplied an area for you to list your Invacare Representatives.

Name	Phone Numbers

© 2022 Invacare Corporation. All rights reserved. Republication, duplication or modification in whole or in part is prohibited without prior written permission from Invacare. Trademarks are identified by ™ and ®. All trademarks are owned by or licensed to Invacare Corporation or its subsidiaries unless otherwise noted.

TABLE of CONTENTS

Description	Page Number
Abbreviations - Wheelchair Dimensions	5
Appearance - Frame Colors and Color Power Set (CPS)	6
Backrest - Adjustable Bands	51
Backrest - Adjustable Height Push Handles	42
Backrest - Adjustable Tension Cover	53
Backrest - Matrix Back	46
Backrest - Premium Cover	50
Backrest - Rear Set Adjustable Height Push Handles	44
Backrest - Stabilization Bar	54
Backrest - Standard Cover	48
Backrest - Standard Push Handles	40
Chassis - Camber Hardware	33
Chassis - Folding Unit	27
Chassis - Frame - Adducted (75 degree)	7
Chassis - Frame - Adducted (85 degree)	12
Chassis - Frame - Painted - Standard (75 Degree)	17
Chassis - Frame - Painted - Standard (85 Degree)	22
Chassis - Rear Wheel Mounting	31
Chassis - Side Support	29
Footrest - Adducted Frame	60
Footrest - Adducted Frame - High Mounted	62

TABLE of CONTENTS

Description	Page Number
Footrest - Calf Straps	64
Footrest - Standard Frame	56
Footrest - Standard Frame - High Mounted	58
Front Wheels - Casters	89
Front Wheels - Forks - Aluminum and Titanium Frame	76
Front Wheels - Forks - Carbon Frame	86
Options and Accessories - Anti-Tipper	123
Options and Accessories - Cane Holder	126
Options and Accessories - Chair Accessories	129
Options and Accessories - Luggage Holder	127
Options and Accessories - Tipper Aids	125
Options and Accessories - Transit Wheels	124
Options and Accessories - Transport Bag	128
Rear Wheels - 22 Inch Wheels	92
Rear Wheels - 24 and 25 Inch Wheels - High Performance	100
Rear Wheels - 24 Inch Wheels	94
Rear Wheels - 26 Inch Wheels	97
Rear Wheels - Handrim Hardware	109
Rear Wheels - Handrims	106
Rear Wheels - Handrims - Carbolife	108
Rear Wheels - Inner Tube and Rim Protection	113

TABLE of CONTENTS

Description	Page Number
Rear Wheels - Quick Release Axles and Bearings	117
Rear Wheels - Spinergy LXX - 24" and 25" Wheels	104
Rear Wheels - Spoke Guards	115
Rear Wheels - Tires	111
Rear Wheels - Wheel Locks	119
Seat Accessories - Fixing Clips	122
Seat Accessories - Seat Bag	121
Seating - Light Seat Covers	37
Seating - Standard Seat Covers	35
Seating - Standard Seat Cushion	39
Siderests - Adjustable Arm	72
Siderests - Clothes Guard	66
Siderests - Fixed Mudguards	68
Siderests - Removable Mudguards	70
Siderests - Tubular Armrest	74


Abbreviations - Wheelchair Dimensions

Dimension Key:

XX = Dimensions in CM


XXX = Dimensions in MM

XX,A = Metric Decimal




SB (SEAT WIDTH)


ST (SEAT DEPTH)




RH (BACK HEIGHT)




UL (KNEE TO HEEL DEPTH)




SHv (FRONT SEAT TO FLOOR)



SHh (REAR SEAT TO FLOOR)



SHD (FRONT TO REAR DIFFERENCE)



SCM_290A

Appearance - Frame Colors and Color Power Set (CPS)

-0184	-0450	-0001
-0185	-0451	-0182
-0186	-0457	-0400
-0187	-3524	-0404
-0411	-5969	-0440
-0419	-9005	-0441
-0445		

Frame Color Code

Color Power Set (CPS) Code

Applies To The Following Parts And Options


- Starec and High Performance Rear Wheel Hubs w/ Quick Release Protection.
- 22" and 24" Starec Rear Wheel Rims (Except Chrome).
- Starec Front Wheels, Side Supporters.
- Adapter Plates.
- Forks.

CHAMPION COLORS

Appearance - Frame Colors and Color Power Set (CPS)

Frame Color Code Number	Frame Color	Color Power Set (CPS) Code Number	Color Power Set (CPS) Color Description
-0184	Mercury Blue	-0001	Chrome
-0185	Dynamite Orange	-0182	Champagne
-0186	Atomic Violet	-0400	Techno-Anthracite
-0187	Nitro Pink	-0404	B-Blue
-0411	Matte Black	-0440	Comet Red
-0419	Silver Grey	-0441	Arctic Blue
-0445	Ivory White	N/A	N/A
-0450	Smokey Grey	N/A	N/A
-0451	Flashy Yellow	N/A	N/A
-0457	Mint Green	N/A	N/A
-3524	Royal Red	N/A	N/A
-5969	Royal Blue	N/A	N/A
-9005	Glossy Black	N/A	N/A

Frame - Adducted (75 degree)



SCM_236A

Frame - Adducted (75 degree)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1-4		See Chart	Frame - Adducted (75 degree, ST340 - ST400) (Pair)	2	1
5-8		See Chart	Frame - Adducted (75 degree, ST420 - ST480) (Pair)	2	1
9-12		See Chart	Frame - Adducted (75 degree, Long ST340 - ST400) (Pair)	2	1
13-16		See Chart	Frame - Adducted (75 degree, Long ST420 - ST480) (Pair)	2	1
17		1529372	Front Frame Protectors (Pair)	2	1

ST340 - ST400 (Pair)

Color	Item #1 ST340	Item #2 ST360	Item #3 ST380	Item #4 ST400
Mercury Blue	1587010-0184	1587011-0184	1587012-0184	1587013-0184
Dynamite Orange	1587010-0185	1587011-0185	1587012-0185	1587013-0185
Atomic Violet	1587010-0186	1587011-0186	1587012-0186	1587013-0186
Nitro Pink	1587010-0187	1587011-0187	1587012-0187	1587013-0187
Matte Black	1587010-0411	1587011-0411	1587012-0411	1587013-0411
Silver Grey	1587010-0419	1587011-0419	1587012-0419	1587013-0419
Ivory White	1587010-0445	1587011-0445	1587012-0445	1587013-0445
Smokey Grey	1587010-0450	1587011-0450	1587012-0450	1587013-0450
Flashy Yellow	1587010-0451	1587011-0451	1587012-0451	1587013-0451
Mint Green	1587010-0457	1587011-0457	1587012-0457	1587013-0457
Royal Red	1587010-3524	1587011-3524	1587012-3524	1587013-3524
Royal Blue	1587010-5969	1587011-5969	1587012-5969	1587013-5969
Glossy Black	1587010-9005	1587011-9005	1587012-9005	1587013-9005

Frame - Adducted (75 degree)

ST420 - ST480 (Pair)

Color	Item #5 ST420	Item #6 ST440	Item #7 ST460	Item #8 ST480
Mercury Blue	1587014-0184	1587015-0184	1587016-0184	1587017-0184
Dynamite Orange	1587014-0185	1587015-0185	1587016-0185	1587017-0185
Atomic Violet	1587014-0186	1587015-0186	1587016-0186	1587017-0186
Nitro Pink	1587014-0187	1587015-0187	1587016-0187	1587017-0187
Matte Black	1587014-0411	1587015-0411	1587016-0411	1587017-0411
Silver Grey	1587014-0419	1587015-0419	1587016-0419	1587017-0419
Ivory White	1587014-0445	1587015-0445	1587016-0445	1587017-0445
Smokey Grey	1587014-0450	1587015-0450	1587016-0450	1587017-0450
Flashy Yellow	1587014-0451	1587015-0184	1587016-0451	1587017-0451
Mint Green	1587014-0457	1587015-0185	1587016-0457	1587017-0457
Royal Red	1587014-3524	1587015-3524	1587016-3524	1587017-3524
Royal Blue	1587014-5969	1587015-5969	1587016-5969	1587017-5969
Glossy Black	1587014-9005	1587015-9005	1587016-9005	1587017-9005

Frame - Adducted (75 degree)

Long ST340 - ST400 (Pair)


Color	Item #9 ST340	Item #10 ST360	Item #11 ST380	Item #12 ST400
Mercury Blue	1587076-0184	1587077-0184	1587078-0184	1587079-0184
Dynamite Orange	1587076-0185	1587077-0185	1587078-0185	1587079-0185
Atomic Violet	1587076-0186	1587077-0186	1587078-0186	1587079-0186
Nitro Pink	1587076-0187	1587077-0187	1587078-0187	1587079-0187
Matte Black	1587076-0411	1587077-0411	1587078-0411	1587079-0411
Silver Grey	1587076-0419	1587077-0419	1587078-0419	1587079-0419
Ivory White	1587076-0445	1587077-0445	1587078-0445	1587079-0445
Smokey Grey	1587076-0450	1587077-0450	1587078-0450	1587079-0450
Flashy Yellow	1587076-0451	1587077-0451	1587078-0451	1587079-0451
Mint Green	1587076-0457	1587077-0457	1587078-0457	1587079-0457
Royal Red	1587076-3524	1587077-3524	1587078-3524	1587079-3524
Royal Blue	1587076-5969	1587077-5969	1587078-5969	1587079-5969
Glossy Black	1587076-9005	1587077-9005	1587078-9005	1587079-9005

Frame - Adducted (75 degree)

Long ST420 - ST480 (Pair)

Color	Item #13 ST420	Item #14 ST440	Item #15 ST460	Item #16 ST480
Mercury Blue	1587080-0184	1587081-0184	1587082-0184	1587083-0184
Dynamite Orange	1587080-0185	1587081-0185	1587082-0185	1587083-0185
Atomic Violet	1587080-0186	1587081-0186	1587082-0186	1587083-0186
Nitro Pink	1587080-0187	1587081-0187	1587082-0187	1587083-0187
Matte Black	1587080-0411	1587081-0411	1587082-0411	1587083-0411
Silver Grey	1587080-0419	1587081-0419	1587082-0419	1587083-0419
Ivory White	1587080-0445	1587081-0445	1587082-0445	1587083-0445
Smokey Grey	1587080-0450	1587081-0450	1587082-0450	1587083-0450
Flashy Yellow	1587080-0451	1587081-0451	1587082-0451	1587083-0451
Mint Green	1587080-0457	1587081-0457	1587082-0457	1587083-0457
Royal Red	1587080-3524	1587081-3524	1587082-3524	1587083-3524
Royal Blue	1587080-5969	1587081-5969	1587082-5969	1587083-5969
Glossy Black	1587080-9005	1587081-9005	1587082-9005	1587083-9005

Frame - Adducted (85 degree)



SCM_234A

Frame - Adducted (85 degree)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1-4		See Chart	Frame - Adducted (85 degree (ST340 - ST400) (Pair)	2	1
5-8		See Chart	Frame - Adducted (85 degree (ST420 - ST480) (Pair)	2	1
9-12		See Chart	Frame - Adducted (85 degree, Long ST340 - ST400) (Pair)	2	1
13-16		See Chart	Frame - Adducted (85 degree, Long ST420 - ST480) (Pair)	2	1
17		1529372	Front Frame Protectors (Pair)	2	1

Adducted (85 degree (ST340 - ST400) (Pair)

Color	Item #1 ST340	Item #2 ST360	Item #3 ST380	Item #4 ST400
Mercury Blue	1587100-0184	1587101-0184	1587102-0184	1587103-0184
Dynamite Orange	1587100-0185	1587101-0185	1587102-0185	1587103-0185
Atomic Violet	1587100-0186	1587101-0186	1587102-0186	1587013-0186
Nitro Pink	1587100-0187	1587101-0187	1587102-0187	1587103-0187
Matte Black	1587100-0411	1587101-0411	1587102-0411	1587103-0411
Silver Grey	1587100-0419	1587101-0419	1587102-0419	1587103-0419
Ivory White	1587100-0445	1587101-0445	1587102-0445	1587103-0445
Smokey Grey	1587100-0450	1587101-0450	1587102-0450	1587103-0450
Flashy Yellow	1587100-0451	1587101-0451	1587102-0451	1587103-0451
Mint Green	1587100-0457	1587101-0457	1587102-0457	1587103-0457
Royal Red	1587100-3524	1587101-3524	1587102-3524	1587103-3524
Royal Blue	1587100-5969	1587101-5969	1587102-5969	1587103-5969
Glossy Black	1587100-9005	1587101-9005	1587102-9005	1587103-9005

Frame - Adducted (85 degree)

ST420 - ST480 (Pair)

Color	Item #5 ST420	Item #6 ST440	Item #7 ST460	Item #8 ST480
Mercury Blue	1587104-0184	1587105-0184	1587106-0184	1587107-0184
Dynamite Orange	1587104-0185	1587105-0185	1587106-0185	1587107-0185
Atomic Violet	1587104-0186	1587105-0186	1587106-0186	1587107-0186
Nitro Pink	1587104-0187	1587105-0187	1587106-0187	1587107-0187
Matte Black	1587104-0411	1587105-0411	1587106-0411	1587107-0411
Silver Grey	1587104-0419	1587105-0419	1587106-0419	1587107-0419
Ivory White	1587104-0445	1587105-0445	1587106-0445	1587107-0445
Smokey Grey	1587104-0450	1587105-0450	1587106-0450	1587107-0450
Flashy Yellow	1587104-0451	1587105-0451	1587106-0451	1587107-0451
Mint Green	1587104-0457	1587105-0457	1587106-0457	1587107-0457
Royal Red	1587104-3524	1587105-3524	1587106-3524	1587107-3524
Royal Blue	1587104-5969	1587105-5969	1587106-5969	1587107-5969
Glossy Black	1587104-9005	1587105-9005	1587106-9005	1587107-9005

Frame - Adducted (85 degree)

Adducted (85 degree, Long ST340 - ST400) (Pair)


Color	Item #9 ST340	Item #10 ST360	Item #11 ST380	Item #12 ST400
Mercury Blue	1587145-0184	1587146-0184	1587147-0184	1587148-0184
Dynamite Orange	1587145-0185	1587146-0185	1587147-0185	1587148-0185
Atomic Violet	1587145-0186	1587146-0186	1587147-0186	1587148-0186
Nitro Pink	1587145-0187	1587146-0187	1587147-0187	1587148-0187
Matte Black	1587145-0411	1587146-0411	1587147-0411	1587148-0411
Silver Grey	1587145-0419	1587146-0419	1587147-0419	1587148-0419
Ivory White	1587145-0445	1587146-0445	1587147-0445	1587148-0445
Smokey Grey	1587145-0450	1587146-0450	1587147-0450	1587148-0450
Flashy Yellow	1587145-0451	1587146-0451	1587147-0451	1587148-0451
Mint Green	1587145-0457	1587146-0457	1587147-0457	1587148-0457
Royal Red	1587145-3524	1587146-3524	1587147-3524	1587148-3524
Royal Blue	1587145-5969	1587146-5969	1587147-5969	1587148-5969
Glossy Black	1587145-9005	1587146-9005	1587147-9005	1587148-9005

Frame - Adducted (85 degree)

Adducted (85 degree, Long ST420 - ST480) (Pair)

Color	Item #13 ST420	Item #14 ST440	Item #15 ST460	Item #16 ST480
Mercury Blue	1587149-0184	1587150-0184	1587151-0184	1587152-0184
Dynamite Orange	1587149-0185	1587150-0185	1587151-0185	1587152-0185
Atomic Violet	1587149-0186	1587150-0186	1587151-0186	1587152-0186
Nitro Pink	1587149-0187	1587150-0187	1587151-0187	1587152-0187
Matte Black	1587149-0411	1587150-0411	1587151-0411	1587152-0411
Silver Grey	1587149-0419	1587150-0419	1587151-0419	1587152-0419
Ivory White	1587149-0445	1587150-0445	1587151-0445	1587152-0445
Smokey Grey	1587149-0450	1587150-0450	1587151-0450	1587152-0450
Flashy Yellow	1587149-0451	1587150-0451	1587151-0451	1587152-0451
Mint Green	1587149-0457	1587150-0457	1587151-0457	1587152-0457
Royal Red	1587149-3524	1587150-3524	1587151-3524	1587152-3524
Royal Blue	1587149-5969	1587150-5969	1587151-5969	1587152-5969
Glossy Black	1587149-9005	1587150-9005	1587151-9005	1587152-9005

Frame - Painted - Standard (75 Degree)



SCM_237A

Frame - Painted - Standard (75 Degree)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1-4		See Chart	Frame - Standard (75 Degree - Short (ST340-ST400) (Pair)	2	1
5-7		See Chart	Frame - Standard (75 Degree - Short (ST420-ST460) (Pair)	2	1
8-11		See Chart	Frame - Standard (75 Degree - Long (ST360-ST400) (Pair)	2	1
12-14		See Chart	Frame - Standard (75 Degree - Long (ST420-ST460) (Pair)	2	1
15		1529372	Front Frame Protectors (Pair)	2	1

Standard (75 Degree - Short (ST340-ST400) (Pair)

Color	Item #1 ST340	Item #2 ST360	Item #3 ST380	Item #4 ST400
Mercury Blue	1585854-0184	1585855-0184	1585856-0184	1585857-0184
Dynamite Orange	1585854-0185	1585855-0185	1585856-0185	1585857-0185
Atomic Violet	1585854-0186	1585855-0186	1585856-0186	1585857-0186
Nitro Pink	1585854-0187	1585855-0187	1585856-0187	1585857-0187
Matte Black	1585854-0411	1585855-0411	1585856-0411	1585857-0411
Silver Grey	1585854-0419	1585855-0419	1585856-0419	1585857-0419
Ivory White	1585854-0445	1585855-0445	1585856-0445	1585857-0445
Smokey Grey	1585854-0450	1585855-0450	1585856-0450	1585857-0450
Flashy Yellow	1585854-0451	1585855-0451	1585856-0451	1585857-0451
Mint Green	1585854-0457	1585855-0457	1585856-0457	1585857-0457
Royal Red	1585854-3524	1585855-3524	1585856-3524	1585857-3524
Royal Blue	1585854-5969	1585855-5969	1585856-5969	1585857-5969
Glossy Black	1585854-9005	1585855-9005	1585856-9005	1585857-9005

Frame - Painted - Standard (75 Degree)

Standard (75 Degree - Short (ST420-ST460) (Pair)

Color	Item #5 ST420	Item #6 ST440	Item #7 ST460
Mercury Blue	1585858-0184	1585859-0184	1585860-0184
Dynamite Orange	1585858-0185	1585859-0185	1585860-0185
Atomic Violet	1585858-0186	1585859-0186	1585860-0186
Nitro Pink	1585858-0187	1585859-0187	1585860-0187
Matte Black	1585858-0411	1585859-0411	1585860-0411
Silver Grey	1585858-0419	1585859-0419	1585860-0419
Ivory White	1585858-0445	1585859-0445	1585860-0445
Smokey Grey	1585858-0450	1585859-0450	1585860-0450
Flashy Yellow	1585858-0451	1585859-0451	1585860-0451
Mint Green	1585858-0457	1585859-0457	1585860-0457
Royal Red	1585858-3524	1585859-3524	1585860-3524
Royal Blue	1585858-5969	1585859-5969	1585860-5969
Glossy Black	1585858-9005	1585859-9005	1585860-9005

Frame - Painted - Standard (75 Degree)

Standard (75 Degree - Long (ST360-ST400) (Pair)


Color	Item #8 ST340	Item #9 ST360	Item #10 ST380	Item #11 ST400
Mercury Blue	1585986-0184	1585987-0184	1585988-0184	1585989-0184
Dynamite Orange	1585986-0185	1585987-0185	1585988-0185	1585989-0185
Atomic Violet	1585986-0186	1585987-0186	1585988-0186	1585989-0186
Nitro Pink	1585986-0187	1585987-0187	1585988-0187	1585989-0187
Matte Black	1585986-0411	1585987-0411	1585988-0411	1585989-0411
Silver Grey	1585986-0419	1585987-0419	1585988-0419	1585989-0419
Ivory White	1585986-0445	1585987-0445	1585988-0445	1585989-0445
Smokey Grey	1585986-0450	1585987-0450	1585988-0450	1585989-0450
Flashy Yellow	1585986-0451	1585987-0451	1585988-0451	1585989-0451
Mint Green	1585986-0457	1585987-0457	1585988-0457	1585989-0457
Royal Red	1585986-3524	1585987-3524	1585988-3524	1585989-3524
Royal Blue	1585986-5969	1585987-5969	1585988-5969	1585989-5969
Glossy Black	1585986-9005	1585987-9005	1585988-9005	1585989-9005

Frame - Painted - Standard (75 Degree)

Standard (75 Degree - Long (ST420-ST460) (Pair)

Color	Item #12 ST420	Item #13 ST440	Item #14 ST460
Mercury Blue	1586087-0184	1586088-0184	1586089-0184
Dynamite Orange	1586087-0185	1586088-0185	1586089-0185
Atomic Violet	1586087-0186	1586088-0186	1586089-0186
Nitro Pink	1586087-0187	1586088-0187	1586089-0187
Matte Black	1586087-0411	1586088-0411	1586089-0411
Silver Grey	1586087-0419	1586088-0419	1586089-0419
Ivory White	1586087-0445	1586088-0445	1586089-0445
Smokey Grey	1586087-0450	1586088-0450	1586089-0450
Flashy Yellow	1586087-0451	1586088-0451	1586089-0451
Mint Green	1586087-0457	1586088-0457	1586089-0457
Royal Red	1586087-3524	1586088-3524	1586089-3524
Royal Blue	1586087-5969	1586088-5969	1586089-5969
Glossy Black	1586087-9005	1586088-9005	1586089-9005

Frame - Painted - Standard (85 Degree)



Frame - Painted - Standard (85 Degree)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1-4		See Chart	Frame - Standard (85 Degree - Long (ST340-ST400) (Pair)	2	1
5-7		See Chart	Frame - Standard (85 Degree - Short (ST420-ST460) (Pair)	2	1
8-11		See Chart	Frame - Standard (85 Degree - Long (ST340-ST400) (Pair)	2	1
12-14		See Chart	Frame - Standard (85 Degree - Long (ST420-ST460) (Pair)	2	1
15		1529372	Front Frame Protectors (Pair)	2	1

Standard (85 Degree - Short (ST340-ST400) (Pair)

Color	Item #1 ST340	Item #2 ST360	Item #3 ST380	Item #4 ST400
Mercury Blue	1586641-0184	1586642-0184	1586643-0184	1586644-0184
Dynamite Orange	1586641-0185	1586642-0185	1586643-0185	1586644-0185
Atomic Violet	1586641-0186	1586642-0186	1586643-0186	1586644-0186
Nitro Pink	1586641-0187	1586642-0187	1586643-0187	1586644-0187
Matte Black	1586641-0411	1586642-0411	1586643-0411	1586644-0411
Silver Grey	1586641-0419	1586642-0419	1586643-0419	1586644-0419
Ivory White	1586641-0445	1586642-0445	1586643-0445	1586644-0445
Smokey Grey	1586641-0450	1586642-0450	1586643-0450	1586644-0450
Flashy Yellow	1586641-0451	1586642-0451	1586643-0451	1586644-0451
Mint Green	1586641-0457	1586642-0457	1586643-0457	1586644-0457
Royal Red	1586641-3524	1586642-3524	1586643-3524	1586644-3524
Royal Blue	1586641-5969	1586642-5969	1586643-5969	1586644-5969
Glossy Black	1586641-9005	1586642-9005	1586643-9005	1586644-9005

Frame - Painted - Standard (85 Degree)

Standard (85 Degree - Short (ST420-ST460) (Pair)

Color	Item #5 ST420	Item #6 ST440	Item #7 ST460
Mercury Blue	1586645-0184	1586646-0184	1586647-0184
Dynamite Orange	1586645-0185	1586646-0185	1586647-0185
Atomic Violet	1586645-0186	1586646-0186	1586647-0186
Nitro Pink	1586645-0187	1586646-0187	1586647-0187
Matte Black	1586645-0411	1586646-0411	1586647-0411
Silver Grey	1586645-0419	1586646-0419	1586647-0419
Ivory White	1586645-0445	1586646-0445	1586647-0445
Smokey Grey	1586645-0450	1586646-0450	1586647-0450
Flashy Yellow	1586645-0451	1586646-0451	1586647-0451
Mint Green	1586645-0457	1586646-0457	1586647-0457
Royal Red	1586645-3524	1586646-3524	1586647-3524
Royal Blue	1586645-5969	1586646-5969	1586647-5969
Glossy Black	1586645-9005	1586646-9005	1586647-9005

Frame - Painted - Standard (85 Degree)

Standard (85 Degree (ST340-ST400) (Pair)


Color	Item #8 ST340	Item #9 ST360	Item #10 ST380	Item #11 ST400
Mercury Blue	1586700-0184	1586701-0184	1586702-0184	1586703-0184
Dynamite Orange	1586700-0185	1586701-0185	1586702-0185	1586703-0185
Atomic Violet	1586700-0186	1586701-0186	1586702-0186	1586703-0186
Nitro Pink	1586700-0187	1586701-0187	1586702-0187	1586703-0187
Matte Black	1586700-0411	1586701-0411	1586702-0411	1586703-0411
Silver Grey	1586700-0419	1586701-0419	1586702-0419	1586703-0419
Ivory White	1586700-0445	1586701-0445	1586702-0445	1586703-0445
Smokey Grey	1586700-0450	1586701-0450	1586702-0450	1586703-0450
Flashy Yellow	1586700-0451	1586701-0451	1586702-0451	1586703-0451
Mint Green	1586700-0457	1586701-0457	1586702-0457	1586703-0457
Royal Red	1586700-3524	1586701-3524	1586702-3524	1586703-3524
Royal Blue	1586700-5969	1586701-5969	1586702-5969	1586703-5969
Glossy Black	1586700-9005	1586701-9005	1586702-9005	1586703-9005

Frame - Painted - Standard (85 Degree)

Standard (85 Degree - Long (ST420-ST460) (Pair)

Color	Item #12 ST420	Item #13 ST440	Item #14 ST460
Mercury Blue	1586704-0184	1586705-0184	1586835-0184
Dynamite Orange	1586704-0185	1586705-0185	1586835-0185
Atomic Violet	1586704-0186	1586705-0186	1586835-0186
Nitro Pink	1586704-0187	1586705-0187	1586835-0187
Matte Black	1586704-0411	1586705-0411	1586835-0411
Silver Grey	1586704-0419	1586705-0419	1586835-0419
Ivory White	1586704-0445	1586705-0445	1586835-0445
Smokey Grey	1586704-0450	1586705-0450	1586835-0450
Flashy Yellow	1586704-0451	1586705-0451	1586835-0451
Mint Green	1586704-0457	1586705-0457	1586835-0457
Royal Red	1586704-3524	1586705-3524	1586835-3524
Royal Blue	1586704-5969	1586705-5969	1586835-5969
Glossy Black	1586704-9005	1586705-9005	1586835-9005

Chassis - Folding Unit




SCM_239A

Chassis - Folding Unit

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1588012	Folding Unit Complete, 364mm.	1	1
2	B	1588013	Folding Unit Complete, 384mm.	1	1
3	C	1588014	Folding Unit Complete, 404mm.	1	1
4	D	1588015	Folding Unit Complete, 424mm.	1	1
5	E	1588016	Folding Unit Complete, 444mm.	1	1
6	F	1588017	Folding Unit Complete, 464mm.	1	1
7	G	1588018	Folding Unit Complete, 484mm.	1	1
8		1588019	Kit, Folding Unit Sleeves	1	1
9		1588007	Kit, Pressure Part, Folding Unit	2	1
10		1588036	Kit, Unlock Handle Protector	1	1
11		1588026	Folding Aid	1	1
12		1588025	Defolding Aid	1	1
NOTE: A - Use with SB360. B - Use with SB380. C - Use with SB400. D - Use with SB420. E - Use with SB440. F - Use with SB460. G - Use with SB480.					

Chassis - Side Support



SCM_240A


Chassis - Side Support

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		See Chart	Chassis - Side Support	1	1
2		1588008	Kit, Hardware, Frame, Side Support	1	1
3		1588010	Vertical Brace, Complete (Pair)	1	1
4		1588011	Kit, Hardware, Vertical Brace	1	1

Chassis - Side Support

Color	Item 1 Pair
Chrome	1588024-0001
Black Anodized	1588024-0004
Champagne	1588024-0182
Techno-Anthrazit	1588024-0400
B-Blue	1588024-0404
Comet Red	1588024-0440
Arctic Blue	1588024-0441

Chassis - Rear Wheel Mounting




Chassis - Rear Wheel Mounting

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		See Chart	Chassis - Rear Wheel Mounting	1	1
2		See Chart	Chassis - Rear Wheel Mounting	1	1
3		1588004	Kit, Hardware, Adapter Plate	1	1
4		1588005	Adapter Plate Clamp Part, Folding Unit Pressure Part	2	1
5		1588006	Adapter Plate Clamp Part	4	1
6		1588007	Kit, Pressure Part, Folding Unit	2	1

Chassis - Rear Wheel Mounting

Color	Item 1 Right	Item 2 Left
Chrome	1588020-0001	1588021-0001
Black Anodized	1588020-0004	1588021-0004
Champagne	1588020-0182	1588021-0182
Techno-Anthrazit	1588020-0400	1588021-0400
B - Blue	1588020-0404	1588021-0404
Comet Red	1588020-0440	1588021-0440
Arctic Blue	1588020-0441	1588021-0441

Chassis - Camber Hardware




SCM_277A

Chassis - Camber Hardware

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A,C	SP1589278	Kit, 0° Adapter Sleeve, Head, Short (Pair)	2	1
2	B,D	SP1589279	Kit, 2° Adapter Sleeve, Head, Short (Pair)	2	1
3	E	SP1589280	Kit, 2° Adapter Sleeve, Head, Long (Pair)	2	1
4	F	SP1589281	Kit, 4° Adapter Sleeve, Head, Long (Pair)	2	1
5	G	SP1589282	Kit, 4° Adapter Sleeve, Head, Short (Pair)	2	1
6		SP1589284	Spacer, Adapter Sleeve (7mm) (Pair)	2	1
7		SP1589285	Spacer, Adapter Sleeve (3mm) (Pair)	2	1
NOTE: A - With Clothes Guard or Clothes Guard and Transit Wheel. B - With Clothes Guard, Clothes Guard and Transit Wheel, or Fixed Mudguard and Transit Wheel. C - With Removable Mudguard, Removable Mudguard and Transit Wheel. D - With Removable Mudguard, and Removable Mudguard and Transit Wheel. E - With Fixed Mudguard. F - With Clothes Guard, Removable Mudguard, and Fixed Mudguard. G - With Clothes Guard, Removable Mudguard, Fixed Mudguard, and Transit Wheel.					

Seating - Standard Seat Covers



Seating - Standard Seat Covers

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
		See Chart	Seating - Seat Covers (SB36-SB48) x (ST34-ST40)	1	1
		See Chart	Seating - Seat Covers (SB36-SB48) x (ST42-ST48)	1	1
8		1588192	Kit, Seat Rail - ST340 (Pair)	2	1
9		1588193	Kit, Seat Rail - ST360 (Pair)	2	1
10		1588194	Kit, Seat Rail - ST380 (Pair)	2	1
11		1588195	Kit, Seat Rail - ST400 (Pair)	2	1
12		1588196	Kit, Seat Rail - ST420 (Pair)	2	1
13		1588197	Kit, Seat Rail - ST440 (Pair)	2	1
14		1588198	Kit, Seat Rail - ST460 (Pair)	2	1
15		1588199	Kit, Seat Rail - ST480 (Pair)	2	1


Seat Covers (SB36-SB48) x (ST34-ST40)

Model	ST34	ST36	ST38	ST40
SB36	1581477	1581478	1581479	1581480
SB38	1581485	1581486	1581593	1581594
SB40	1581899	1581900	1581901	1581902
SB42	1581970	1581971	1581972	1581973
SB44	1581978	1581979	1581980	1581981
SB46	1581986	1581987	1581988	1581989
SB48	1581994	1581995	1581996	1582018

Seat Covers (SB36-SB48) x (ST42-ST48)

Model	ST42	ST44	ST46	ST48
SB36	1581481	1581482	1581483	1581484
SB38	1581718	1581748	1581749	1581750
SB40	1581903	1581967	1581968	1581969
SB42	1581974	1581975	1581976	1581977
SB44	1581982	1581983	1581984	1581985
SB46	1581990	1581991	1581992	1581993
SB48	1582019	1582020	1582021	1582022

Seating - Light Seat Covers



SCM_243A


Seating - Light Seat Covers

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
		See Chart	Seating - Light Seat Covers (SB36-SB42) x (ST36-ST42)	1	1
5		1588193	Kit, Seat Rail - ST360 (Pair)	2	1
6		1588194	Kit, Seat Rail - ST380 (Pair)	2	1
7		1588195	Kit, Seat Rail - ST400 (Pair)	2	1
8		1588196	Kit, Seat Rail - ST420 (Pair)	2	1

Light Seat Covers (SB36-SB42) x (ST36-ST42)


Model	ST36	ST38	ST40	ST42
SB36	1583419	1583420	1583421	1583422
SB38	1583423	1583424	1583425	1583426
SB40	1583427	1583428	1583429	1583430
SB42	1583431	1583432	1583433	1583434

Seating - Standard Seat Cushion



Sizes	SB36	SB38	SB40	SB42	SB44	SB46	SB48
ST34	1556326	1556327	1556328	1556329	1556330	1556331	1556332
ST36	1556335	1556336	1556337	1556338	1556339	1556340	1556341
ST38	1556336	1556344	1556345	1556346	1556347	1556348	1556349
ST40	1556337	1556345	1556352	1556353	1556354	1556355	1556356
ST42	1556338	1556346	1556353	1556359	1556360	1556361	1556362
ST44	1556339	1556347	1556354	1556360	1556365	1556366	1556367
ST46	1556340	1556348	1556355	1556361	1556366	1556370	1556371
ST48	1556341	1556349	1556356	1556362	1556367	1556371	1556374

Backrest - Standard Push Handles




SCM_233A

Backrest - Standard Push Handles

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1588221	Latch Bolt w/ Hardware , Clothes Guard	1	1
2		1588222	Latch Bolt w/ Hardware , Mudguard Fixed	1	1
3		1588223	Latch Bolt w/ Hardware , Mudguard Removable or Siderest	1	1
4		1588224	Kit, Hardware, Latch Bolt, Clothes Guard	1	1
5		1588225	Kit, Hardware, Latch Bolt, Mudguard Fixed	1	1
6		1588226	Kit, Hardware, Latch Bolt, Mudguard Removable or Siderest	1	1
7		1588218	Kit, Intermediate Tube, Short	1	1
8		1588219	Kit, Intermediate Tube, Long	1	1
9		1588220	Kit, Hardware, Intermediate Tube	1	1
10		1588217	Kit, Screw, Backrest Tube	1	1
11		1525122	Backrest Tube Standard, Short (Pair)	2	1
12		1525123	Backrest Tube Standard, Long (Pair)	2	1
13		1525124	Backrest Tube Mini, Short (Pair)	2	1
14		1525125	Backrest Tube Mini, Long (Pair)	2	1
15		1507734	Handle, Black (Pair)	2	1
16		1525126	Handle Black, Mini (Pair)	2	1
17		1524956	Telescopic Tube Short, Foldable Push Handle (Pair)	2	1
18		1525121	Telescopic Tube Medium, Foldable Push Handle (Pair)	2	1
19		1507733	Foldable Push Handles (Pair)	2	1
20		1546913	Telescopic Back Tube, Standard - Short (Pair)	2	1
21		1546914	Telescopic Back Tube, Standard - Medium (Pair)	2	1
22		1546946	Telescopic Back Tube, Standard - Long (Pair)	2	1
23		1432083	Kit, Hardware, Backrest Band	1	1

Backrest - Adjustable Height Push Handles




SCM_232A

Backrest - Adjustable Height Push Handles

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1438203	Kit, Adjustable Height, Integrated Push Handles, Standard, Short	1	1
2		1549985	Adjusting Lever, Adjustable Height, Integrated Push Handles (Pair)	1	1
3		1507734	Handle, Black (Pair)	2	1
4		1507752	Backrest Tube, Straight (265mm), Integrated Push Handles (Pair)	2	1
5		1507753	Backrest Tube, Straight (280mm), Integrated Push Handles (Pair)	2	1
6		1507862	Backrest Tube, Straight (295mm), Integrated Push Handles (Pair)	2	1
7		1508438	Backrest Tube, Straight (310mm), Integrated Push Handles (Pair)	2	1
8		1508466	Backrest Tube, Straight (325mm), Integrated Push Handles (Pair)	2	1
9		1508498	Backrest Tube, Straight (340mm), Integrated Push Handles (Pair)	2	1
10		1509103	Backrest Tube, Straight (355mm), Integrated Push Handles (Pair)	2	1
11		1509104	Backrest Tube, Straight (370mm), Integrated Push Handles (Pair)	2	1
12		1509105	Backrest Tube, Straight (385mm), Integrated Push Handles (Pair)	2	1
13		1509106	Backrest Tube, Straight (400mm), Integrated Push Handles (Pair)	2	1
14		1509107	Backrest Tube, Straight (415mm), Integrated Push Handles (Pair)	2	1
15		1509108	Backrest Tube, Straight (430mm), Integrated Push Handles (Pair)	2	1
16		1588217	Kit, Screw, Backrest Tube	1	1

Backrest - Rear Set Adjustable Height Push Handles




SCM_231A

Backrest - Rear Set Adjustable Height Push Handles

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1517876	Push Handles, Adjustable Height, Rear Set (Pair)	2	1
2		1438938	Push Handle Tube Complete, Rearset	1	1
3		1438939	Mounting Hardware - Right Complete, Rearset Push Handles	1	1
4		1438940	Mounting Hardware - Left Complete, Rearset Push Handles	1	1
5		1529799	Adjusting Lever (Pair)	2	1
6		1581293	Sleeves, Mounting Hardware, Rearset Push Handles (Pair)	2	1
7		1546913	Telescopic Back Tube, Standard - Short (Pair)	2	1
8		1546914	Telescopic Back Tube, Standard - Medium (Pair)	2	1
9		1546946	Telescopic Back Tube, Standard - Long (Pair)	2	1
10		1588217	Kit, Screw, Backrest Tube	1	1

Backrest - Matrix Back




SCM_230A

Backrest - Matrix Back

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1438203	Kit, Adjustable Height, Integrated Push Handles, Standard, Short	1	1
2		1549985	Adjusting Lever, Adjustable Height, Integrated Push Handles (Pair)	1	1
3		1507734	Handle, Black (Pair)	2	1
4		1561018	Tube, Plug Ø24mm (Pair)	2	1
5		1588219	Kit, Intermediate Tube, Long	1	1
6		1540065	Kit, Screw, Telescopic Tube	1	1

Backrest - Standard Cover




Backrest - Standard Cover

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1-3		See Chart	Backrest - Standard Cover	1	1
4		1432083	Kit, Hardware, Backrest Band	1	1

Back Width	Item #1 RH30-34,5	Item #2 RH36-42	Item #3 RH43,5-46,5
SB36	1552072	1552084	1552096
SB38	1552073	1552085	1552097
SB40	1552074	1552086	1552098
SB42	1552075	1552087	1552099
SB44	1552076	1552088	1552100
SB46	1552077	1552089	1552101
SB48	1552078	1552090	1552102


Backrest - Premium Cover



SCM_430A

Size	Short	Middle	Long	Extra Short
SB28-30	SP1605999	SP1606005	SP1606011	SP1606368
SB32-34	SP1606000	SP1606006	SP1606012	SP1606369
SB36-38	SP1606001	SP1606007	SP1606013	SP1606370
SB40-42	SP1606002	SP1606008	SP1606014	SP1606407
SB44-46	SP1606003	SP1606009	SP1606015	SP1606408
SB48-50	SP1606004	SP1606010	SP1606016	SP1606409

Backrest - Adjustable Bands




SCM_246A

Backrest - Adjustable Bands

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1513167	Adjustment Strap w/ One Strap (Pair)	2	1
2		1513168	Adjustment Strap w/ Two Straps (Pair)	2	1
3		1513169	Adjustment Strap w/ Four Straps (Pair)	2	1
4		1513170	Backrest Band, Two Part, S, Width (5cm)	1	1
5		1513171	Backrest Band, Two Part, Width (10cm)	1	1
6		1513172	Adjustable End Part (5 cm)	1	1
7		1513173	Adjustable End Part (10cm)	1	1
8		1432083	Kit, Hardware, Backrest Band	1	1


Backrest - Adjustable Tension Cover



Backrest Cover, Adjustable Tension (Aeroback)

Aeroback	Item #1 RH30-33	Item #2 RH34,5-40,5	Item #3 RH42-46,5
SB36-38	1512157	1512162	1512167
SB40-42	1512158	1512163	1512168
SB44-46	1512159	1512164	1512169
SB48	1512160	1512165	1512170

Backrest - Stabilization Bar




SCM_248A

Backrest - Stabilization Bar

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1517980	Stabilizing Bar, Complete, Backrest Tube - SB36	1	1
2		1517981	Stabilizing Bar, Complete, Backrest Tube - SB38	1	1
3		1517982	Stabilizing Bar, Complete, Backrest Tube - SB40	1	1
4		1517983	Stabilizing Bar, Complete, Backrest Tube - SB42	1	1
5		1517984	Stabilizing Bar, Complete, Backrest Tube - SB44	1	1
6		1517985	Stabilizing Bar, Complete, Backrest Tube - SB46	1	1
7		1517986	Stabilizing Bar, Complete, Backrest Tube - SB48	1	1
8		1518000	Stabilizing Bar - SB36	1	1
9		1518001	Stabilizing Bar - SB38	1	1
10		1518002	Stabilizing Bar - SB40	1	1
11		1518003	Stabilizing Bar - SB42	1	1
12		1518004	Stabilizing Bar - SB44	1	1
13		1518005	Stabilizing Bar - SB46	1	1
14		1518006	Stabilizing Bar - SB48	1	1
15		1512093	Kit, Hardware, Stabilizing Bar, Left	1	1
16		1512094	Kit, Hardware, Stabilizing Bar, Right	1	1

Footrest - Standard Frame




Footrest - Standard Frame

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1-3		See Chart	Footrest - Standard	1	1
4		1512177	Footrest Narrow	1	1
5		1512178	Footrest Medium	1	1
6		1512179	Footrest Wide	1	1
7		1588357	Kit, Telescoic Footplate, Aluminum Frame	1	1
8		1588358	Kit, Telescoic Footplate, Titanium Frame	1	1
9		1588359	Kit, Telescoic Footplate, Carbon Frame	1	1
10		1588355	Lateral Hardware, Conical, Short	1	1
11		1588356	Lateral Hardware, Conical, Long	2	1
12		1588360	Footrest Sleeve, Aluminum Frame	1	1
13		1588361	Footrest Sleeve, Titanium Frame	1	1
14		1588362	Footrest Sleeve, Carbon Frame	1	1
15		1579076	Tube, Short, Footplate Joint (Pair)	2	1
16		1579077	Tube, Medium, Footplate Joint (Pair)	2	1
17		1579078	Tube, Long, Footplate Joint (Pair)	2	1
18		1432123	Kit, Hardware, Footplate	1	1
19		1513175	Spring, Folding Aid, Champion 2003	1	1

Footrest - Standard

Option	Item #1 Aluminum	Item #2 Titanium	Item #3 Carbon
SB36	1588328	1588334	1588349
SB38	1588329	1588335	1588350
SB40	1588330	1588336	1588351
SB42	1588331	1588346	1588352
SB44	1588332	1588347	1588353
SB46	1588333	1588348	1588354

Footrest - Standard Frame - High Mounted



SCM_251A


Footrest - Standard Frame - High Mounted

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1-2		See Chart	Footrest - Standard - High Mounted	1	1
3	A	1512177	Footrest Narrow	1	1
4	B	1512178	Footrest Medium	1	1
5	C	1512179	Footrest Wide	1	1
6	D	1588355	Lateral Hardware, Conical, Short	1	1
7	E	1588356	Lateral Hardware, Conical, Long	1	1
8		1588391	Kit, Telescopic Tube, Footrest, High Mounted (Pair)	2	1
9	F	1544613	Kit, Clamping, Footrest, High Mounted (Pair)	2	1
10	G	1588392	Kit, Clamping, Footrest, High Mounted, Frame, S Carbon (Pair)	2	1
11	A	1579076	Tube, Short, Footplate Joint (Pair)	2	1
12	B	1579077	Tube, Medium, Footplate Joint (Pair)	2	1
13	C	1579078	Tube, Long, Footplate Joint (Pair)	2	1
14		1432123	Kit, Hardware, Footplate	1	1
15		1513175	Spring, Folding Aid, Champion 2003	1	1
NOTE: A - Used with SB38-SB40. B - Used with SB42-SB44. C - Used with SB46. D - Used with SB38, SB42, and SB46. E - Used with SB40 and SB44. F - Aluminum and Titanium. G - Carbon.					

Footrest - Standard - High Mounted

Option	Item #1 Aluminum or Titanium	Item #2 Carbon
SB38	1588369	1588386
SB40	1588370	1588387
SB42	1588371	1588388
SB44	1588384	1588389
SB46	1588385	1588390

Footrest - Adducted Frame




Footrest - Adducted Frame

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		See Chart	Footrest - Adducted Frame	1	1
2	A	1512177	Footrest Narrow	1	1
3	B	1512178	Footrest Medium	1	1
4		1588355	Lateral Hardware, Conical, Short	1	1
5		1588356	Lateral Hardware, Conical, Long	1	1
6		1588391	Kit, Telescopic Tube, Footrest, High Mounted (Pair)	2	1
7		1544613	Kit, Clamping, Footrest, High Mounted (Pair)	2	1
8	A	1579076	Tube, Short, Footplate Joint (Pair)	2	1
9	B	1579077	Tube, Medium, Footplate Joint (Pair)	2	1
10		1432123	Kit, Hardware, Footplate	1	1
11		1513175	Spring, Folding Aid, Champion 2003	1	1
NOTE: A - Used with SB42-SB44. B - Used with SB46-SB48.					

Footrest - Adducted Frame

Option	Item #1
SB42	1588393
SB44	1588394
SB46	1588395
SB48	1588396

Footrest - Adducted Frame - High Mounted




Footrest - Adducted Frame - High Mounted

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		See Chart	Footrest - Adducted Frame - High Mounted	1	1
2	A	1512177	Footrest Narrow	1	1
3	B	1512178	Footrest Medium	1	1
4	C	1588355	Lateral Hardware, Conical, Short	1	1
5	D	1588356	Lateral Hardware, Conical, Long	1	1
6		1588391	Kit, Telescopic Tube, Footrest, High Mounted (Pair)	2	1
7		1544613	Kit, Clamping, Footrest, High Mounted (Pair)	2	1
8	A	1579076	Tube, Short, Footplate Joint (Pair)	2	1
9	B	1579077	Tube, Medium, Footplate Joint (Pair)	2	1
10		1432123	Kit, Hardware, Footplate	1	1
11		1513175	Spring, Folding Aid, Champion 2003	1	1
NOTE: A - Used with SB42-SB44. B - Used with SB46-SB48. C - Used with SB42-SB46. D - Used with SB44-SB48.					

Footrest - Adducted Frame - High Mounted

Option	Item #1
SB42	1588393
SB44	1588394
SB46	1588395
SB48	1588396

Footrest - Calf Straps




SCM_253A

Footrest - Calf Straps

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A,B	1561529	Calf Strap, Complete (940mm)	1	1
2	C,D	1588397	Calf Strap, Assembled Supporter (660mm)	1	1
3	E,F	1588398	Calf Strap, Assembled Supporter (780mm)	1	1
4	G	1588399	Calf Strap, Assembled Supporter (940mm)	1	1
NOTE: A - Used with SB36 to SB46 Standard Frames. B - Used with SB40 to SB48 Adducted Frames. C - Used with SB36 Standard Frames. D - Used with SB40 Adducted Frames. E - Used with SB38 to SB44 Standard Frames. F - Used with SB44 to SB48 Adducted Frames. G - Used with SB48 Standard Frames.					


Siderests - Clothes Guard



Siderests - Clothes Guard

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1588250	Kit, Clothes Guard, Carbon, Small w/ Hardware (Pair)	2	1
2		1588251	Kit, Clothes Guard, Carbon, Medium w/ Hardware (Pair)	2	1
3		1588252	Kit, Clothes Guard, Carbon, Large w/ Hardware (Pair)	2	1
4		1588253	Kit, Clothes Guard, Carbon, Extra Large w/ Hardware (Pair)	2	1
5		1588254	Clothes Guard, Carbon, Small	1	2
6		1588255	Clothes Guard, Carbon, Medium	1	2
7		1588256	Clothes Guard, Carbon, Large	1	2
8		1588257	Clothes Guard, Carbon, Extra Large	1	2
9		1588258	Kit, Hardware, Rear, Clothes Guard	1	1
10		1588259	Kit, Hardware, Front, Clothes Guard	1	1

Siderests - Fixed Mudguards




SCM_255A

Siderests - Fixed Mudguards

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1588323	Mudguard, Carbon, Small, Complete (Pair)	2	1
2		1588324	Mudguard, Carbon, Large, Complete (Pair)	2	1
3		1529782	Mudguard, Carbon, Small, Short - Right	1	2
4		1529784	Mudguard, Carbon, Large, Long - Right	1	2
5		1529781	Mudguard, Carbon, Small, Short - Left	1	2
6		1529783	Mudguard, Carbon, Large, Long - Left	1	2
7		1588325	Kit, Hardware, Rear, Mudguard	1	1
8		1588326	Kit, Hardware, Front, Mudguard	1	1


Siderests - Removable Mudguards



Siderests - Removable Mudguards

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1588229	Mudguard, Carbon, Extra Small, Long - Left	1	1
2		1588231	Mudguard, Carbon, Small, Long - Left	1	1
3		1588233	Mudguard, Carbon, Medium, Long - Left	1	1
4		1588235	Mudguard, Carbon, Large, Long - Left	1	1
5		1588237	Mudguard, Carbon, Extra Large, :Long - Left	1	1
6		1588230	Mudguard, Carbon, Extra Small, Long - Right	1	1
7		1588232	Mudguard, Carbon, Small, Long - Right	1	1
8		1588234	Mudguard, Carbon, Medium, Long - Right	1	1
9		1588236	Mudguard, Carbon, Large, Long - Right	1	1
10		1588238	Mudguard, Carbon, Extra Long, Long - Right	1	1
11		1588239	Mudguard, Carbon, Extra Small - Left	1	1
12		1588241	Mudguard, Carbon, Small - Left	1	1
13		1588243	Mudguard, Carbon, Medium - Left	1	1
14		1588245	Mudguard, Carbon, Large - Left	1	1
15		1588247	Mudguard, Carbon, Extra Large - Left	1	1
16		1588240	Mudguard, Carbon, Extra Small - Right	1	1
17		1588242	Mudguard, Carbon, Small - Right	1	1
18		1588244	Mudguard, Carbon, Medium - Right	1	1
19		1588246	Mudguard, Carbon, Large - Right	1	1
20		1588248	Mudguard, Carbon, Extra Large - Right	1	1
21		1547532	Mudguard, Carbon Holder, Long (Pair)	2	1
22		1547537	Kit, Siderest, Height Clamp (Pair)	2	1
23		1588249	Kit, Hardware Mudguard, Vertical Brace (Pair)	2	1

Siderests - Adjustable Arm




SCM_257A

Siderests - Adjustable Arm

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1517677	Siderest, Adjustable Height, Short - Right	1	1
2		1517678	Siderest, Adjustable Height, Short - Left	1	1
3		1517682	Siderest, Adjustable Height, Short (Pair)	2	1
4		1439497	Arm Pad w/ Screws, Short (270mm) (Pair)	2	1
5		1547537	Kit, Siderest, Height Clamp (Pair)	2	1
6		1588327	Hardware, Siderest, Frame (Pair)	2	1

Siderests - Tubular Armrest




SCM_258A

Siderests - Tubular Armrest

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		SP1589262	Tubular Arm w/ Standard Pad, Complete w/ Hardware (Pair)	2	1
2		SP1589263	Tubular Arm w/ Standard Pad, Complete w/ Hardware, Right	1	1
3		SP1589264	Tubular Arm w/ Standard Pad, Complete w/ Hardware, Left	1	1
4		SP1589265	Tubular Arm w/ Standard Pad (Pair)	2	1
5		SP1589266	Tubular Arm w/ Standard Pad, Right	1	1
6		SP1589267	Tubular Arm w/ Standard Pad, Left	1	1
7		SP1589268	Tubular Arm Holder w/ Hardware (Pair)	2	1
8		SP1589269	Tubular Arm Holder w/ Hardware, Right	1	1
9		SP1589270	Tubular Arm Holder w/ Hardware, Left	1	1
10		SP1589271	Tubular Arm (Pair)	2	1
11		SP1589272	Tubular Arm, Right	1	1
12		SP1589273	Tubular Arm, Left	1	1
13		SP1589274	Kit, Screws, Armrest Holder (Pair)	2	1
14		SP1589275	Kit, Screws, Tubular Arm (Pair)	2	1
15		SP1589276	Standard Pad	1	2

Front Wheels - Forks - Aluminum and Titanium Frame



SCM_259A

Front Wheels - Forks - Aluminum and Titanium Frame

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1562560-0004	Supporter, Front Wheel (Pair), Black	2	1
2		1562561-0004	Supporter, Front Wheel - Left, Black	1	1
3		1562562-0004	Supporter, Front Wheel - Right, Black	1	1
4		SP1590122	Kit, Fork, Front Pin Bent (Pair)	2	1
5		SP1590123	Kit, Fork, Front Pin, Straight (Pair)	2	1
6		See Chart	Front Wheels - Fork - Aluminum and Titanium Frame	2	1
7		See Chart	Front Wheels - Fork - Aluminum and Titanium Frame	2	1
8		See Chart	Front Wheels - Fork - Aluminum and Titanium Frame	2	1
9		1510219	Kit, Bearing	1	1
		See Chart	Configuration Table, Front Fork Mounted, 75 Degree Frame - Short	1	1
		See Chart	Configuration Table, Front Fork Mounted, 75 Degree Frame - Long	1	1
		See Chart	Configuration Table, Front Fork Mounted, 85 Degree Frame - Short	1	1
		See Chart	Configuration Table, Front Fork Mounted, 85 Degree Frame - Long	1	1

Front Wheels - Fork - Aluminum and Titanium Frame

Color	Item #6 (Small)	Item #7 (Medium)	Item #8 (Large)
Chrome	1562565-0001	1562566-0001	1562567-0001
Black Anodized	1562565-0004	1562566-0004	1562567-0004
Champagne	1562565-0182	1562566-0182	1562567-0182
Techno-Anthrazit	1562565-0400	1562566-0400	1562567-0400
B - Blue	1562565-0404	1562566-0404	1562567-0404
Comet Red	1562565-0440	1562566-0440	1562567-0440
Arctic Blue	1562565-0441	1562566-0441	1562567-0441

Front Wheels - Forks - Aluminum and Titanium Frame

Configuration Table, Front Fork Mounted, 75 Degree Frame - Short

Seat Height Front	Seat Height Rear	ST340-380 Caster 3"	ST340-380 Caster 4"	ST340-380 Caster 5"	ST400-480 Caster 3"	ST400-480 Caster 4"	ST400-480 Caster 5"
450	390	L	M	M	M	S	S
450	400	M	M	S	M	S	S
450	410	M	M	S	M	S	S
450	420	M	S	S	M	S	S
450	430	M	S	S	M	S	S
450	440	M	S	S	M	S	S
450	450	S	S	S	S	S	S
460	390	L	L	M	L	M	M
460	400	L	L	M	M	M	S
460	410	L	M	M	M	M	S
460	420	L	M	M	M	M	S
460	430	M	M	S	M	M	S
460	440	M	M	S	M	S	S
460	450	M	M	S	M	S	S
460	460	M	S	S	M	S	S
470	400	L	L	L	L	M	M
470	410	L	L	M	L	M	M
470	420	L	L	M	L	M	M
470	430	L	M	M	L	M	M
470	440	L	M	M	M	M	S
470	450	L	M	M	M	M	S
470	460	M	M	S	M	M	S
470	470	M	M	S	M	M	S
480	410	L	L	L	L	L	M
480	420	L	L	L	L	L	M
480	430	L	L	L	L	L	M
480	440	L	L	M	L	M	M
480	450	L	L	M	L	M	M
480	460	L	L	M	L	M	M
480	470	L	M	M	L	M	M
480	480	L	M	M	L	M	M
490	420	-	L	L	L	L	L
490	430	-	L	L	L	L	L
490	440	-	L	L	L	L	M
490	450	-	L	L	L	L	M
490	460	-	L	L	L	L	M
490	470	L	L	M	L	L	M
490	480	L	L	M	L	L	M
490	490	L	L	M	L	L	M
500	430	-	-	L	-	L	L
500	440	-	-	L	-	L	L

Front Wheels - Forks - Aluminum and Titanium Frame

Configuration Table, Front Fork Mounted, 75 Degree Frame - Short

Seat Height Front	Seat Height Rear	ST340-380 Caster 3"	ST340-380 Caster 4"	ST340-380 Caster 5"	ST400-480 Caster 3"	ST400-480 Caster 4"	ST400-480 Caster 5"
500	450	-	-	L	-	L	L
500	460	-	-	L	-	L	L
500	470	-	L	L	-	L	L
500	480	-	L	L	-	L	L
500	490	-	L	L	-	L	L
510	440	-	-	-	-	-	L
510	450	-	-	-	-	-	L
510	460	-	-	-	-	-	L
510	470	-	-	-	-	-	L
510	480	-	-	L	-	-	L
510	490	-	-	L	-	-	L
520	450	-	-	-	-	-	-
520	460	-	-	-	-	-	-
520	470	-	-	-	-	-	-
520	480	-	-	-	-	-	-
520	490	-	-	-	-	-	-
530	460	-	-	-	-	-	-
530	470	-	-	-	-	-	-
530	480	-	-	-	-	-	-
530	490	-	-	-	-	-	-
540	470	-	-	-	-	-	-
540	480	-	-	-	-	-	-
540	490	-	-	-	-	-	-

Front Wheels - Forks - Aluminum and Titanium Frame

Configuration Table, Front Fork Mounted, 75 Degree Frame - Long

Seat Height Front	Seat Height Rear	ST340-380 Caster 3"	ST340-380 Caster 4"	ST340-380 Caster 5"	ST400-480 Caster 3"	ST400-480 Caster 4"	ST400-480 Caster 5"
450	390	S	S	S	S	S	S
450	400	S	S	S	S	S	S
450	410	S	S	S	S	S	S
450	420	S	S	S	S	S	S
450	430	S	S	S	S	S	S
450	440	S	S	S	S	S	S
450	450	S	S	S	S	S	S
460	390	M	S	S	S	S	S
460	400	S	S	S	S	S	S
460	410	S	S	S	S	S	S
460	420	S	S	S	S	S	S
460	430	S	S	S	S	S	S
460	440	S	S	S	S	S	S
460	450	S	S	S	S	S	S
460	460	S	S	S	S	S	S
470	400	M	S	S	S	S	S
470	410	M	S	S	S	S	S
470	420	S	S	S	S	S	S
470	430	S	S	S	S	S	S
470	440	S	S	S	S	S	S
470	450	S	S	S	S	S	S
470	460	S	S	S	S	S	S
470	470	S	S	S	S	S	S
480	410	M	M	S	M	S	S
480	420	M	M	S	S	S	S
480	430	M	S	S	S	S	S
480	440	M	S	S	S	S	S
480	450	S	S	S	S	S	S
480	460	S	S	S	S	S	S
480	470	S	S	S	S	S	S
480	480	S	S	S	S	S	S
490	420	L	M	M	M	S	S
490	430	M	M	S	M	S	S
490	440	M	M	S	M	S	S
490	450	M	S	S	M	S	S
490	460	M	S	S	S	S	S
490	470	M	S	S	S	S	S
490	480	S	S	S	S	S	S
490	490	S	S	S	S	S	S
500	430	L	L	M	M	M	S
500	440	L	M	M	M	M	S

Front Wheels - Forks - Aluminum and Titanium Frame

Configuration Table, Front Fork Mounted, 75 Degree Frame - Long

Seat Height Front	Seat Height Rear	ST340-380 Caster 3"	ST340-380 Caster 4"	ST340-380 Caster 5"	ST400-480 Caster 3"	ST400-480 Caster 4"	ST400-480 Caster 5"
500	450	L	M	M	M	M	S
500	460	M	M	S	M	S	S
500	470	M	M	S	M	S	S
500	480	M	S	S	M	S	S
500	490	M	S	S	M	S	S
510	440	L	L	L	L	M	M
510	450	L	L	M	L	M	M
510	460	L	L	M	M	M	S
510	470	L	M	M	M	M	S
510	480	L	M	M	M	M	S
510	490	M	M	S	M	M	S
520	450	L	L	L	L	L	M
520	460	L	L	L	L	M	M
520	470	L	L	M	L	M	M
520	480	L	L	M	L	M	M
520	490	L	M	M	L	M	M
530	460	-	L	L	L	L	L
530	470	-	L	L	L	L	M
530	480	-	L	L	L	L	M
530	490	-	L	L	L	L	M
540	470	-	-	L	-	L	L
540	480	-	-	L	-	L	L
540	490	-	-	L	-	L	L

Front Wheels - Forks - Aluminum and Titanium Frame

Configuration Table, Front Fork Mounted, 85 Degree Frame - Short

Seat Height Front	Seat Height Rear	ST340-380 Caster 3"	ST340-380 Caster 4"	ST340-380 Caster 5"	ST400-480 Caster 3"	ST400-480 Caster 4"	ST400-480 Caster 5"
450	390	M	M	S	M	S	S
450	400	M	M	S	M	S	S
450	410	M	S	S	S	S	S
450	420	S	S	S	S	S	S
460	390	L	M	M	M	M	S
460	400	L	M	M	M	M	S
460	410	M	M	S	M	S	S
460	420	M	M	S	M	S	S
460	430	M	M	S	M	S	S
470	390	L	L	M	L	M	M
470	400	L	L	M	L	M	M
470	410	L	M	M	M	M	S
470	420	L	M	M	M	M	S
470	430	L	M	M	M	M	S
470	440	L	M	M	M	M	S
480	390	L	L	L	L	L	M
480	400	L	L	L	L	L	M
480	410	L	L	L	L	M	M
480	420	L	L	M	L	M	M
480	430	L	L	M	L	M	M
480	440	L	M	M	L	M	M
480	450	L	M	M	L	M	M
490	390	-	L	L	L	L	L
490	400	-	L	L	L	L	L
490	410	-	L	L	L	L	M
490	420	-	L	L	L	L	M
490	430	-	L	L	L	L	M
490	440	-	L	L	L	L	M
490	450	L	L	M	L	L	M
490	460	L	L	M	L	L	M
500	400	-	-	L	-	L	L
500	410	-	-	L	-	L	L
500	420	-	-	L	-	L	L
500	430	-	-	L	-	L	L
500	440	-	-	L	-	L	L
500	450	-	L	L	-	L	L
500	460	-	L	L	L	L	L
500	470	-	L	L	L	L	L
510	410	-	-	-	-	-	L
510	420	-	-	-	-	-	L
510	430	-	-	-	-	-	L

Front Wheels - Forks - Aluminum and Titanium Frame

Configuration Table, Front Fork Mounted, 85 Degree Frame - Short

Seat Height Front	Seat Height Rear	ST340-380 Caster 3"	ST340-380 Caster 4"	ST340-380 Caster 5"	ST400-480 Caster 3"	ST400-480 Caster 4"	ST400-480 Caster 5"
510	440	-	-	-	-	-	L
510	450	-	-	-	-	-	L
510	460	-	-	L	-	L	L
510	470	-	-	L	-	L	L
510	480	-	-	L	-	L	L
520	420	-	-	-	-	-	-
520	430	-	-	-	-	-	-
520	440	-	-	-	-	-	-
520	450	-	-	-	-	-	-
520	460	-	-	-	-	-	L
520	470	-	-	-	-	-	L
520	480	-	-	-	-	-	L
520	490	-	-	-	-	-	L
530	430	-	-	-	-	-	-
530	440	-	-	-	-	-	-
530	450	-	-	-	-	-	-
530	460	-	-	-	-	-	-
530	470	-	-	-	-	-	-
530	480	-	-	-	-	-	-
530	490	-	-	-	-	-	-
540	440	-	-	-	-	-	-
540	450	-	-	-	-	-	-
540	460	-	-	-	-	-	-
540	470	-	-	-	-	-	-
540	480	-	-	-	-	-	-
540	490	-	-	-	-	-	-

Front Wheels - Forks - Aluminum and Titanium Frame

Configuration Table, Front Fork Mounted, 85 Degree Frame - Long


Seat Height Front	Seat Height Rear	ST340-380 Caster 3"	ST340-380 Caster 4"	ST340-380 Caster 5"	ST400-480 Caster 3"	ST400-480 Caster 4"	ST400-480 Caster 5"
450	390	S	S	S	S	S	S
450	400	S	S	S	S	S	S
450	410	S	S	S	S	S	S
450	420	S	S	S	S	S	S
460	390	S	S	S	S	S	S
460	400	S	S	S	S	S	S
460	410	S	S	S	S	S	S
460	420	S	S	S	S	S	S
460	430	S	S	S	S	S	S
470	390	M	S	S	S	S	S
470	400	S	S	S	S	S	S
470	410	S	S	S	S	S	S
470	420	S	S	S	S	S	S
470	430	S	S	S	S	S	S
470	440	S	S	S	S	S	S
480	390	M	M	S	S	S	S
480	400	M	M	S	S	S	S
480	410	M	S	S	S	S	S
480	420	M	S	S	S	S	S
480	430	S	S	S	S	S	S
480	440	S	S	S	S	S	S
480	450	S	S	S	S	S	S
490	390	L	M	M	M	S	S
490	400	L	M	M	M	S	S
490	410	L	M	M	M	S	S
490	420	M	M	S	M	S	S
490	430	M	M	S	S	S	S
490	440	M	S	S	S	S	S
490	450	M	S	S	S	S	S
490	460	M	S	S	S	S	S
500	400	L	L	M	M	M	S
500	410	L	L	M	M	M	S
500	420	L	M	M	M	M	S
500	430	L	M	M	M	S	S
500	440	M	M	S	M	S	S
500	450	M	M	S	M	S	S
500	460	M	M	S	M	S	S
500	470	M	S	S	M	S	S
510	410	L	L	L	L	M	M
510	420	L	L	L	L	M	M
510	430	L	L	M	M	M	S

Front Wheels - Forks - Aluminum and Titanium Frame

Configuration Table, Front Fork Mounted, 85 Degree Frame - Long

Seat Height Front	Seat Height Rear	ST340-380 Caster 3"	ST340-380 Caster 4"	ST340-380 Caster 5"	ST400-480 Caster 3"	ST400-480 Caster 4"	ST400-480 Caster 5"
510	440	L	L	M	M	M	S
510	450	L	M	M	M	M	S
510	460	L	M	M	M	M	S
510	470	M	M	S	M	M	S
510	480	M	M	S	M	M	S
520	420	L	L	L	L	L	M
520	430	L	L	L	L	L	M
520	440	L	L	L	L	M	M
520	450	L	L	M	L	M	M
520	460	L	L	M	L	M	M
520	470	L	L	M	L	M	M
520	480	L	M	M	L	M	M
520	490	L	M	M	M	M	S
530	430	-	L	L	L	L	L
530	440	-	L	L	L	L	M
530	450	-	L	L	L	L	M
530	460	-	L	L	L	L	M
530	470	-	L	L	L	L	M
530	480	L	L	M	L	M	M
530	490	L	L	M	L	M	M
540	440	-	-	L	-	L	L
540	450	-	-	L	-	L	L
540	460	-	-	L	-	L	L
540	470	-	-	L	-	L	L
540	480	-	-	L	L	L	M
540	490	-	L	L	L	L	M

Front Wheels - Fork - Carbon Frame



SCM_260A

Front Wheels - Fork - Carbon Frame

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1-3		See Chart	Items #1-3. Front Wheels - Fork Support - Carbon Frame	1	2
4		See Chart	Item #4. Front Wheels - Fork Support - Carbon Frame	1	1
5		See Chart	Item #5. Front Wheels - Fork Support - Carbon Frame	1	1
6		1451718	Kit, Bearing	1	1
7		See Chart	Item #7. Bearing Block Cap	1	1
8		1512124	Kit, Excenter, Carbon Frame	1	1

Items #1-3. Front Wheels - Fork - Carbon Frame

Color	Item #1 5 in. Fork	Item #2 6 in. Fork	Item #3 8 in. Fork
Chrome	1432004-0001	1432005-0004	1432006-0001
Black Anodized	1432004-0004	1432005-0001	1432006-0004
Champagne	1432004-0182	1432005-0182	1432006-0182
Techno-Anthraxit	1432004-0400	1432005-0400	1432006-0400
B - Blue	1432004-0404	1432005-0404	1432006-0404
Comet Red	1432004-0440	1432005-0440	1432006-0440
Arctic Blue	1432004-0441	1432005-0441	1432006-0441

Item #4. Bearing Block A w/ Bearing

Option	Item #4 Bearing Block "B"
Chrome	1511607-0001
Black Anodized	1511607-0004
Champagne	1511607-0182
Techno-Anthraxit	1511607-0400
B - Blue	1511607-0404
Comet Red	1511607-0440
Arctic Blue	1511607-0441

Front Wheels - Fork - Carbon Frame




Item #5. Bearing Block B w/ Bearing

Option	Item #5 Bearing Block "B"
Chrome	1511608-0001
Black Anodized	1511608-0004
Champagne	1511608-0182
Techno-Anthrazit	1511608-0400
B - Blue	1511608-0404
Comet Red	1511608-0440
Arctic Blue	1511608-0441

Item #7. Bearing Block Cap

Option	Item #7. Bearing Block Cap
Chrome	1487961-0001
Black Anodized	1487961-0004
Champagne	1487961-0182
Techno-Anthrazit	1487961-0400
B - Blue	1487961-0404
Comet Red	1487961-0440
Arctic Blue	1487961-0441

Front Wheels - Casters



SCM_261A

Front Wheels - Casters

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1513180	3" Sports Wheel, Transparent	1	2
2		1513181	4" Skater Wheel, Transparent	1	2
3		SP1638032-NLA	NLA - 4" Caster, Black, Aluminum Rim, Frog Legs	1	2
3		SP1638083	4" Caster, Silver, Aluminum Rim, Frog Legs	1	2
4		SP1638084	5" Caster, Black, Aluminum Rim, Frog Legs	1	2
4		SP1638085	5" Caster, Silver, Aluminum Rim, Frog Legs	1	2
5		See Chart	Item #5. 4 in. Starec (Kuschall) Casters	1	2
6		See Chart	Item #6. 5 in. Starec (Kuschall) Casters	1	2
7		1562626-NLA	NLA - 4" Everyday Wheel, Three Spokes, Black	1	2
7		SP1638072	4" Comfort Wheel, Black	1	2
8		SP1638071	5" Everyday Wheel, 3 Spokes, Black (Pair)	2	1
8		SP1638073	5" Comfort Wheel, Black	1	2
9		1588188-NLA	NLA - Axle Set Complete 3" Sport / 4" and 5" Starec / 4" Skater Wheel	1	2
10		1570108	Axle Set w/out Spacer	1	2
11		1570104	Ball Bearing 608-2RS 8/22*7 (Pair)	2	1
12		1570105	Ball Bearing 6000-2RS (Pair)	2	1

Item #5. 4 in. Starec (Kuschall) Casters


Color	Part Number
Chrome	1513183-0001
Champagne	1513183-0182
Techno-Anthrazit	1513183-0400
B - Blue	1513183-0404
Black Anodized	1513183-0004
Comet Red	1513183-0440
Arctic Blue	1513183-0441

Item #6. 5 in. Starec (Kuschall) Casters

Color	Part Number
Chrome	1513187-0001
Champagne	1513187-0182
Techno-Anthracite	1513187-0400
B - Blue	1513187-0404
Black Anodized	1513187-0004
Comet Red	1513187-0440
Arctic Blue	1513187-0441

Front Wheels - Casters

Rear Wheels - 22 Inch Wheels



SCM_262A



Rear Wheels - 22 Inch Wheels

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		See Chart	Item #1. Rear Wheel w/o Tire	1	2
2		SP1601741	Rear Wheel (22"), Starec Hub X, No Handrim, Rightrun Black/Grey	1	2
3		904.1139	Rear Wheel (22"), Starec Hub X, No Handrim, Treaded	1	2
4		SP1601860	Rear Wheel (22"), Starec Hub X, Handrim Aluminum Black, Rightrun Black/Grey	1	2
5		1562656	Rear Wheel (22"), Starec Hub X, Handrim Aluminum Black, Treaded	1	2
6		1571073	Spoke Reflector (Pair)	2	1
7		1579082	Kit, Spokes, Complete (228mm)	10	1

Item #1. Rear Wheel w/o Tire

Color	Part Number
Colorless	1559083-0407
Chrome	1559083-0001
Champagne	1559083-0182
Techno-Anthrazit	1559083-0400
B - Blue	1559083-0404
Comet Red	1559083-0440
Arctic Blue	1559083-0441

Rear Wheels - 24 Inch Wheels



SCM_263A

Rear Wheels - 24 Inch Wheels

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		See Chart	Item #1. Rear Wheel, 24 in., w/o Handrim, w/o Tire	1	2
2		SP1601745	Rear Wheel (24"), Starec Hub X, No Handrim, Rightrun Black/Grey Tire	1	2
3		904.1141	Rear Wheel (24"), Starec Hub X, No Handrim, Treaded	1	2
4		See Chart	Item #4. Rear Wheel, 24 in., w/o Handrim	1	2
5		SP1601864	Rear Wheel (24"), Starec Hub X, Aluminum Handrim, Black, Rightrun Black/Grey	1	2
6		1562660	Rear Wheel (24"), Starec Hub X, Aluminum Handrim, Black, Treaded	1	2
7		See Chart	Item #7. Rear Wheel, 24 in., w/ Handrim	1	2
8		1571073	Spoke Reflector (Pair)	2	1
9		1579086	Kit, Spoke, Complete (253mm)	1	1

Item #1. Rear Wheel, 24 in., w/o Handrim, w/o Tire

Hub Color	Part Number
Colorless	1567104-0407-NLA
Chrome	1567104-0001
Champagne	1567104-0182-NLA
Techno-Antrazit	1567104-0400
B - Blue	1567104-0404
Comet Red	1567104-0440
Arctic Blue	1567104-0441

Item #4. Rear Wheel, 24 in., w/o Handrim


Hub Color	Part Number
Colorless	1554399-0407
Chrome	1554399-0001
Champagne	1554399-0182
Techno-Antrazit	1554399-0400-NLA
B - Blue	1554399-0404
Comet Red	1554399-0440
Arctic Blue	1554399-0441

Rear Wheels - 24 Inch Wheels

Item #7. Rear Wheel, 24 in., w/ Handrim

Hub Color	Part Number
Colorless	1554400-0407
Chrome	1554400-0001
Champagne	1554400-0182
Techno-Antrazit	1554400-0400
B - Blue	1554400-0404
Comet Red	1554400-0440
Arctic Blue	1554400-0441

Rear Wheels - 26 Inch Wheels



SCM_264A

Rear Wheels - 26 Inch Wheels

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		See Chart	Item #1. Rear Wheel, 26 Inches, w/o Handrim, w/o Tire	1	2
2		SP1601749	Rear Wheel 26", Starec Hub X, w/out Handrim, Rightrun Black/Grey	1	2
3		See Chart	Item #3. Rear Wheel, 26 Inches, w/o Handrim, w/ Tire	1	2
4		SP1601866	Rear Wheel 26", Starec Hub X, Handrim Aluminum Black, Rightrun Black/Grey	1	2
5		See Chart	Item #5. Rear Wheel, 26 Inches, w/ Handrim, and Tire	1	2
6		1571073	Spoke Reflector (Pair)	2	1
7		1579090	Kit, Spokes, Complete (278mm)	10	1

Item #1. Rear Wheel, 26 Inches, w/o Handrim, w/o Tire

Colors	Part Number
Colorless	1567099-0407
Chrome	1567099-0001
Champagne	1567099-0182
Techno-Antracite	1567099-0400
B - Blue	1567099-0404
Comet Red	1567099-0440
Arctic Blue	1567099-0441

Item #3. Rear Wheel, 26 Inches, w/o Handrim, w/ Tire


Colors	Part Number
Colorless	1554403-0407
Chrome	1554403-0001
Champagne	1554403-0182
Techno-Antracite	1554403-0400
B - Blue	1554403-0404
Comet Red	1554403-0440
Arctic Blue	1554403-0441

Rear Wheels - 26 Inch Wheels

Item #5. Rear Wheel, 26 Inches, w/ Handrim, and Tire

Colors	Part Number
Clear	1554417-0407
Chrome	1554417-0001
Champagne	1554417-0182
Techno-Antracite	1554417-0400
B - Blue	1554417-0404
Comet Red	1554417-0440
Arctic Blue	1554417-0441

Rear Wheels - 24 and 25 Inch Wheels - High Performance



SCM_265A

Rear Wheels - 24 and 25 Inch Wheels - High Performance

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		See Chart	Item #1. Rear Wheel, 24 Inches w/ Handrim and Rightrun Tire	1	2
2		See Chart	Item #2. Rear Wheel, 24 Inches, w/ Handrim and Evo Tire	1	2
3		See Chart	Item #3. Rear Wheel, 25 Inches, w/ Handrim and Evo Tire	1	2
4		See Chart	Item #4. Rear Wheel, 24 Inches, w/o Handrim and w/o Tire	1	2
5		See Chart	Item #5. Rear Wheel, 25 Inches, w/o Handrim and w/o Tire	1	2
6		1571073	Spoke Reflector (Pair)	2	1
7		1578867	Kit, Spokes, High Performance Wheel (24")	10	1
8		1578868	Kit, Spokes, High Performance Wheel (25")	10	1

Item #1. Rear Wheel, 24 Inches, w/ Handrim and Rightrun Tire

Colors	Part Number
Black	SP1601752-0004
Chrome	SP1601752-0001
Champagne	SP1601752-0182
Techno-Anthracite	SP1601752-0400
B - Blue	SP1601752-0404
Comet Red	SP1601752-0440
Arctic Blue	SP1601752-0441

Item #2. Rear Wheel, 24 Inches, w/ Handrim and Evo Tire

Colors	Part Number
Clear	1562670-0004
Chrome	1562670-0001
Champagne	1562670-0182
Techno-Antracite	1562670-0400
B - Blue	1562670-0404
Comet Red	1562670-0440
Arctic Blue	1562670-0441

Rear Wheels - 24 and 25 Inch Wheels - High Performance

Item #3. Rear Wheel, 25 Inches, w/ Handrim and Evo Tire

Colors	Part Number
Champagne	1562669-0182
Techno-Anthracite	1562669-0400
B - Blue	1562669-0404
Comet Red	1562669-0440
Arctic Blue	1562669-0441

Item #4. Rear Wheel, 24 Inches, w/o Handrim, w/o Tire


Colors	Part Number
Black	1517723-0004
Chrome	1517723-0001
Champagne	1517723-0182
Techno-Antracite	1517723-0400
B - Blue	1517723-0404
Comet Red	1517723-0440
Arctic Blue	1517723-0441

Rear Wheels - 24 and 25 Inch Wheels - High Performance


Item #5. Rear Wheel, 25 Inches, w/o Handrim, w/o Tire

Colors	Part Number
Black	1562937-0004
Chrome	1562937-0001
Champagne	1562937-0182
Techno-Antracite	1562937-0400
B - Blue	1562937-0404
Comet Red	1562937-0440
Arctic Blue	1562937-0441

Rear Wheels - Spinergy LXX - 24" and 25" Wheels




HANDRIM	SCHWALBE		HANDRIM	SCHWALBE	
		24"			25"
		24"			25"



3

24"

25"



6

4

24" INSIDE

25" INSIDE

5

24" OUTSIDE

25" OUTSIDE

SCM_266B

Rear Wheels - Spinergy LXX - 24" and 25" Wheels

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1579101	Rear Wheel (24"), Spinergy LXX, Aluminum Handrim, Black, Schwalbe One Tire, Complete	1	2
2		1579105	Rear Wheel (24"), Spinergy LXX, Aluminum Handrim, Black, Marathon Evo Tire, Complete	1	2
3		1579102	Rear Wheel (24"), Spinergy LXX, No Handrim, No Tire	1	2
4		1579103	(24") Spoke Inside, Spinergy LXX Wheel (236mm)	1	9
5		1579104	(24") Spoke Outside, Spinergy LXX Wheel (249mm)	1	9
6		1571088	Spoke Reflector (Pair)	2	1

Rear Wheels - Handrims

HANDRIMS		ALUMINUM				BLACK ALUMINUM				TITANIUM		POLISHED STEEL		SUPER GRIP			SURGE	
																		
		22"	24"	25"	26"	22"	24"	25"	26"	24"	24"	26"	24"	25"	26"	24"	26"	
		1	3	5	7	9	10	11	12	13	15	16	17	18	19	20	21	
		2	4	6	8					14			17	18	19	20	21	
			3	5			10	11		13	15			18		20	21	
			4	6						14				18		20	21	
			3	5			10			13	15			18		20	21	
			4	6						14				18		20	21	


SCM 267A

SCM_267A

Rear Wheels - Handrims

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1568713-0005	Handrim, Aluminum, Distance (22" - 3cm)	1	2
2		1568700	Handrim, Aluminum (22" - 4cm)	1	2
3		1568714-0005	Handrim, Aluminum, Distance (24" - 3cm)	1	2
4		1568702-NLA	NLA - Handrim, Aluminum Light, Clear Anodized (24" - 4cm)	1	2
5		1568730-0005	Handrim 25" Aluminum, Distance 3cm	1	2
6		1568708	Handrim, Aluminum Light, High Performance (25" - 4cm)	1	2
7		1568731-0005	Handrim, Aluminum, Distance (26" - 3cm)	1	2
8		1568709	Handrim, Aluminum (26" - 4cm)	1	2
9		1568713-0411	Handrim, Aluminum, Black, Distance (22" - 3cm)	1	2
10		1568714-0411	Handrim, Aluminum, Black, Distance (24" - 3cm)	1	2
11		1568730-0411	Handrim, Aluminum, Black, Distance (25" - 3cm)	1	2
12		1568731-0411	Handrim, Aluminum, Black, Distance (26" - 3cm)	1	2
13		1568716	Handrim, Titanium, Distance (24" - 3cm)	1	2
14		1568705	Handrim 24" Titanium, 4cm	1	2
15		1568715	Handrim, Polished Steel, Distance (24" - 3cm)	1	2
16		1568777	Handrim, Polished Steel, Distance (26" - 3cm)	1	2
17		1568701	Handrim, Super-Grip (22")	1	2
18		1568703	Handrim, Super-Grip (24")	1	2
19		1568710	Handrim, Super-Grip (26")	1	2
20	A	SP60122339	Handrim 24" Carbolife Gekko S	1	2
21		1556917	Handrim, LT Surge (25")	1	2
NOTE: A - 1556916 is no longer available. Use part number SP60122339 as a replacement.					



Rear Wheels - Handrims - Carbolife



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		SP1638081	Handrim, Carbolife, Curve (24")	1	2
2		SP1638080	Handrim, Carbolife, Gekko (24")	1	2



Rear Wheels - Handrim Hardware

ALUMINUM		BLACK ALUMINUM		TITANIUM		STEEL		SUPERGRIP		SURGE		GEKKO		CURVE	
															
SPOKE GUARD	NO SPOKE GUARD	SPOKE GUARD	NO SPOKE GUARD	SPOKE GUARD	NO SPOKE GUARD	SPOKE GUARD	NO SPOKE GUARD	SPOKE GUARD	NO SPOKE GUARD	SPOKE GUARD	NO SPOKE GUARD	SPOKE GUARD	NO SPOKE GUARD	SPOKE GUARD	NO SPOKE GUARD
1	1	1	1	1	1	1	1	1	1	3	3	1	1	1	1
4	4	4	4	4	4	4	4	4	4	4	6	4	4	4	4
											7				
															KUSCHALL WHEEL
															HIGH PERF. WHEEL
															SPINERGY WHEEL

3cm

ALUMINUM		BLACK ALUMINUM		TITANIUM		STEEL		SUPERGRIP		SURGE	
											
SPOKE GUARD	NO SPOKE GUARD	SPOKE GUARD	NO SPOKE GUARD	SPOKE GUARD	NO SPOKE GUARD	SPOKE GUARD	NO SPOKE GUARD	SPOKE GUARD	NO SPOKE GUARD	SPOKE GUARD	NO SPOKE GUARD
3	2			3	2	3	2	3	2		
6	5			6	5	6	5	6	5		
	5				5		5		5		
											KUSCHALL WHEEL
											HIGH PERF. WHEEL
											SPINERGY WHEEL

4cm

SCM_282A

Rear Wheels - Handrim Hardware

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1495958	Kit, Hardware, Handrim (3 cm)	1	2
2		1495959	Kit, Hardware, Handrim (4 cm)	1	2
3		1495961	Kit, Hardware, Handrim w/ Spokes Guard (4cm)	1	2
4		1496019	Kit, Hardware, Handrim, High Performance (3 cm)	1	2
5		1496020	Kit, Hardware, Handrim, High Performance (4cm)	1	2
6		1497111	Kit, Hardware, Handrim, High Performance w/ Spokes Guard (4cm)	1	2
7	A	1496015	Kit, Hardware, Handrim, Mountain Bike (3cm)	1	2
NOTE: A - GEKKO and CURVE Handrims use this Handrim Hardware.					

Rear Wheels - Tires


		RIGHTRUN TIRE				TREADED TIRE		MARATHON EVO			SCHWALBE TIRE		SOLID TIRE	
														
		22"	24"	25"	26"	22"	24"	24"	25"	26"	24"	25"	22"	24"
TIRES	STAREC WHEEL	1	4		12	2	5	6		13	7		3	8
	HIGH PERFORMANCE		4	9			5	6	10		7	11		
	SPINERGY LXX		4				5	6			7			

SCM_268A

Rear Wheels - Tires

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		SP1590117	Tire, Rightrun, Black/Grey (22")	1	2
2		1570836	Tire, Treaded (22")	1	2
3		1570837	Tire, Solid, Light, Grey (22")	1	2
4		SP1590118	Tire, Rightrun, Black/Grey (24")	1	2
5		1570829	Tire, Treaded (24")	1	2
6		1570838	Tire, Marathon Plus Evo (24")	1	2
7		1578802	Tire, Schwalbe One (24")	1	2
8		1570839	Tire, Solid, Light, Grey (24")	1	2
9		SP1590412	Tire, Rightrun Black/Grey (25")	1	2
10		1570841	Tire, Marathon Plus Evo (25")	1	2
11		1578803	Tire, Schwalbe One (25")	1	2
12		SP1590119	Tire, Rightrun, Black/Grey (26")	1	2
13		1570842	Tire, Marathon Plus Evo (26")	1	2

Rear Wheels - Inner Tube and Rim Protection




SCM_283A

Rear Wheels - Inner Tube and Rim Protection

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1570848	Inner Tube (22" x 1")	1	2
2		1570849	Inner Tube (24" x 1")	1	2
3		1570851	Inner Tube, High-Performance Wheel (26")	1	2
4		1570850	Inner Tube (26" x 1")	1	2
5		1570844	Rim Protection (22")	1	2
6		1570845	Rim Protection (24" Large)	1	2
7		1495962	Rim Protection, High Performance / Spinergy Wheel (24")	1	2
8		1570846	Rim Protection (26")	1	2
9		1571072	Rear Wheel Insert Flat Free (20, 22, 24, 26" Tire)	1	2

Rear Wheels - Spoke Guards




SCM_269A

Rear Wheels - Spoke Guards

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1564453	Spoke Guard, Transparent (22")	1	2
2	A,B	1564457	Spoke Guard, Transparent (24")	1	2
3	A	1564492	Spoke Guard, Transparent (26")	1	2
4	A,B	1564456	Spoke Guard Kuschall, Black (24")	1	2
NOTE: A - Used with Starec Wheels. B - Used with Starec and High Performance Wheels.					

Rear Wheels - Quick Release Axles and Bearings



SCM_270A


Rear Wheels - Quick Release Axles and Bearings

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		See Chart	Item #1. Quick Release Axle - Long (Pair)	2	1
2		1507627	Quick Release Axle, Long w/ Cap, Spinergy Wheel (Pair)	2	1
3		1493755	Tetra Clip w/ Quick Release Axle (Pair), Long	2	1
4		1571069	Bearing R8-2RS 1/2"/1.1/8"x5/16" (Pair)	2	1

Item #1. Quick Release Axle - Long (Pair)

Color	Part Number
Arctic Blue	1509137-0441
B - Blue	1509137-0404
Black	1509137-0004
Chrome	1509137-0001
Champagne	1509137-0182
Colorless	1509137-0407
Comet Red	1509137-0440
Techno-Anthracite	1509137-0400


Rear Wheels - Wheel Locks



Rear Wheels - Wheel Locks

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1512127	Brake, Straight Lever, Left	1	1
2		1512128	Brake, Straight Lever, Right	1	1
3		1517658	Brake, Straight Lever, Long Shoe, Left	1	1
4		1517657	Brake, Straight Lever, Long Shoe, Right	1	1
5		1512125	Brake, Bent Lever, Left	1	1
6		1512126	Brake, Bent Lever, Right	1	1
7		1517656	Brake, Lever Bent, Long Shoe, Left	1	1
8		1517655	Brake, Lever Bent, Long Shoe, Right	1	1
9		SP1512131	Performance Brake Left	1	1
10		SP1512132	Performance Brake Right	1	1
11		SP1512133	Active Brake Left	1	1
12		SP1512134	Active Brake Right	1	1
13		1512129	Pull Brake Left	1	1
14		1512130	Pull Brake Right	1	1
15		1517660	Pull Brake, Long Shoe Left	1	1
16		1517659	Pull Brake, Long Shoe Right	1	1
17		1513176	Clamping piece , Brake, set	1	1
18		1513177	Brake Handle Standard (Pair)	2	1
19		1513178	Brake Handle, Performance Brake (Pair)	2	1
20		1513179	Brake Lever Extension (Pair)	2	1
21		1544424	Sleeves, Brake Clamping, Titanium Frame (Pair)	2	1
22		1579328	Kit, Brake Shoe Standard (40mm)	1	1
23		1579329	Kit, Brake Shoe Long (50mm), 4 Pcs.	1	1



Seat Accessories - Seat Bag



SCM_272A

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		SP1589933	Seat Bag SB36	1	1
2		SP1589934	Seat Bag SB38	1	1
3		SP1589935	Seat Bag SB40	1	1
4		SP1589936	Seat Bag SB42	1	1
5		SP1589937	Seat Bag SB44	1	1
6		SP1589938	Seat Bag SB46	1	1
7		SP1589942	Seat Bag SB48	1	1
8		SP1589943	Kit, Fixing Clip	1	1


Seat Accessories - Fixing Clips



SCM_273A


Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1589115	Kit, Fixing Clips, Short	1	1
2		1589116	Kit, Fixing Clips, Long	1	1
3		1589078	Kit, Hardware Tape	1	1

Options and Accessories - Anti-Tipper




Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1588417	Anti-Tipper, Complete - Right	1	1
2		1588418	Anti-Tipper, Complete - Left	1	1
3		1588419	Anti-Tipper Holder, Complete - Right	1	1
4		1588776	Anti-Tipper Holder, Complete - Left	1	1
5		1588777	Anti-Tipper Tube, Complete	1	1
6		1588798	Kit, Hardware, Anti-Tipper	1	1

Options and Accessories - Transit Wheels




Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1588799	Transit Wheels Complete (Pair)	2	1
2		1588800	Transit Wheels Holder (Pair)	2	1
3		1588777	Anti-Tipper Tube, Complete	1	2
4		1588801	Transit Wheels Parts (Pair)	2	1
5		1588798	Kit, Hardware, Anti-Tipper	1	1

Options and Accessories - Tipper Aids




Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1588804	Tipper Aid, Complete	1	1
2		1588803	Tipper Aid and Cane Holder Combination, Complete	1	1
3		1588802	Kit, Adaptation Holder	1	1

Options and Accessories - Cane Holder



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1588805	Cane Holder, Complete	1	1
2		1588803	Tipper Aid and Cane Holder Combination, Complete	1	1
3		1588802	Kit, Adaptation Holder	1	1

Options and Accessories - Luggage Holder



SCM_286A


Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1547626	Luggage Holder Complete (Pair)	2	1

Options and Accessories - Transport Bag



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1512059	Transport Bag, Standard	1	1
2		1512060	Transport Bag, Large	1	1

Options and Accessories - Chair Accessories



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1518285	Tool Kit	1	1
2		1517924	Air Pump	1	1
3		1518286	Air Man	1	1
4		1518287	Wheel Bag 24" - 26"	1	1
5		1547625	Bottle Holder	1	1
6		1547622	Mobile Phone Holder	1	1
7		1589117	Kit, Kuschall Watch Mounted, Bearing Block	1	1
8		SP1601739	Kit, Kuschall Watch Mounted, Supporter, Front Wheel	1	1

Invacare Corporation

www.invacare.com



Yes, you can.

One Invacare Way
Elyria, Oh USA
44036-2125
800-333-6900

Technical Support
800-832-4707
Customer Service
844-695-8724

Fax Orders
866-577-3507

© 2022 Invacare Corporation. All rights reserved. Republication, duplication or modification in whole or in part is prohibited without prior written permission from Invacare. Trademarks are identified by ™ and ®. All trademarks are owned by or licensed to Invacare Corporation or its subsidiaries unless otherwise noted. (Revised 8/18/22)

Form No. 16-249