

ICC - Allto Bed
Model BH2380/76



Yes, you can.

Parts Catalog Usage Guide

The information contained in this document is subject to change without notice. This catalog contains service parts for Form No. 16-253. If you have any questions regarding this catalog, please contact one of our Customer Service Representatives.

For Customer Service Support Call:

1-800-333-6900

For Faxed In Orders Call:

1-800-678-4682 or

1-440-366-0709

For Technical Service Call:

1-800-832-4707

Literature Can Be Found at:

www.invacare.com or

1-800-828-6272

SYMBOL DEFINITION

Throughout this catalog you will notice letters to draw your attention to specific footnotes which are located at the bottom of the page or below specified charts.

Examples:

A,B,C etc. = Special Comments

NOTE: = General information pertaining to that specific page

NLA = No Longer Available

HARDWARE MEASUREMENTS

The measurements for the hardware will be listed as follows:

Screws, Bolts Diameter X Thread Count/Pitch X Length

Washers Inside Diameter X Outside Diameter X Thickness

Spacers Inside Diameter X Outside Diameter X Length

For your convenience, we have supplied an area for you to list your Invacare Representatives.

Name	Phone Numbers


© 2020 Invacare Corporation. All rights reserved. Reproduction, duplication or modification in whole or in part is prohibited without prior written permission from Invacare. Trademarks are identified by ™ and ®. All trademarks are owned by or licensed to Invacare Corporation or its subsidiaries unless otherwise noted.

TABLE of CONTENTS


Description	Page Number
Electronics & Pendants	3

Electronics & Pendants


NOTE: This Bed is obsolete and only parts listed are available to order.
NOTE: All Frame Parts are obsolete and No Longer Available. For Deck Replacements, see Echo 2380/76 Deck Replacement Guide.




Shaft Approx. 3 Ft. - Red Plug



Shaft Approx. 1 Ft. - Blue Plug



Shaft Approx. 1-1/2 Ft. - White Plug



Comes w/ Groundwire and Clips



DB9 Connector



DB9 Connector

SICC_110A

Electronics & Pendants

NOTE: This Bed is obsolete and only parts listed are available to order.

NOTE: All Frame Parts are obsolete and No Longer Available. For Deck Replacements, see Echo 2380/76 Deck Replacement Guide.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	MTRM1703	Head Actuator (Red Plug)	1	1
2	B	MTRM1702	Foot Actuator (Blue Plug)	1	1
3	C	MTRM1747	NLA - Hilo Actuator (White Plug)	1	1
4	D	JNCTNBOXKIT1 -NLA	NLA - Junction Control Box (Comes w/ Ground Wire & Clips)	1	1
5	E	1190233	Kit, Pendant, Full/Partial Rails, Linak # PNDNT-6FNC-1-2 (With DB9 Connector)	1	1
6	E	PNDNTKIT3-NL A	NLA - Kit, Pendant, Assist Rails (Old Square), Linak # PNDNT-6FNC-1-2 (With DB9 Connector) (May 01 - July 02)	1	1
<p>NOTE: A - Will work on old Allto Beds. VO7508AW76U B - Will work on old Allto Beds. VO7508AX76U C - Reference 6183M-1905220-01 D - J Box Prior to 11/29/01. NAX5105-3A29MPZ E - These Pendants are not designed for the Allto, but these will work. However the Pendant cannot be fastened to the Control Box using the screws.</p>					

Invacare Corporation

www.invacare.com



Yes, you can.

One Invacare Way
Elyria, Oh USA
44036-2125
800-333-6900

Technical Support
800-832-4707
Customer Service
800-333-6900

Literature Orders
800-333-6900

© 2020 Invacare Corporation. All rights reserved. Republication, duplication or modification in whole or in part is prohibited without prior written permission from Invacare. Trademarks are identified by TM and ©. All trademarks are owned by or licensed to Invacare Corporation or its subsidiaries unless otherwise noted. (Revised 03/23/19)

Form No. 16-253