

Battery Boxes and Trays

NOTE: Light Duty, 250-S Models use UI batteries with 22NF battery boxes.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	G	1097007	Battery Box Assembly, Rear, 22NF	1	2
2	A,B,G	1097003	Battery Box Top, Rear, 22NF	1	2
3	H	1057944	Battery Box Bottom with Strap, 22NF	1	2
3		1057954-NLA	NLA - Strap, Battery	1	2
3		1041433-NLA	NLA - Pin, Handle	1	2
4	I	1071532-NLA	NLA - Battery Tray Assembly with Hardware	1	1
5		1071490-NLA	NLA - Tray with Strap and Buckle	1	2
6		1025752-NLA	NLA - Screw, Hex Head (1/4-20 x 1-1/2")	1	4
7		1071566-NLA	NLA - Pin, Tray	1	4
8		15680X012	Collar, Anti-Rattle (7/8")	1	4
9		7930333	Locknut (1/4-20)	25	1
10		1071531-NLA	NLA - Support, Tray	1	1
11		1021406	Screw, Hex Head (1/4-20 x 1-1/4")	1	4
12		1074127-NLA	NLA - Decal, Tray	1	1
13		1071543-NLA	NLA - Bracket, Tray, Black	1	2
14		1071505	Screw, Socket w/ Patch (M6 x 1 x 20mm)	1	3
15	J	1071504-NLA	NLA - Screw, Socket w/ Patch (M6 x 1 x 35mm)	1	3
16		1079193-NLA	NLA - Tray, Battery with Strap, Weight Shifting Tilt	1	1
	A,E	50000X166-NLA	NLA - Fuse - Buss 2-1/4" L (60 AMP)	1	2
	E,F	1089385	Fuse, MIDI (Non-G/B & G/B) (80 AMP)	1	2
	C,E,K	1060318-NLA	NLA - Hardware, Battery Terminal, UI Batteries (Pair)	2	1
	D,E	1058618-NLA	NLA - Hardware, Battery Terminal, 22NF Batteries (Pair)	2	1

NOTE: A - As of 11/5/99 the battery harness fuse changed from Bussman 60 AMP to a MIDI 80 AMP fuse. This battery box top is compatible for all chairs manufactured before 11/5/99 and must be ordered whenever a replacement Bussman 60 AMP fuse is requested.
 B - It is recommended that new terminal hardware be ordered when a new battery box top is ordered.
 C - Hardware includes small battery terminal caps (gray and orange), 1/4-20 x 7/8" hex flange screws, 1/4-20 hex flange nuts, cable ties, locking tab washers and instruction sheet for two batteries.
 D - Hardware includes battery clamps, battery clamp caps (red and black), small battery terminal caps (red and black), 1/4-20 x 7/8" hex flange screws, 1/4-20 hex flange nuts, cable ties and locking tab washers for two batteries.
 E - Not Shown.
 F - After 11/4/99 - Dual Battery Box or After 11/9/00 - Single Battery Box.
 G - 22NF and UI 250-S Models.
 H - #1057944 includes: (2) each #1057954 Battery Strap and (2) each #1041433 Handle Pin.
 I - Includes items 5-13.
 J - (3) each required for Non-Weight Shifting Tilt and (2) each required for Weight Shifting Tilt.
 K - 250-S Models.