


# Shroud Covers - Before 3/1/00



**NOTE:** When converting an integrated sling seat to accommodate shrouds, order (1070573) charger connector/trunk stop bracket and attaching bolt (1045267).

# Shroud Covers - Before 3/1/00

## Item #'s 1-5 Shroud Covers Before 3/1/00

Item Number	Description	Bright Plum Metallic	Bright Blue Metallic	Bright Silver Metallic	Blue Green Chameleon	Candy Apple Red
1	Side Panel (Captain) - Right	1074942-NLA See Note(s) B, D	1074943-NLA See Note(s) B, D	1074945-NLA See Note(s) B, D	1074944-NLA See Note(s) B, D	1085094-NLA See Note(s) B, D
1	Side Panel (Captain) - Left	1074946-NLA See Note(s) B, D	1074947-NLA See Note(s) B, A, D	1074949-NLA See Note(s) B, D	1074948-NLA See Note(s) B, D	1085095-NLA See Note(s) B, D
1	Side Panel (Sling Seat) - Right	1074975-NLA See Note(s) B, D	1074976-NLA See Note(s) B, A, D	1074978-NLA See Note(s) B, D	1074977-NLA See Note(s) D, B	1085092-NLA See Note(s) A, B, D
1	Side Panel (Sling Seat) - Left	1074979-NLA See Note(s) B, E	1074980-NLA See Note(s) B, E	1074982-NLA See Note(s) E, B	1074981-NLA See Note(s) E, B	1085093-NLA See Note(s) B, D, E
2	Battery Cover - 16" Wide	1074970-NLA See Note(s) C	1074971-NLA See Note(s) C	1074973-NLA See Note(s) C	1074972-NLA See Note(s) A, C	1085088-NLA See Note(s) C
2	Battery Cover - 18" Wide	1074966-NLA See Note(s) A, C	1074967-NLA See Note(s) C	1074969-NLA See Note(s) C	1074968-NLA See Note(s) C	1085087-NLA See Note(s) C
2	Battery Cover - 20" Wide	1074962-NLA See Note(s) C	1074963-NLA See Note(s) A	1074965-NLA See Note(s) C	1074964-NLA See Note(s) C	1085086-NLA See Note(s) C
3	Trunk Lid - 16" Wide	1074958-NLA See Note(s) C, A	1074959-NLA See Note(s) C	1074961-NLA See Note(s) C, A	1074960-NLA See Note(s) C	1085091-NLA See Note(s) C
3	Trunk Lid - 18" Wide	1074954-NLA See Note(s) C	1074955-NLA See Note(s) C	1074957-NLA See Note(s) C	1074956-NLA See Note(s) C	1085090-NLA See Note(s) C
3	Trunk Lid - 20" Wide	1074950-NLA See Note(s) C, A	1074951-NLA See Note(s) C	1074953-NLA See Note(s) C	1074952-NLA See Note(s) C	1085089-NLA See Note(s) C
4	Rear Storage Compartment, - 16" Wide	1070562-NLA See Note(s) C	1070562-NLA See Note(s) C	1070562-NLA See Note(s) C	1070562-NLA See Note(s) C	1070562-NLA See Note(s) C
4	Rear Storage Compartment, - 18" Wide	1070561-NLA See Note(s) C	1070561-NLA See Note(s) C	1070561-NLA See Note(s) C	1070561-NLA See Note(s) C	1070561-NLA See Note(s) C
4	Rear Storage Compartment, - 20" W See Note(s) C	1070560-NLA	1070560-NLA See Note(s) C	1070560-NLA See Note(s) C	1070560-NLA See Note(s) C	1070560-NLA See Note(s) C
5	Storage Shelf, Rear	1071655-NLA	1071655-NLA	1071655-NLA	1071655-NLA	1071655-NLA

NOTE: A - Shrouds made before 3/1/00 are NLA. See page titled, "Shroud Covers - After 2/29/00" and select a color conversion package.  
 B - Decals not included. See page titled, "Packaging, Literature and Decals."  
 C - Captain seat shrouds on 18" wide only.  
 D - Side Panel (Captain Seat) includes rubber bushing.  
 E - Side Panel (Sling Seat) includes rubber bushing.

## Shroud Covers - Before 3/1/00

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1-5		See Chart	Item #'s 1-5 Shroud Covers - Before 3/1/00	1	1
6		1070565-NLA	NLA - Bracket, Shroud, Black - Right	1	1
6		1070566-NLA	NLA - Bracket, Shroud, Black - Left	1	1
7		1045267	Screw, Hex Head w/ Patch (1/4-20 x 3/4")	1	4
8		1028516	Washer, Steel (3/16 x 7/16 x 1/16")	1	2
9		7930332-NLA	NLA - Locknut, High Profile (#10-32)	25	1
10		1070575-NLA	NLA - Screw, SEMS w/ Washer (#10-32 x 7/8")	1	2
11		1070588-NLA	NLA - Screw, Phillips Truss Head (#10-32 x 1/2")	1	6
12		1028755	Screw, Phillips Pan Head (#10-32 x 5/8")	1	4
13		1070587	Washer, Flat (3/16 x 5/8 x 1/16")	1	4
	A,B	1070564-NLA	NLA - Hardware, Side Panel Attaching	1	1
14		1070577-NLA	NLA - Plate, Trunk Clip	1	2
15		1070583-NLA	NLA - Stud, Wing Head (3/4" Long)	1	2
16		1070581-NLA	NLA - Washer, Nylon Cupped (9/32 x 19/32 x 3/32")	1	4
17		1070580-NLA	NLA - Washer, Nylon Retainer (1/4 x 1/2 x 1/32")	1	4
18		1013497-NLA	NLA - Fastener, Turn, Clip-On (1/4")	1	4
19	C	1070585	Stud, Oval Slotted (1" Long)	1	2
19	D	1070574-NLA	NLA - Stud, Wing Head (1" Long)	1	2
20		1070556-NLA	NLA - Fastener, Dual Lock (1" W x 2" L)	1	2
21		1070555-NLA	NLA - Fastener, Dual Lock (5/8" W x 2" L)	1	2
22		1070594-NLA	NLA - Bushing, Rubber	1	2
23		1070563-NLA	NLA - Plate, Trunk	1	2
24		1070582-NLA	NLA - Rivet	1	4
25		1062096-NLA	NLA - Fastener, Loop (1/2" W x 1" L)	1	4
NOTE: A - Not Shown. B - Includes items 6-13. C - Before 2/11/98. D - After 2/10/98.					