







ISP Polaris EX

ISP3000 & ISP4000



ISP Polaris EX

ISP3000 & ISP4000

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A,G	ISP4000S-NLA	NLA - Humidifier, Heated Polaris EX, Silver	1	1
	C	1132517	Cap, Reservoir	1	1
	C	1132522	Reservoir, Water (Only)	1	1
	C	1117424-NLA	NLA - Elbow, Humidifier	1	1
	C,J	1117423-NLA	NLA - Plate, Bottom, Assembly	1	1
	C,K	1132523-NLA	NLA - Bottom Plate Assembly	1	1
	C	1132521-NLA	NLA - O-Ring	1	1
	B,I	1136074-NLA	NLA - Humidifier Reservoir Assembly	1	1
	C	1120165-NLA	NLA - Coupling, Hose	1	1
	C	11008X001-NLA	NLA - Screw, Phillips Pan Head (#10-32 x 1-3/8")	1	1
	B	1124211-NLA	NLA - Package, Replacement Filter - ISP3000	3	1
	B	1124213-NLA	NLA - Door, Access with Screw	1	1
	B	1023281-NLA	NLA - Screw, Phillips Pan Head Tap (#8-18 x 1/2")	1	1
14	B	1119677-NLA	NLA - Case, Carrying	1	1
	A,D	1085420-NLA	NLA - Tube with Cuff, Smooth Bore, Clear, Packaged (6 FT. L)	1	1
	A,C	1119768-NLA	NLA - Insert, Bottom Humidifier	1	2
	A,C	1119769	Carton, Shipper	1	1
	A,B	1119765	Insert, Bottom	1	1
	A,B	1119766-NLA	NLA - Insert, Top	1	1
	A,B	1119767	Carton, Shipper	1	1
	A	1110559	Owner's Manual, ISP3000	1	1
	A	1118358	Owner's Manual, ISP4000	1	1
	A,B,L	1141773-NLA	NLA - Kit, Filter Ultra Fine	1	1
			Filter, Ultra Fine CPAP	3	1
			Instructions, Ultra Fine Filter	1	1

NOTE: A - Not Shown.
B - For use on ISP3000.
C - For use on ISP4000.
D - Includes Part #1075695 Smooth Bore Tube with Cuff enclosed in Polybag.
G - Includes Silver Humidifier Base, Silver Humidifier Shroud, Hot Plate Assembly, Hose Coupling Kit, Shipper Carton and Owner's Manual.
I - Includes Reservoir, Cap, Bottom Plate, O-Ring and Elbow.
J - Plate used with old style Water Reservoir, Before 8/17/05.
K - Used with new Water Reservoir, after 8/16/05.
L - Can be used in addition to CPAP Filter, for additional filtration.