


Cantilever Arms - AT99N

Non-Locking Bolt On

Available on Spree GT chairs only



Cantilever Arms - AT99N

Non-Locking Bolt On

Available on Spree GT chairs only

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1118594	Cantilever Arm Assembly - Right	1	1
2			Spacer, Coved (1/4 x 11/16 x 1/2")	1	4
3			NLA - Cam, Cantilever Arm	1	1
4			Locknut (1/4-28)	1	2
5			Washer (1/4 x 1/2 x 1/16")	1	1
6			Washer, Nylon (1/4 x 1/2 x 1/32")	1	1
7			NLA - Tube, Cantilever Arm, Black	1	1
8			NLA - Screw, Hex Head (1/4-28 x 2-3/4")	1	1
9			NLA - Screw, Hex Head (1/4-28 x 2-1/2")	1	1
10			Grip, Foam Hand (6-7/8" L)	1	1
11			Plug Button, Black (1")	1	2
1	A	1118595	Cantilever Arm Assembly - Left	1	1
2			Spacer, Coved (1/4 x 11/16 x 1/2")	1	4
3			NLA - Cam, Cantilever Arm	1	1
4			Locknut (1/4-28)	1	2
5			Washer (1/4 x 1/2 x 1/16")	1	1
6			Washer, Nylon (1/4 x 1/2 x 1/32")	1	1
7			NLA - Tube, Cantilever Arm, Black	1	1
8			NLA - Screw, Hex Head (1/4-28 x 2-3/4")	1	1
9			NLA - Screw, Hex Head (1/4-28 x 2-1/2")	1	1
10			Grip, Foam Hand (6-7/8" L)	1	1
11			Plug Button, Black (1")	1	2

NOTE: A - Includes items 2-11