

Power Articulating Legrest Kits

ELRPW (After 7/26/04)

NOTE: 2G Tilt Only Systems originally ordered without a T.R.C.M.2G controller must include a T.R.C.M.2G controller and 4-way switch with hardware. Refer to the pages titled "Electronics" for MKIV T.R.C.M. 2G and MKIV T.A.C. of this book for a T.R.C.M. upgrade kit. Therefore, add the power leg kit dated after 7/26/04. Power Articulating Legs use sector style rigging hangers only and do not come packaged with the kits.

For 2GT, 2GR, 2GT/R Seating Systems When a T.R.C.M. 2G Is Present

ELRPW After 7/26/04

Quantity	Part Number	Description
1	1126095-NLA See Note(s) A	Kit, ELRPW
1	1126091	NLA - Legrest Assembly, Power Articulating, Sector, Black (13"-20" Range) - Right
1	1126092	NLA - Legrest Assembly, Power Articulating, Sector, Black (13"-20" Range) - Left
1	1101520-NLA	NLA - Kit, Power Leg Harness and Mounting Hardware
1	1098337-NLA	NLA - Instruction Sheet, Power Leg Installation

NOTE: A - Kit includes the 4 part numbers listed below the kit part number