


IVC™ I-Fit™ & CareGuard™ Shower Chairs

Model 90-2, 95-2, 9780 and 9781



IVC™ I-Fit™ & CareGuard™ Shower Chairs

Model 90-2, 95-2, 9780 and 9781

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1116880	Backrest	1	1
2	B	1116872	Leg Assembly, Composite	1	4
3	B	1121530-NLA	NLA - Tip, Rubber, Grey	1	4
4	C	1094945	Button, Snap	1	4
5	C	1063124	Anti-Rattle, Grey (1-1/8")	1	4
6	C	1110289-NLA	NLA - Tip, Rubber, Grey	1	4
NOTE: A - For use with 96-2 (Aluminum) and 9781 (Composite) models only with Backrests. B - For use with Shower Chairs with Composite Legs (9780 and 9781). C - For use with Shower Chairs with Aluminum Legs (90-2 and 95-2).					