

SX5RC Conversion Kit

NOTE: SX5RC Conversion Kit will convert a SX5 Wheelchair into a SX5RC Wheelchair.

Nylon Upholstery

Chair Width	Part Number
14" Wide	1110201-NLA See Note(s) A
16" Wide	1110202-NLA See Note(s) B
18" Wide	1110203-NLA See Note(s) C
20" Wide	1110204-NLA See Note(s) D
22" Wide	1110205-NLA See Note(s) E

NOTE: A - No Longer Available. Order through the Configurator using RCL14
 B - No Longer Available. Order through the Configurator using RCL16
 C - No Longer Available. Order through the Configurator using RCL18
 D - No Longer Available. Order through the Configurator using RCL20
 E - No Longer Available. Order through the Configurator using RCL22

NOTE: SX5RC Conversion Kit includes: Right/Left Back Assemblies, Back Upholstery, Headrest Upholstery, Spreader Bar, Hardware and Instruction Sheet.

NOTE: Rear Wheel Types U2244BK, 2244BK and 4244BK are not compatible with SX5RC Conversion Kits.

NOTE: If the wheelchair has 20-inch rear wheels, clamp-on wheel locks (Right 8881181914; Left 8881181915) and Mounting Hardware (1030987) (1ea.) must be used with the short wheel lock shoe extension (1048656) (2 ea.).

NOTE: If the wheelchair has 22-inch rear wheels, clamp-on wheel locks (Right 8881181914; Left 8881181915) and Mounting Hardware (1030987) (1ea.) must be used (without a wheel lock shoe extension).

NOTE: If the wheelchair has 24-inch rear wheels, bolt-on wheel locks must be used with the large wheel lock shoe extension (1037421) (1 ea.).

Vinyl Upholstery

Chair Width	Part Number
14" Wide	8881110195-NLA See Note(s) B, A
16" Wide	8881110196-NLA See Note(s) A, C
18" Wide	8881110197-NLA See Note(s) A, D
20" Wide	8881110198-NLA See Note(s) E, A
22" Wide	8881110199-NLA See Note(s) A, F

NOTE: A - Specify upholstery color when ordering, see category titled, "Upholstery Colors."
 B - No Longer Available. Order through the Configurator using RCL14
 C - No Longer Available. Order through the Configurator using RCL16
 D - No Longer Available. Order through the Configurator using RCL18
 E - No Longer Available. Order through the Configurator using RCL20
 F - No Longer Available. Order through the Configurator using RCL22

NOTE: SX5RC Conversion Kit includes: Right/Left Back Assemblies, Back Upholstery, Headrest Upholstery, Spreader Bar, Hardware and Instruction Sheet.

NOTE: Rear Wheel Types U2244BK, 2244BK and 4244BK are not compatible with SX5RC Conversion Kits.

NOTE: If the wheelchair has 20-inch rear wheels, clamp-on wheel locks (Right 8881181914; Left 8881181915) and Mounting Hardware (1030987) (1ea.) must be used with the short wheel lock shoe extension (1048656) (2 ea.).

NOTE: If the wheelchair has 22-inch rear wheels, clamp-on wheel locks (Right 8881181914; Left 8881181915) and Mounting Hardware (1030987) (1ea.) must be used (without a wheel lock shoe extension).

NOTE: If the wheelchair has 24-inch rear wheels, bolt-on wheel locks must be used with the large wheel lock shoe extension (1037421) (1 ea.).