



Frames and Hardware
(Before 6/24/99)



Before 6/95



6/1/95 thru 6/23/99

Frames and Hardware (Before 6/24/99)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		15680X013	Anti-Rattle (1")	1	2
2		1037351-NLA	NLA - Guide, Seat	1	1
3		1059251-NLA	NLA - Screw, Socket Head (1/4-28 x 1-1/2")	1	2
4		1038129-NLA	NLA - Socket, Arm, Front	1	1
5		1052957	Screw, Hex Head w/ Patch (1/4-20 x 2")	1	1
6		1026907	Washer (1/4 x 1/2 x 1/16")	1	1
7	B,A	1054571	Kit, Sector Block - Right	1	1
7	B,A	1054572	Kit, Sector Block - Left	1	1
8		1025765	Locknut (1/4-28)	1	2
9	B,C	8881038582-NL A	NLA - Frame w/o Hardware, Color - Right (16" D)	1	1
9	B,C	8881038583-NL A	NLA - Frame w/o Hardware, Color - Left (16" D)	1	1
10	B,D	8881053371-NL A	NLA - Frame w/o Hardware, Color - Right (16" D)	1	1
10	B,D	8881053372-NL A	NLA - Frame w/o Hardware, Color - Left (16" D)	1	1
10	B,D	8881053336-NL A	NLA - Frame w/o Hardware, Color - Right (18" D)	1	1
10	B,D	8881053337-NL A	NLA - Frame w/o Hardware, Color - Left (18" D)	1	1
NOTE: A - Includes items 5 and 6 B - Specify frame finish, see page titled "Frame Finishes and Touch-up Paint" C - Before 6/95 D - 6/1/95 thru 6/23/99					