



Frames & Hardware (After 6/23/99)



Adjustable Back Height Frame



Adjustable Angle Back Frame

Frames & Hardware (After 6/23/99)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1056948	Guide, Seat	1	1
2		1055347-NLA	NLA - Screw, Hex Head (1/4-28 x 1-1/2")	1	2
3		1056405-NLA	NLA - H-Block with Socket Cup	1	1
4		1052957	Screw, Hex Head w/ Patch (1/4-20 x 2")	1	1
5		1026907	Washer (1/4 x 1/2 x 1/16")	1	1
6	A,B	1054571	Kit, Sector Block - Right	1	1
6	A,B	1054572	Kit, Sector Block - Left	1	1
7		1025765	Locknut (1/4-28)	1	2
8		1013869	Cap, End, Tube, Black (1")	1	1
9	B,C	8881086832	Frame w/o Hardware, Color - Right (16" D)	1	1
9	B,C	8881086833	Frame w/o Hardware, Color - Left (16" D)	1	1
9	B,C	8881086834	Frame w/o Hardware, Color - Right (18" D)	1	1
9	B,C	8881086835	Frame w/o Hardware, Color - Left (18" D)	1	1
10	B,D	8881127076	Frame w/o Hardware, Color - Right (16" D)	1	1
10	B,D	8881127077	Frame w/o Hardware, Color - Left (16" D)	1	1
10	B,D	8881127080-NL A	NLA - Frame w/o Hardware, Color - Right (18" D)	1	1
10	B,D	8881127081	NLA - Frame w/o Hardware, Color - Left (18" D)	1	1
NOTE: A - Includes items 4 and 5 B - Specify frame finish, see page titled "Frame Finishes and Touch-up Paint" C - Adjustable Back Height Frame D - Adjustable Angle Back Frame					