


Shrouds and Mounting Hardware



Shrouds and Mounting Hardware

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		See Chart	#1 Side Shrouds with Decals	1	1
2		1104797	Decal, Drive/Push	1	2
3		1121907-NLA	NLA - Decal, Shroud - Pronto M71 Jr.	1	2
4		1120033-NLA	NLA - Decal, Warning	1	1
5		1025451	Screw, Slotted Hex Washer Head (#10-32 x 5/8")	1	8
6		1070587	Washer, Flat (3/16 x 5/8 x 1/16")	1	12
7		1028755	Screw, Phillips Pan Head (#10-32 x 5/8")	1	4
8		1122958-NLA	NLA - Plate, Shroud - Right	1	1
8		1122959-NLA	NLA - Plate, Shroud - Left	1	1
9	A	1122269-NLA	NLA - Kit, Side Shroud Hardware	2	1
10		1109035-NLA	NLA - Clamp, Shroud	1	1
11		1025713-NLA	NLA - Screw, Hex Head (1/4-20 x 1")	1	1
12		1109036-NLA	NLA - Bracket, Clamp Shroud Mounting, Black - Right	1	1
12		1109037-NLA	NLA - Bracket, Clamp Shroud Mounting, Black - Left	1	1
13		1026907	Washer (1/4 x 1/2 x 1/16")	1	1
14		1025713-NLA	NLA - Screw, Hex Head (1/4-20 x 1")	1	1
15	B	1108464-NLA	NLA - Kit, Shroud Clamp Hardware - Right	1	1
15	B	1108465-NLA	NLA - Kit, Shroud Clamp Hardware - Left	1	1
16		1123240-NLA	NLA - Battery Shroud w/ Decal, Black - Rear	1	1
17		1095489	Decal, Invacare 3-D Logo (2-1/2" L)	1	2
18		1028755	Screw, Phillips Pan Head (#10-32 x 5/8")	1	4
19		1070587	Washer, Flat (3/16 x 5/8 x 1/16")	1	4
20		1123238-NLA	NLA - Battery Shroud with Decals, Black - Front	1	1
21		1118381-NLA	NLA - Label, Warning Battery Terminal	1	1
NOTE: A - Includes items 5-8 B - Includes items 10-14					

#1 Side Shrouds with Decals

Side	Magic Black	Deep Blue Pearl	Candy Red Pearl
Right	1122272-NLA	1122274-NLA	1122270-NLA
Left	1122273-NLA	1122275-NLA	See Note See Note(s) A
NOTE: A - No longer available as a service part.			