




# MK5 Controller Cables




**MK5 EX Controller Option Cable**




**Battery Extension Cable**



**Power Option Cable**



**Tilt Sensor Cable**



**Tilt Sensor Splitter Cable**

## MK5 Controller Cables

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1118231-NLA	NLA - Cable, MK5 EX Controller Option	1	1
2		1122199	Cable, Battery Extension, Red Connector	1	1
3	B	1123923-NLA	NLA - Cable, Power Controller Option, MK5 EX	1	1
4	C	1126122-NLA	NLA - Cable, Tilt Sensor	1	1
5	C	1126123-NLA	NLA - Cable, Splitter Tilt Sensor	1	1
NOTE: A - The 5 pin connector attaches to MKIV joystick options. B - Power Option Cable comes standard on chairs with MK5 EX controller. C - Items 4 and 5 are used in conjunction with each other.					