



Adjustable Height Arms



Arm Attaching Hardware

Adjustable Height Arms

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	8881053578	Kit, Arm Pad w/ Screws, Desk Length (9" L x 2" W) (Pair)	2	1
1	A	8881053577	Kit, Arm Pad w/ Screws, Full Length (14" L x 2" W) (Pair)	2	1
2		1027095	Plug Button, Black (7/8")	1	2
3		1111714	Package, Screw Phillips Pan Head (#10-32 x 1-1/8")	2	1
4		1070859	Arm Tube Assembly w/ Plug Buttons, Desk Length, Black - Right/Left	1	1
4		1041522	Arm Tube Assembly w/ Plug Buttons, Full Length, Black - Right/Left	1	1
5	B	1035892	Adjustable Height Arm Base Assembly - Right	1	1
			Base, Armrest, Adjustable Height	1	1
6			Spring, Arm Lock	1	1
7			Pin, Arm Lock, Black	1	1
8			Rivet	1	6
9			NLA - Guard, Clothing	1	1
10			Kit, Lock Handle w/Rivets (Pair)	1	1
5	B	1035893	Adjustable Height Arm Base Assembly - Left	1	1
			Base, Armrest, Adjustable Height	1	1
6			Spring, Arm Lock	1	1
7			Pin, Arm Lock, Black	1	1
8			Rivet	1	6
9			NLA - Guard, Clothing	1	1
10			Kit, Lock Handle w/Rivets (Pair)	1	1
11	C	1083083	Kit, Arm Attaching Hardware (Pair)	1	1
12			Socket, Arm, Front - Right	1	1
12			Socket, Arm, Front - Left	1	1
13			Socket, Arm, Rear - Right/Left	1	2
14			Screw, Hex Head (1/4-28 x 2-1/4")	1	4
15			Washer (1/4 x 1/2 x 1/16")	1	4
16			Nut, Hex (1/4-28)	1	6
12		1083081	Socket, Arm, Front - Right	1	1
12		1083082	Socket, Arm, Front - Left	1	1
13		1048835	Socket, Arm, Rear - Right/Left	1	2
14		1055350	Screw, Hex Head (1/4-28 x 2-1/4")	1	4
15		1026907	Washer (1/4 x 1/2 x 1/16")	1	4
16		1085111	Nut, Hex (1/4-28)	1	6

NOTE: A - Specify upholstery color U67 Black when ordering
 B - Includes items 6-10
 C - Includes items 12-16