


Contracture Platform Footrest



Contracture Platform Footrest

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1127694-NLA	NLA - Kit, Contracture Platform Hardware	1	1
2			Locknut (1/4-28)	1	2
3			NLA - Screw, Hex Head (5/16-24 x 3")	1	1
4			Screw, Socket Button Head (1/4-28 x 1-1/2")	1	2
5			Bolt, Lug (5/16-24 x 1-3/8")	1	2
6			Locknut (5/16-24)	1	3
7			Locknut (1/4-20)	1	2
8			NLA - Screw, Hex Head (1/4-20 x 1-1/2")	1	2
9	B	1127831-NLA	NLA - Bracket, Contracture Support, Black	1	1
10		1052436-NLA	NLA - Support, Contracture Platform, Black, Long	1	1
11		1052441	Tube, Support, Lower, 110 Degree, Black, Long	1	1
12		1127695	Kit, Contracture Footboard Hardware	1	1
13		Screw, Socket Flat Head (1/4-20 x 3/4")	1	4	
14		NLA - Half Clamp, Threaded	1	2	
15		Screw, Socket Button Head (1/4-20 x 1")	1	8	
16		NLA - Half Clamp, Non-Threaded	1	2	
17		Plug Button, Black (3/4")	1	2	
18		1044076	Footboard	1	1
NOTE: A - Includes items 2-8 B - Includes items 13-17					