




Accessories for Model 7750 (Square Tube)




**1011225
UPPER BED BRACKET**



**1011224
LOWER BED BRACKET
with WALL BUMPER**



**1011229
GRAB BAR, CHAIN, and CLAMP**



**1011223
CLAMP**

Accessories for Model 7750 (Square Tube)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	10301X000-NLA	NLA - Hook, "S"	1	2
2		17573X000-NLA	NLA - Key	1	1
3		17572X000	Rod, Clamp	1	1
4		16324X000	Nut, Wing (3/8"-16)	1	2
5		1046727	Cap, Plug, Grey	1	2
6		1011223-NLA	NLA - Clamp Only	1	1
7		17077X000-NLA	NLA - Hook Rod	1	2
8		40186X019-NLA	NLA - Tube, End Cap	1	1
9		17101X000-NLA	NLA - Grab Bar, Trapeze, White	1	1
10		95103X000-NLA	NLA - Chain, 8 Link	1	1
NOTE: A - As of 8/98 Hook Rod is NLA, replace with Clamp Rod - item #3.					