




Battery Harnesses



Battery Cable Extension Harness



Elevate Actuator and Battery Charger Harness

Battery Harnesses

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A,D	1127304-NLA	NLA - Kit, Battery Harnesses with Terminal Hardware	1	1
			Tie Wrap (11-1/2" L)	1	6
			NLA - Instruction Sheet, Battery Installation	1	1
2			NLA - Harness, Battery - Front	1	1
3			Harness, Battery w/ 75 AMP Fuse and Black Connector	1	1
4			NLA - Harness, Battery - Rear	1	1
5			Locknut (1/4-20)	1	4
6			Washer, Locking Tab	1	2
7			NLA - Screw, Hex Flange Whiz Lock(1/4-20 x 5/8")	1	4
8			Harness, Battery w/ 75 AMP Fuse and Red Connector	1	1
9	B	1116381	Cable, Battery Extension, Black Connector	1	1
10	C	1127260	Harness, Actuator Sensor and Battery Charger	1	1

NOTE: A - Includes items 2-8, Cables Ties (not shown), and Instruction Sheet (not shown).
 B - Used on M50 and M51 models
 C - Used on M61 models
 D - Part number 1125093 is no longer available from the Invacare Service Parts Department. Refer to the Invacare Website under "Product Literature Search". Enter the part number in the "Keywords" field and download a free PDF version copy.