


# Rear Maintenance Preventative Kit

1

Brake Release Lever  
With Nut




Seat Locking Pin




Shock Absorbers With Nut




Seat Post




Locking Pin With  
Screw and Nut



Brake Release Pivot  
With Screw and Nut



2



3



## Rear Maintenance Preventative Kit

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1128212-NLA	NLA - Kit, Rear Preventative Maintenance Includes: 1 ea. Locking Pin with Screw & Nut, 1 ea. Seat Post, 2 ea. Seat Post Locking Pins, 2 ea. Brake Release Lever with Nut, 1 ea. Brake Release Pivot with Screw and Nut, and 4 ea. Shock Absorber with Nut	1	1
2		1132322-NLA	NLA - Pin, Seat Post Locking (Pkg of 10)	10	1
3		1132323-NLA	NLA - Package, Brake Release Lever w/ Nut	10	1