



Battery Box Kit



Battery Box Kit

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1128213-NLA	NLA - Kit, Battery Box Includes: 1 ea. Battery Box Top with Wiring, 1 ea. Battery Box Bottom, 1 ea. Battery Box Hardware, and 1 ea. Jumper Wire	1	1
2		1132327-NLA	NLA - Circuit Breakers	10	1
3		1132328-NLA	NLA - Harness, Charger Port with Cover	1	1
	A	12AH-2	Batteries, 12 volt -12 Amp Hour (Pair)	2	1
NOTE: A - Not Shown					