



Swing Arm Caster



Swing Arm Caster

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1128214-NLA	NLA - Kit, Swing Arm Caster Includes: 2 ea. Solid Rear Wheel Casters with Bearings, 2 ea. Rear Caster Forks with Nuts, 4 ea. Head Tube Bearings, 2 ea. Head Tube Spacers, 2 ea. Axle Bolts with Nuts, 4 ea. Rear Caster Spacers, and 2 ea. Head Tube Caps	1	1
2		1132329-NLA	NLA - Kit, Caps, End, Black, Head Tube	10	1