

Tiller Assembly with Harness

(Manufactured After 12/01)


Lynx SX3P (After S/N 00F623034)

Lynx LX3 (After S/N 00F622314)


Panther LX4 (After S/N 00F622474)

Panther MX4 (After S/N 00F623494)


Lynx SX3P




Lynx SX3




Lynx LX3



Lynx LX3^{Plus}
Panther LX4/MX4



Lynx SX3P/SX3/LX3/LX3^{Plus}
Panther LX4/MX4



Tiller Assembly with Harness
(Manufactured After 12/01)
Lynx SX3P (After S/N 00F623034)
Lynx LX3 (After S/N 00F622314)
Panther LX4 (After S/N 00F622474)
Panther MX4 (After S/N 00F623494)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1102248-NLA	NLA - Tiller Assembly w/ Harness - Lynx SX3P, Taupe	1	1
2	B	1113581-NLA	NLA - Tiller Assembly w/ Harness - Lynx SX3, Candy Red	1	1
2	B	1113582-NLA	NLA - Tiller Assembly w/ Harness - Lynx SX3, Deep Blue	1	1
3	B	1113579-NLA	NLA - Tiller Assembly w/ Harness - Lynx LX3, Candy Red	1	1
3	B	1113580-NLA	NLA - Tiller Assembly w/ Harness - Lynx LX3, Deep Blue	1	1
4	B	1113586-NLA	NLA - Delta Tiller Assembly w/ Harness, Candy Red - Lynx LX3 Plus	1	1
4	B	1113587-NLA	NLA - Delta Tiller Assembly with Harness, Deep Blue - Lynx LX3 Plus	1	1
4	B	1113588-NLA	NLA - Delta Tiller Assembly w/ Harness, Magic Black - Lynx LX3 Plus	1	1
5	B	1113583-NLA	NLA - Delta Tiller Assembly w/ Harness, Candy Red - Panther LX4	1	1
5	B	1113584-NLA	NLA - Delta Tiller Assembly w/ Harness, Deep Blue - Panther LX4	1	1
5	B	1113585-NLA	NLA - Delta Tiller Assembly w/ Harness, Magic Black - Panther LX4	1	1
6		1101278-NLA	NLA - Delta Tiller Assembly w/ Harness, Champagne - Panther MX4	1	1
7		1102927-NLA	NLA - Key Switch with Keys	1	1
7		1088704-NLA	NLA - Keys Only (Pair)	2	1
8		1091863-NLA	NLA - Charger Port - Lynx SX3P, Panther LX4 and Zoom-3	1	1
9		1091862-NLA	NLA - Cover, Charger Port - Lynx SX3/LX3/LX3 Plus and Panther MX4	1	1
10		1088706-NLA	NLA - Grips, Hand (Pair)	2	1
11	E	1088705-NLA	NLA - Knob, Tiller Speed	1	1
12		1088707-NLA	NLA - Canopy, LED	1	1
13	F	1101276-NLA	NLA - Potentiometer	1	1
14		1088708-NLA	NLA - Lever, Throttle	1	1
15		1088709-NLA	NLA - Handle, Latch with Spring and Pin	1	1
16		1099708-NLA	NLA - Swivel Bracket Assembly	1	1
17		1101275-NLA	NLA - Handle Bar w/ Hardware - LX3Plus, LX4 and MX4	1	1
18	D	1101271-NLA	NLA - Cover Assembly, Tiller - LX3 Plus (Woodgrain Finish)	1	1
18	D	1102929-NLA	NLA - Cover Assembly, Tiller - LX3 Plus (Snake Skin Finish)	1	1

Tiller Assembly with Harness
(Manufactured After 12/01)
Lynx SX3P (After S/N 00F623034)
Lynx LX3 (After S/N 00F622314)
Panther LX4 (After S/N 00F622474)
Panther MX4 (After S/N 00F623494)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
18	D	1102393	Cover Assembly, Tiller - LX4 (Woodgrain Finish)	1	1
18	D	1102930-NLA	NLA - Cover Assembly, Tiller - LX4 (Snake Skin Finish)	1	1
18	D	1101277-NLA	NLA - Cover Assembly, Tiller - MX4 (Woodgrain Finish)	1	1
	C	P755	Extension, Throttle Lever (Pair)	2	1
	C	1097017-NLA	NLA - Kit, Decal w/ Instruction Sheet, Forward/Reverse - Left Hand	1	1
	C	1088794-NLA	NLA - Tiller Hardware - All Models Includes: 2 each Bolt (M6 x 1.0 x 35) 2 each Nut, Self-Lock (6 x 10 x 8) 1 each Washer (8 x 39 x 2) 1 each Washer (8 x 39 x 1.6) 1 each Bolt (M8 x 1.25 x 50) 1 each Washer (8 x 22 x 2) 1 each Nut, Self-Lock (8 x 1.25 x 13)	1	1
NOTE: B - After 12/01 C - Not Shown D - Old style has charger port on left side with warning label above; new style has black plug button siliconed inside hole. E - This Knob is used on all tiller models - it is not shown on some. F - The Throttle Lever is included with this item. If the tiller for this part is for is not a "Delta" tiller, please contact the technical services department for instructions on replacing this item.					