



Aerosol Compressor - Stratos Portable

Model - IRC1720

Compressor
Shown for
Reference Only



Aerosol Compressor - Stratos Portable

Model - IRC1720

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1133559-NLA	NLA - Adaptor, AC/DC Switch - Portable	1	1
2		1133561-NLA	NLA - Adaptor, Car - Portable	1	1
3		1133560-NLA	NLA - Cord, Power AC - Portable	1	1
4	B	IRC1725-NLA	NLA - Kit, Nebulizer - Portable	1	1
	A	1133551	Filter, Disk Replacement	5	1
	A	1133562-NLA	NLA - Cover, Filter	1	1
	A	1133563-NLA	NLA - Cover, Power Input	1	1
	A	IRC1722	Case, Carrying	1	1
	A	1130202-NLA	NLA - Owner's Manual - Portable	1	1
	A	1130218	Instuction Sheet, Portable Battery	1	1
	A	IRC1726	Mask - Pediatric	1	1
	A,C	IRC1726-2-NLA	NLA - Kit, Masks, Pediatric	1	1
			Mask - Pediatric	1	2
	A	IRC1727	Mask - Adult	1	1
	A,C	IRC1727-2-NLA	NLA - Kit, Masks, Adult	1	1
			Mask - Adult	1	2

NOTE: A - Not Shown.
 B - Kit includes Top and Bottom Nebulizer, Nozzle, Expiratory Valve, Inhalation Valve, Holder, Mouth Piece and Tubing.
 C - Kit includes (2) Masks.