


Contoura Seating Shroud Covers



Contoura Seating Shroud Covers

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1125705	NLA - Kit, Controller and Seat Frame Shroud Covers with Hardware	1	1
2			NLA - Cover, Back, Controller Shroud	1	1
4			Screw, Button Head w/ Patch (1/4-20 x 1/2")	1	4
5			T-Nut (1/4-20)	1	2
6			NLA - Screw, Phillips Pan Head (#6-32 x 1-1/8")	1	2
7			NLA - Nut, Hex Lock (#6-32)	1	2
8			NLA - Bracket, Shroud, Black - Right	1	1
8			NLA - Bracket, Shroud, Black - Left	1	1
9			Fastener, Dual Lock (1/2" W x 7/8" L)	1	8
10			NLA - Cover Assembly, Seat Frame Shroud	1	1
3	B	1125728-NLA	Kit, Controller Shroud Mounting Hardware	1	1
4			Screw, Button Head w/ Patch (1/4-20 x 1/2")	1	4
5			T-Nut (1/4-20)	1	2
6			NLA - Screw, Phillips Pan Head (#6-32 x 1-1/8")	1	2
7			NLA - Nut, Hex Lock (#6-32)	1	2
8			NLA - Bracket, Shroud, Black - Right	1	1
8			NLA - Bracket, Shroud, Black - Left	1	1
9			Fastener, Dual Lock (1/2" W x 7/8" L)	1	4
10	C	1125704-NLA	NLA - Cover Assembly, Seat Frame Shroud	1	1
9		1090197	Fastener, Dual Lock (1/2" W x 7/8" L)	1	4
11		1114308	Screw, Phillips Pan Head Tap (#10-32 x 1/2")	1	2
NOTE: A - Includes items 2, 4-10 B - Includes items 4-9 C - Includes Shroud Cover, items 9 and 11					