


Footrest Support Assembly

93M2 / H93M2

NOTE: Due to design change, for chairs built prior to 10/1/00 it is suggested that Release Handles be replaced in pairs. See category titled, "Front Riggings" for complete Pivot/Slide Tube with Footplate Assembly.



Footrest Support Assembly

93M2 / H93M2

NOTE: Due to design change, for chairs built prior to 10/1/00 it is suggested that Release Handles be replaced in pairs. See category titled, "Front Riggings" for complete Pivot/Slide Tube with Footplate Assembly.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1035539	Handle, Release, Black - Right	1	1
1		1035540	Handle, Release, Black - Left	1	1
2		47565X057	Bushing (7/32 x 3/8 x 11/32")	1	1
3		7930360-NLA	NLA - Package, Screw, Phillips Pan Head (#10-32 x 11/16")	25	1
4		1008935	Spring, Release - Right	1	1
4		1008936	Spring, Release - Left	1	1
5		7930359-NLA	NLA - Package, Screw, Phillips Pan Head (#10-32 x 1/2")	25	1
6		1027095	Plug Button, Black (7/8")	1	2
7	A,G	1035888-NLA	NLA - Tube, Footrest Support, Adult, Chrome - Right	1	1
7	A,G	1035889-NLA	NLA - Tube, Footrest Support, Adult, Chrome - Left	1	1
7	A,H	1149413-NLA	NLA - Tube, Footrest Support, Silver Vein - Right	1	1
7	A,H	1149414-NLA	NLA - Tube, Footrest Support, Silver Vein - Left	1	1
7	B,G	1035491-NLA	NLA - Tube, Footrest Support, Hemi, Chrome - Right	1	1
7	B,G	1035490-NLA	NLA - Tube, Footrest Support, Hemi, Chrome - Left	1	1
7	B,H	1149410-NLA	NLA - Tube, Footrest Support, Silver Vein - Right	1	1
7	B,H	1149409-NLA	NLA - Tube, Footrest Support, Silver Vein - Left	1	1
8		11232X009	Bolt, Lug (5/16-24 x 1-3/8")	1	1
9		1107321	Nut, Low Cap (5/16"-24)	1	1
10		1030979	Screw, Phillips Pan Head Tap (#10-32 x 3/4")	1	1
11		1022874	Bumper	1	1
12		67005X000	Rivnut (#10-32)	1	1
13		7930365-NLA	NLA - Locknut (#10-32)	25	1
14	C	1006490	Tip, Plastic, Black	1	1
15	A,D,E, G	1035890-NLA	NLA - Footrest Support Assembly, Adult, Chrome - Right	1	1
15	A,D,E, G	1035891-NLA	NLA - Footrest Support Assembly, Adult, Chrome - Left	1	1
15	A,D,E, H	1149415-NLA	NLA - Footrest Support Assembly, Silver Vein - Right	1	1
15	A,D,E, H	1149416-NLA	NLA - Footrest Support Assembly, Silver Vein - Left	1	1
15	B,D,E, G	1035493-NLA	NLA - Footrest Support Assembly, Hemi, Chrome - Right	1	1
15	B,D,E, G	1035492-NLA	NLA - Footrest Support Assembly, Hemi, Chrome - Left	1	1
15	B,D,E, H	1149412-NLA	NLA - Footrest Support Assembly, Silver Vein - Right	1	1
15	B,D,E, H	1149411-NLA	NLA - Footrest Support Assembly, Silver Vein - Left	1	1

Footrest Support Assembly

93M2 / H93M2

NOTE: Due to design change, for chairs built prior to 10/1/00 it is suggested that Release Handles be replaced in pairs. See category titled, "Front Riggings" for complete Pivot/Slide Tube with Footplate Assembly.

NOTE: A - Front Rigging 93M2 - Adult Chairs.
B - Front Rigging H93M2 - Hemi Chairs.
C - Before 10/1/00 Only.
D - 14-22" Chair Width.
E - Includes items 1-14.
G - Before 7/6/07.
H - After 7/5/07.