


ASBA Seat Frame and Hardware (16"-24" Wide) (Before 7/5/05)



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	See Chart	#1 Seat Frame with Hardware - Non-Recliner	1	1
2		1029501	Plug Button, Black (1")	1	2
3		1090768	Kit, Lock Handle w/Rivets (Pair)	2	1
4		47471X074	Spring, Arm Lock	1	2
5		1028594	Pin, Arm Lock, Black	1	2

NOTE: A - Includes seat frame and items 2-5

#1 Seat Frame with Hardware

Non-Recliner

Chair Width	16"-17" Depth	18" - 19" Depth	20" - 21" Depth	22" Depth
16" Width	1118313-NLA	1118322-NLA	1118331-NLA	1118340-NLA
17" Width	1118314-NLA	1118323-NLA	1118332-NLA	1118341-NLA
18" Width	1118315-NLA	1118324-NLA	1118333-NLA	1118342-NLA
19" Width	1118316-NLA	1118325-NLA	1118334-NLA	1118343-NLA
20" Width	1118317-NLA	1118326-NLA	1118335-NLA	1118344-NLA
21" Width	1118318-NLA	1118327-NLA	1118336-NLA	1118345-NLA
22" Width	1118319-NLA	1118328-NLA	1118337-NLA	1118346-NLA
23" Width	1118320-NLA	1118329-NLA	1118338-NLA	1118347-NLA

ASBA Seat Frame and Hardware (16"-24" Wide) (Before 7/5/05)

#1 Seat Frame with Hardware

Non-Recliner

Chair Width	16"-17" Depth	18" - 19" Depth	20" - 21" Depth	22" Depth
24" Width	1118321-NLA	1118330-NLA	1118339-NLA	1118348-NLA

NOTE: Seat frame includes items 2-5