


Seat Frame Mounting Hardware

Elevate Only For ASBA Seating (Before 7/5/05)



Seat Frame Mounting Hardware

Elevate Only For ASBA Seating (Before 7/5/05)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1125777-NLA	Kit, A.S.B.A. Seat Frame Mounting Hardware	1	1
2			NLA - Screw, Hex Head w/ Patch (5/16-18 x 1-1/2")	1	4
3			NLA - Screw, Hex Head w/ Patch (5/16-18 x 7/8")	1	6
4			Washer, Flat (21/64 x 53/64 x 7/64")	1	6
5			NLA - Bracket, Seat Frame Mounting, Black - Right	1	1
5			NLA - Bracket, Seat Frame Mounting, Black - Left	1	1
6			Locknut (5/16-18)	1	10
7			Screw, Hex Head w/ Patch (5/16-18 x 1-1/4")	1	4
8			NLA - Plate, A.S.B.A. Seat Frame Mounting	1	1
9			NLA - Nut, Strap	1	2
10			Decal, Seat Adjustment Warning	1	2

NOTE: A - Includes items 2-10