


# Back Upholstery, Nylon

## Elevate Only For ASBA/ADJASBA Seating



**Single Liners**

Width	16"-17" High	18"-19" High	20" High	21"-22" High	23"-24" High
16" W	1045175	1045207	1068947	1045159	1054529
17" W	1045176	1045208	1068948	1045160	1054530
18" W	1045177	1045209	1068949	1045161	1054531

# Back Upholstery, Nylon

## Elevate Only For ASBA/ADJASBA Seating

### Dual Liners

Width	16"-17" High	18" - 19" High	20" High	21" - 22" High	23" - 24" High
16" W	1045180	1045212	1068952	1045164	1054534
17" W	1045181	1045213	1068953	1045165	1054535
18" W	1045182	1045214	1068954	1045166	1054536
19" W	1045183	1045215	1068955	1045167	1054537
20" W	1045184	1045216	1068956	1045168	1054538
21" W	1045185	1045217	1068957	1045169	1054539
22" W	1045186	1045218	1068958	1045170	1054540

NOTE: Dual liner back upholstery is used with chairs that have a seat deeper than 17" or back wider than 18".

### Back Upholstery Hardware

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
	A	See Chart	Single Liners	1	1
	A	See Chart	Dual Liners - Nylon	1	1
	A	1026298	Package, Screw with Washer (#10 x 7/8") (Pair)	2	1

NOTE: A - Not Shown