


Adjustable Angle Stroller Handle (After 8/16/05)

Stroller Handle Hardware (item 1) includes hardware only (items 2-8). If a complete Stroller Handle is needed then items 9 & 11 must be ordered as well.



Adjustable Angle Stroller Handle (After 8/16/05)

Stroller Handle Hardware (item 1) includes hardware only (items 2-8). If a complete Stroller Handle is needed then items 9 & 11 must be ordered as well.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1127791	Kit, Adjustable Angle Stroller Handle Hardware	1	1
2			Screw, Allen Head (1/4-28 x 1-1/4")	1	2
3			Locknut (1/4-28)	1	2
4			Locknut (1/4-28)	1	2
5			Screw, Socket Button Head (1/4-28 x 1-1/4")	1	2
6			Folding Hinge Assembly	1	2
7			Pin, Quick Release with Lanyard	1	2
8			NLA - Rivet (3/16 x 7/16")	1	2
9		See Chart	#9 Crossbrace Tube	1	1
10	B	1055398	Grip, Hand, Foam (4")	1	2
10	C	1127767	Grip, Foam (9")	1	2
11	B,D	1134900-NLA	NLA - Stroller Handle Assembly - Right/Left	1	2
11	C,D	1127789	Stroller Handle Assembly - Right/Left	1	2

NOTE: A - Includes items 2-8
 B - 10"-11" Wide chairs
 C - 12"-18" Wide chairs
 D - Includes item 10, Foam Grip

#9 Crossbrace Tube

10-11" Wide	12-13" Wide	14-15" Wide	16-18" Wide
1134896-NLA	1129473	1129474	1129475