


Adjustable Angle Stroller Handle (8/16/05 thru 5/29/08)

Stroller Handle Hardware (item 1) includes hardware only (items 2-8). If a complete Stroller Handle is needed then items 9 & 11 must be ordered as well.



Adjustable Angle Stroller Handle

(8/16/05 thru 5/29/08)

Stroller Handle Hardware (item 1) includes hardware only (items 2-8). If a complete Stroller Handle is needed then items 9 & 11 must be ordered as well.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1127791	Kit, Adjustable Angle Stroller Handle Hardware	1	1
2			Screw, Allen Head (1/4-28 x 1-1/4")	1	2
3			Locknut (1/4-28)	1	2
4			Locknut (1/4-28)	1	2
5			Screw, Socket Button Head (1/4-28 x 1-1/4")	1	2
6			Folding Hinge Assembly	1	2
7			Pin, Quick Release with Lanyard	1	2
8			NLA - Rivet (3/16 x 7/16")	1	2
9		See Chart	#9 Crossbrace Tube	1	1
10		1127767	Grip, Foam (9")	1	2
11	B	1127789	Stroller Handle Assembly - Right/Left	1	2

NOTE: A - Includes items 2-8
B - Includes item 10, Foam Grip

#9 Crossbrace Tube

14"-15" Wide	16"-17" Wide	18"-22"
1129474	1129475	1129476