



Non-Removable Fold Down Side Guards



Non-Removable Fold Down Side Guards

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1114099	NLA - Screw, Hex Head (1/4-28 x 2")	1	8
2		1039232-NLA	NLA - Cam, Adjustment	1	2
3		1071496-NLA	NLA - Pin, Roll (3/16 x 3/4")	1	2
4		1114084-NLA	NLA - Bracket, Support Back, Black -Right	1	1
4		1114085-NLA	NLA - Bracket, Support Back, Black -Left	1	1
5		1058636-NLA	NLA - Washer, Nylon (25/64 x 7/8 x 1/32")	1	2
6		1108477-NLA	NLA - Washer, Nylon (1/4 x 3/4 x 5/8")	1	4
7		1114087-NLA	NLA - Guard, Side, Folding, Black	1	2
8		1026162	Spacer (1/4 x 3/8 x 11/32")	1	4
9		1040776-NLA	NLA - Spacer, Nylon (1/4 x 1 x 1/4")	1	4
10		40157X020	Spacer, Coved (1/4 x 1 x 9/16")	1	4
11		1026907	Washer (1/4 x 1/2 x 1/16")	1	14
12		1116005-NLA	NLA - Bracket, Latch, Black	1	2
13		1039233	Spacer (1/4 x 7/16 x 5/8")	1	2
14		1114077-NLA	NLA - Latch, Torsion Spring - Right	1	1
14		1114078-NLA	NLA - Latch, Torsion Spring - Left	1	1
15		1025728	Washer (1/4 x 1/2 x 1/32")	1	2
16		1025765	Locknut (1/4-28)	1	10
17		1141633	Cord, Nylon (6-1/2 ft.)	1	1
18		1012371	Washer, Nylon (3/8 x 7/8 x 1/16")	1	8
19	A	1114097-NLA	NLA - Kit, Fold Down Back Mounting Hardware with Side Guards	2	1
20		1114099	NLA - Screw, Hex Head (1/4-28 x 2")	1	4
21		1108477-NLA	NLA - Washer, Nylon (1/4 x 3/4 x 5/8")	1	8
22		1114087-NLA	NLA - Guard, Side, Folding, Black	1	2
23		40157X020	Spacer, Coved (1/4 x 1 x 9/16")	1	4
24		1026907	Washer (1/4 x 1/2 x 1/16")	1	4
25		1025765	Locknut (1/4-28)	1	4
26	B	1114149-NLA	NLA - Kit, Fixed Back Mounting Hardware with Side Guards	2	1
NOTE: A - Includes items 1-18 B - Includes items 20-25					