



Pivot Tube Assemblies w/ Footplates - 1651 Composite



**Used on all front riggings except
Wedge Style Elevating Legrests AT5044**



Used on Wedge Style Elevating Legrests AT5044

Pivot Tube Assemblies w/ Footplates - 1651 Composite

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		See Chart	#1 60° Footrest (13-19")	1	1
1		See Chart	#1 70° Non Tapered Footrest (14-19")	1	1
1		See Chart	#1 70° Tapered Footrest (14-1/2-19-1/2")	1	1
1		See Chart	#1 P30E Elevating Legrest - Ratchet Style (16-20")	1	1
1		See Chart	#8 AT5044 Elevating Legrest - Wedge Style (16-20")	1	1
2	A	1120382	Tube, Pivot/Slide, Curved End, Black - Long (9-7/8")	1	1
2	A	1120386	Tube, Pivot/Slide, Curved, 2" Longer, Black (11-7/8")	1	1
2	B	1120384-NLA	NLA - Tube, Pivot/Slide, Straight - Alum. (9-7/8" L)	1	1
2	B	1120385	Tube, Pivot/Slide, Straight - Black (11-7/8" L)	1	1
3		1134426	Package, Plug Button, Black (3/4")	10	1
4		1026879	Footplate, Composite, Small, Black (6" W x 6" L)	1	1
4		1026880	Footplate, Composite, Medium, Black (7" W x 6" L)	1	1
4		1026881	Footplate, Composite, Large, Black (8" W x 6" L)	1	1
5		1025195	Locknut (1/4-20)	1	1
6		60285X006	Screw, Socket Head (1/4-20 x 2-1/4")	1	1
7		1025457	Spacer, Coved (1/4 x 3/4 x 11/32")	1	2
8	E	See Chart	#8 AT5044 Elevating Legrest - Wedge Style (16-20")	1	1
9		1134427	Package, Plug Button, Black, "D" Shaped (3/4")	10	1
10	C	1134428	Pivot/Slide Tube w/ D Insert, Black	1	1
10	C	1134429	Pivot/Slide Tube w/ D Insert, Black (2" Longer)	1	1
11		1026879	Footplate, Composite, Small, Black (6" W x 6" L)	1	1
11		1026880	Footplate, Composite, Medium, Black (7" W x 6" L)	1	1
11		1026881	Footplate, Composite, Large, Black (8" W x 6" L)	1	1
12		1025195	Locknut (1/4-20)	1	1
13		60285X006	Screw, Socket Head (1/4-20 x 2-1/4")	1	1
14		1025457	Spacer, Coved (1/4 x 3/4 x 11/32")	1	2
NOTE: A - 60° Footrests, and P30E Ratchet Style Legrests B - 70° Footrests C - AT5044 Wedge Style Legrests E - Includes items 9-14					

#1 60° Footrest 13-19" Curved End

Right/Left	Small Long 14" wide	Small 2" Longer 14" wide	Medium Long 16" wide	Medium 2" Longer 16" wide	Large Long 17-20" wide	Large 2" Longer 17-20" wide
Right	1028391	1032627-NLA	1028387	1032625-NLA	1028383	1032623-NLA
Left	1028392	1032628-NLA	1028388	1032626-NLA	1028384	1032624-NLA

Pivot Tube Assemblies w/ Footplates - 1651 Composite

#1 70° Non Tapered Footrest

14-19"

Straight

Right/Left	Small Short 14" wide	Small 2" Longer 14" wide	Medium Short 16" wide	Medium 2" Longer 16" wide	Large Short 17-20" wide	Large 2" Longer 17-20" wide
Right	1042784-NLA	1042790	1042786	1042792	1042788	1042794
Left	1042785-NLA	1042791	1042787	1042793	1042789	1042795

#1 70° Tapered Footrest

14-1/2 - 19-1/2"

Straight

Right/Left	Small Short 16" wide	Small 2" Longer 16" wide	Medium Short 17-18" wide	Medium 2" Longer 17-18" wide	Large Short 20" wide	Large 2" Longer 20" wide
Right	1042784-NLA	1042790	1042786	1042792	1042788	1042794
Left	1042785-NLA	1042791	1042787	1042793	1042789	1042795

#1 P30E Elevating Legrest - Ratchet Style

16-20"

Curved End

Right/Left	Small Long 14" wide	Small 2" Longer 14" wide	Medium Long 16" wide	Medium 2" Longer 16" wide	Large Long 17-20" wide	Large 2" Longer 17-20" wide
Right	1028391	1032627-NLA	1028387	1042792	1028383	1042794
Left	1028392	1032628-NLA	1028388	1042793	1028384	1042795

Pivot Tube Assemblies w/ Footplates - 1651 Composite

#8 AT5044 Elevating Legrest - Wedge Style

16-20"

Curved End

Right/Left	Small Long 14" wide	Small 2" Longer 14" wide	Medium Long 16" wide	Medium 2" Longer 16" wide	Large Long 17-20" wide	Large 2" Longer 17-20" wide
Right	1132078	1132107	1132080	1132109	1132082	1132111
Left	1132079	1132108	1132081	1132110	1132083	1132112