



Footplate & Pivot Tube Assemblies - 1651 Composite



Used on all front riggings except
Wedge Style Elevating Legrests AT5044



Used on Wedge Style Elevating Legrests AT5044

Footplate & Pivot Tube Assemblies - 1651 Composite

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	E	See Chart	#1 60° Footrest (13-19")	1	1
1	E	See Chart	#1 60° MFX Footrest (9-3/4-14")	1	1
1	E	See Chart	#1 70° Non Tapered Footrest (14-19")	1	1
1	E	See Chart	#1 AT5044 Elevating Legrest - Ratchet Style (16-20")	1	1
2	A	1120383	Tube, Pivot/Slide, Curved End, Black - Short (8-1/4")	1	1
2	B	1120382	Tube, Pivot/Slide, Curved End, Black - Long (9-7/8")	1	1
2	C	1120384-NLA	NLA - Tube, Pivot/Slide, Straight - Alum. (9-7/8" L)	1	1
3		1134426	Package, Plug Button, Black (3/4")	10	1
4		1026879	Footplate, Composite, Small, Black (6" W x 6" L)	1	1
4		1026880	Footplate, Composite, Medium, Black (7" W x 6" L)	1	1
4		1026881	Footplate, Composite, Large, Black (8" W x 6" L)	1	1
5		1025195	Locknut (1/4-20)	1	1
6		60285X006	Screw, Socket Head (1/4-20 x 2-1/4")	1	1
7		1025457	Spacer, Coved (1/4 x 3/4 x 11/32")	1	2
8	F	See Chart	#8 AT5044 Elevating Legrest - Wedge Style (16-20")	1	1
9		1134427	Package, Plug Button, Black, "D" Shaped (3/4")	10	1
10	D	1134428	Pivot/Slide Tube w/ D Insert, Black	1	1
10	D	1134429	Pivot/Slide Tube w/ D Insert, Black (2" Longer)	1	1
11		1026879	Footplate, Composite, Small, Black (6" W x 6" L)	1	1
11		1026880	Footplate, Composite, Medium, Black (7" W x 6" L)	1	1
11		1026881	Footplate, Composite, Large, Black (8" W x 6" L)	1	1
12		1025195	Locknut (1/4-20)	1	1
13		60285X006	Screw, Socket Head (1/4-20 x 2-1/4")	1	1
14		1025457	Spacer, Coved (1/4 x 3/4 x 11/32")	1	2

NOTE: A - 60° MFX Footrests
 B - 60° Footrests, and AT5044 Ratchet Style Legrests
 C - 70° Footrests
 D - AT5044 Wedge Style Legrests
 E - Includes items 2-7
 F - Includes items 9-14

#1 60° Footrest 13-19" Curved End

Right/Left	Small Long 14" wide	Medium Long 15-16" wide
Right	1028391	1028387
Left	1028392	1028388

Footplate & Pivot Tube Assemblies - 1651 Composite

**#1 60° MFX Footrest
9-3/4 - 14"
Curved End**

Right/Left	Small Short 14" wide	Medium Short 15-16" wide
Right	1028389-NLA	1028385-NLA
Left	1028390-NLA	1028386-NLA

**#1 70° Non Tapered
14-19"
Straight**

Right/Left	Small Short 13-14" wide	Medium Short 15-16" wide
Right	1042784-NLA	1042786
Left	1042785-NLA	1042787

**#1 AT5044 Elevating Legrest - Ratchet Style
16-20"
Curved End**

Right/Left	Small Long 14" wide	Small 2" Longer 14" wide	Medium Long 15-16" wide	Medium 2" Longer 15-16" wide	Large Long 17-22" wide	Large 2" Longer 17-22" wide
Right	1028391	1120416	1028387	1120418	1028383	1120420
Left	1028392	1120417	1028388	1120419	1028384	1120421

Footplate & Pivot Tube Assemblies - 1651 Composite

#8 AT5044 Elevating Legrest - Wedge Style

16-20"

Curved End

Right/Left	Small Long 14" wide	Medium Long 15-16" wide
Right	1132078	1132080
Left	1132079	1132081