



# Footplate and Pivot Tube Assemblies - 1651 Composite



**Used on all front riggings except  
Wedge Style Elevating Legrests AT5044**



**Used on Wedge Style Elevating Legrests AT5044**

# Footplate and Pivot Tube Assemblies - 1651 Composite

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	F	See Chart	#1 60° Footrest (13-19")	1	1
1	F	See Chart	#1 60° MFX Footrest (9-3/4-14")	1	1
1	F	See Chart	#1 70° Non Tapered Footrest (14-19")	1	1
1	F	See Chart	#1 70° Tapered Footrest (14-1/2-19-1/2")	1	1
1	F	See Chart	#1 AT5044 Elevating Legrest - Ratchet Style (16-20")	1	1
2	A	1120382	Tube, Pivot/Slide, Curved End, Black - Long (9-7/8")	1	1
2	B	1120389	Tube, Pivot/Slide, Straight, Black - Long (9-7/8")	1	1
2	C	1120390	Tube, Pivot/Slide, Straight, 2" Longer ALPT, Black (11-7/8")	1	1
2	D	1120391	Tube, Pivot/Slide, Curved End, 2" Longer ALPT, Black (11-7/8")	1	1
3		1134426	Package, Plug Button, Black (3/4")	10	1
4		1026879	Footplate, Composite, Small, Black (6" W x 6" L)	1	1
4		1026880	Footplate, Composite, Medium, Black (7" W x 6" L)	1	1
4		1026881	Footplate, Composite, Large, Black (8" W x 6" L)	1	1
5		1025195	Locknut (1/4-20)	1	1
6		60285X006	Screw, Socket Head (1/4-20 x 2-1/4")	1	1
7		1025457	Spacer, Coved (1/4 x 3/4 x 11/32")	1	2
8	G	See Chart	#8 AT5044 Elevating Legrest - Wedge Style (16-20")	1	1
9		1134427	Package, Plug Button, Black, "D" Shaped (3/4")	10	1
10	E	1134428	Pivot/Slide Tube w/ D Insert, Black	1	1
10	E	1134429	Pivot/Slide Tube w/ D Insert, Black (2" Longer)	1	1
11		1026879	Footplate, Composite, Small, Black (6" W x 6" L)	1	1
11		1026880	Footplate, Composite, Medium, Black (7" W x 6" L)	1	1
11		1026881	Footplate, Composite, Large, Black (8" W x 6" L)	1	1
12		1025195	Locknut (1/4-20)	1	1
13		60285X006	Screw, Socket Head (1/4-20 x 2-1/4")	1	1
14		1025457	Spacer, Coved (1/4 x 3/4 x 11/32")	1	2

NOTE: A - 60° Footrests, and AT5044 Ratchet Style Legrests  
 B - 70° Footrests  
 C - 70° Footrests, and AT5044 Ratchet Style Legrests  
 D - 60° Footrest  
 E - AT5044 Wedge Style Legrests  
 F - Includes items 2-7  
 G - Includes items 9-14

## #1 60° Footrest 13-19" Curved End

Right/Left	Small Long 13-14" wide	Small 2" Longer 13-14" wide	Medium Long 15-16" wide	Medium 2" Longer 15-16" wide	Large Long 17-22" wide	Large 2" Longer 17-22" wide
Right	1028391	1120416	1028387	1120418	1028383	1120420
Left	1028392	1120417	1028388	1120419	1028384	1120421

# Footplate and Pivot Tube Assemblies - 1651 Composite

## #1 60° MFX Footrest

9-3/4 - 14"

Curved End

Right/Left	Small Short 13-14" wide	Medium Short 15-16" wide	Large Short 17-22" wide
Right	1028391	1028387	1028383
Left	1028392	1028388	1028384

## #1 70° Non Tapered; 70LIFTNT Footrest

14-19"

Straight

Right/Left	Small Short 13-14" wide	Small 2" Longer 13-14" wide	Medium Short 15-16" wide	Medium 2" Longer 15-16" wide	Large Short 17-22" wide	Large 2" Longer 17-22" wide
Right	1120404	1120410	1120406	1120412	1120408	1120414
Left	1120405	1120411	1120407	1120413	1120409	1120415

## #1 70° Tapered Footrest

14-1/2 - 19-1/2"

Straight

Right/Left	Small Short 16" wide	Small 2" Longer 16" wide	Medium Short 17-18" wide	Medium 2" Longer 17-18" wide	Large Short 19-22" wide	Large 2" Longer 19-22" wide
Right	1120404	1120410	1120406	1120412	1120408	1120414
Left	1120405	1120411	1120407	1120413	1120409	1120415

# Footplate and Pivot Tube Assemblies - 1651 Composite

## #1 AT5044 Elevating Legrest - Ratchet Style

16-20"

Curved End

Right/Left	Small Long 14" wide	Small 2" Longer 14" wide	Medium Long 15-16" wide	Medium 2" Longer 15-16" wide	Large Long 17-22" wide	Large 2" Longer 17-22" wide
Right	1028391	1120416	1028387	1120418	1028383	1120420
Left	1028392	1120417	1028388	1120419	1028384	1120421

## #8 AT5044 Elevating Legrest - Wedge Style

16-20"

Curved End

Right/Left	Small Long 14" wide	Small 2" Longer 14" wide	Medium Long 15-16" wide	Medium 2" Longer 15-16" wide	Large Long 17-22" wide	Large 2" Longer 17-22" wide
Right	1132078	1132107	1132080	1132109	1132082	1132111
Left	1132079	1132108	1132081	1132110	1132083	1132112