



# Growth Kits



**Growth Kits for chairs with 12" Wheels**



**Growth Kits for chairs with 20", 22", or 24" Wheels**

# Growth Kits

## #1 Growth Kits - Flush Mounted Seat Pans (15"-17" Deep only) (chairs with 12" Wheels)

Kits include Backrest Crossmember, Insert Tubes, Spacer Tubes, Seat Pan, Seat Pan Supports, Push Handle Assembly, and Foot Release Assembly for wheel lock.

Desired Chair Width	Chairs with Chrome Spacer Tubes	Chairs with Color Spacer Tubes
16" W	1140326-NLA See Note(s) B	8881140327 See Note(s) B, A
18" W	1140330-NLA See Note(s) B	8881140331 See Note(s) B, A
20" W	1140322-NLA See Note(s) B	8881140323 See Note(s) B, A
NOTE: A - Specify frame finish when ordering B - Kits may be used for more than 1 seat width size depending upon Seat Pan type.		

## #1 Growth Kits - Drop Mounted Seat Pans (15"-17" Deep only) (chairs with 12" Wheels)

Kits include Backrest Crossmember, Insert Tubes, Spacer Tubes, Seat Pan, Push Handle Assembly, and Foot Release Assembly for wheel lock (Seat Pan Supports are provided but NOT required).

Desired Chair Width	Chairs with Chrome Spacer Tubes	Chairs with Color Spacer Tubes
18" W	1140326-NLA See Note(s) B	8881140327 See Note(s) B, A
20" W	1140330-NLA See Note(s) B	8881140331 See Note(s) A, B
22" W	1140322-NLA See Note(s) B	8881140323 See Note(s) B, A
NOTE: A - Specify frame finish when ordering B - Kits may be used for more than 1 seat width size depending upon Seat Pan type.		

## #2 Growth Kits - Flush Mounted Seat Pans (15"-17" Deep only) (chairs with 20", 22", & 24" Wheels)

Kits include Backrest Crossmember, Insert Tubes, Spacer Tubes, Seat Pan, Seat Pan Supports, Spacer and Slide Tube for stroller handle, and Foam Grip.

Desired Chair Width	Chairs with Chrome Spacer Tubes	Chairs with Color Spacer Tubes
16" W	1140328-NLA See Note(s) B	8881140329 See Note(s) B, A
18" W	1140332-NLA See Note(s) B	8881140333 See Note(s) A, B
20" W	1140324-NLA See Note(s) B	8881140325 See Note(s) B, A
NOTE: A - Specify frame finish when ordering B - Kits may be used for more than 1 seat width size depending upon Seat Pan type.		

# Growth Kits

## #2 Growth Kits - Drop Mounted Seat Pans (15"-17" Deep only) (chairs with 20", 22", & 24" Wheels)

Kits include Backrest Crossmember, Insert Tubes, Spacer Tubes, Seat Pan, Spacer and Slide Tube for stroller handle, and Foam Grip (Seat Pan Supports are provided but NOT required).

Desired Chair Width	Chairs with Chrome Spacer Tubes	Chairs with Color Spacer Tubes
18" W	1140328-NLA See Note(s) B	8881140329 See Note(s) B, A
20" W	1140332-NLA See Note(s) B	8881140333 See Note(s) A, B
22" W	1140324-NLA See Note(s) B	8881140325 See Note(s) B, A
NOTE: A - Specify frame finish when ordering B - Kits may be used for more than 1 seat width size depending upon Seat Pan type.		