


MNS500 Low Air Loss System

Mattress Covers are not shown.



MNS500 Low Air Loss System

Mattress Covers are not shown.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1135278-NLA	NLA - Cover, Mattress Top MNS500 - Narrow (STZ)	1	1
1	A	1135279-NLA	NLA - Cover, Mattress Top MNS500 - Standard (STZ)	1	1
1	A	1135280-NLA	NLA - Cover, Mattress Top MNS500 - Bariatric (STZ)	1	1
2	A	1132042-NLA	NLA - Cover, Mattress Top MNS500 - Narrow (LTZ)	1	1
2	A	1132043-NLA	NLA - Cover, Mattress Top MNS500 - Standard (LTZ)	1	1
2	A	1132044-NLA	NLA - Cover, Mattress Top MNS500 - Bariatric (LTZ)	1	1
3		1141540-NLA	NLA - Connector, Male CPC (1/4" NPT)	1	1
4		1141539-NLA	NLA - Connector, Female CPC 1/4" NPT	1	1
5		1141543-NLA	NLA - Fuse (750 M AMP)	5	1
6		1141587-NLA	NLA - Package, Vibration Mount	4	1
7	B	1141522-NLA	NLA - Package, Thomas Filter	10	1
8		1141542-NLA	NLA - Muffler, Brass (1/8" NPT)	1	1
9		1141533-NLA	NLA - Package, Solenoid Valve	2	1
10		1141561-NLA	NLA - Insert, Coupling (1/2" NPT)	1	1
11	C	See Note	See Footnote	1	1
12		1146109-NLA	NLA - Bladder, Air Sack, MNS5000 Narrow (32" Long)	1	1
12		1143681-NLA	NLA - Bladder, Air Sack, MNS5000 Standard (35" Long)	1	1
12		1143682-NLA	NLA - Bladder, Air Sack, MNS5000 Bariatric (42" Long)	1	1

NOTE: A - Not shown
 B - Photo shows 1(one) filter, part # contains a quantity of 10 (ten) Filters.
 C - No longer available as a service part in the field.