


MNS600 System

Mattress Covers are not shown.



MNS600 System

Mattress Covers are not shown.

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1135281-NLA	NLA - Cover, Mattress Top MNS600 - Narrow (STZ)	1	1
1	A	1135282-NLA	NLA - Cover, Mattress Top MNS600 - Standard (STZ)	1	1
1	A	1135283-NLA	NLA - Cover, Mattress Top MNS600 - Bariatric (STZ)	1	1
2	A	1132045-NLA	NLA - Cover, Mattress Top MNS600 - Narrow (LTZ)	1	1
2	A	1132046-NLA	NLA - Cover, Mattress Top MNS600 - Standard (LTZ)	1	1
2	A	1132047-NLA	NLA - Cover, Mattress Top MNS600 - Bariatric (LTZ)	1	1
3		1141540-NLA	NLA - Connector, Male CPC (1/4" NPT)	1	1
4		1141539-NLA	NLA - Connector, Female CPC 1/4" NPT	1	1
5		1141543-NLA	NLA - Fuse (750 M AMP)	5	1
6		1141587-NLA	NLA - Package, Vibration Mount	4	1
7		1141542-NLA	NLA - Muffler, Brass (1/8" NPT)	1	1
8		1141533-NLA	NLA - Package, Solenoid Valve	2	1
9		1141561-NLA	NLA - Insert, Coupling (1/2" NPT)	1	1
10		1141541-NLA	NLA - Fitting, Elbow (1/4" Barb x 1/8" NPT)	2	1
11		1141572-NLA	NLA - Solenoid Manifold	1	1
12		1141574-NLA	NLA - Compressor - Thomas	1	1
13		1143684-NLA	NLA - Tubing, Vinyl Double (25' L)	1	1
14		1146106-NLA	NLA - Bladder, Air Sack, MNS600 Narrow (32" Long)	1	1
14		1146107-NLA	NLA - Bladder, Air Sack, MNS600 Standard (35" Long)	1	1
14		1146108-NLA	NLA - Bladder, Air Sack, MNS600 Bariatric (42" Long)	1	1

NOTE: A - Not shown