





# Miscellaneous Accessories



**TRRO Pelvic Restraint Belt  
Mounting Location**



**Auto Style Push Button Buckle  
Positioning Strap for TRBKTS Option**



**Seat Positioning Strap  
Mounting Hardware**

## Miscellaneous Accessories

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1122818-NLA	NLA - Strap, Positioning, Seat, Airline Style Buckle w/ Hardware (37" Long Eyelet-to-Eyelet)	1	1
1	B	1036527	Strap, Positioning, Seat/Chest w/ Auto Style Buckle (Eyelet to Eyelet Length 33")	1	1
2		1118569-NLA	NLA - Pin, Threaded	1	2
3		1118568-NLA	NLA - Bushing, Nylon (5/16 x 7/16 x 5/32")	1	2
4		1118579-NLA	NLA - Nut, Hex Jam (3/8-16)	1	4
5	C	1124385	Kit, Seat Positioning Strap Mounting Hardware (both LT & RT sides)	1	1
6			Screw, Hex Head w/ Patch (1/4-20 x 3/4")	1	2
7			Washer (1/4 x 23/32 x 1/16")	1	2
8			Locknut (1/4-20)	1	2
	D	40176X010-NLA	NLA - Pump, Tire	1	1

NOTE: A - Pelvic Restraint Belt for TRRO option. Includes items 2-4  
 B - Positioning Strap for TRBKTS option.  
 C - Hardware for auto style positioning strap. Includes items 6-8  
 D - Not Shown