


ADJASBA Seat Frame Mounting Hardware (16"-24" Wide) Low S-T-F TRRO/TRBKTS Option



ADJASBA Seat Frame Mounting Hardware (16"-24" Wide) Low S-T-F TRRO/TRBKTS Option

| Item Number | Note | Part Number | Description | Quantity Per | |
|-------------|------|-------------|---|--------------|----------|
| | | | | Package | Assembly |
| 1 | A | 1139542 | Kit, Seat Frame Mtg. Hdwr. 16"-24" W, Low S-T-F with TRRO/TRBKTS Option | 1 | 1 |
| 2 | | | Screw, Socket Head w/ Patch (5/16-18 x 7/8") | 1 | 4 |
| 3 | | | NLA - Washer, Ground (11/32 x 13/16 x 1/8") | 1 | 6 |
| 4 | | | NLA - Screw, Socket Head (5/16-18 x 3/4") | 1 | 2 |
| 5 | | | Plate, Low Seat-To-Floor Mounting, Rear TRRO/TRBKTS - Right | 1 | 1 |
| 5 | | | Plate, Low Seat-To-Floor Mounting, Rear TRRO/TRBKTS - Left | 1 | 1 |
| 6 | | | Plate, Nut (5/16-18) | 1 | 2 |
| 7 | | | Washer, Flat (21/64 x 53/64 x 7/64") | 1 | 4 |
| 8 | | | Screw, Hex Head w/ Patch (5/16-18 x 3/4") | 1 | 8 |
| 9 | | | NLA - Nut, Strap | 1 | 2 |
| 10 | | | NLA - Bracket, Upper Seat, Front, Black | 1 | 2 |
| 11 | | | NLA - Bracket, Lower Seat, Front, Black | 1 | 2 |

NOTE: A - Includes items 2-11