

Invacare® Pegasus METRO/ PRO **Canadian Power Scooter**



Yes, you can.

Parts Catalogue Usage Guide

The information contained in this document is subject to change without notice. This catalogue contains service parts for Form No. 18-075C. If you have any questions regarding this catalogue, please contact one of our Customer Service Representatives.

For Customer Service Support Call:

1-800-668-5324 or 1-905-890-8300

For Faxed In Orders Call:

1-800-668-5478 or

1-905-890-5244

For Technical Service Call:

1-800-668-5324 x2655

For Orders:

canadaorders@invacare.com

Literature: www.invacare.ca or 1-800-668-5324

SYMBOL DEFINITION

Throughout this catalog you will notice letters to draw your attention to specific footnotes which are located at the bottom of the page or below specified charts.

Examples:

A,B,C etc. = Special Comments

NOTE: = General information pertaining to that specific page

NLA = No Longer Available

HARDWARE MEASUREMENTS

The measurements for the hardware will be listed as follows:

Screws, Bolts Diameter X Thread Count/Pitch X Length

Washers Inside Diameter X Outside Diameter X Thickness

Spacers Inside Diameter X Outside Diameter X Length

For your convenience, we have supplied an area for you to list your Invacare Representatives.

Name	Phone Numbers

© 2022 Invacare Corporation. All rights reserved. Reproduction, duplication or modification in whole or in part is prohibited without prior written permission from Invacare. Trademarks are identified by ™ and ®. All trademarks are owned by or licensed to Invacare Corporation or its subsidiaries unless otherwise noted.

TABLE of CONTENTS

Description	Page Number
Anti-tipper Wheel	40
Armrests	34
Batteries and Cables	63
Battery Charger	65
Bumpers	80
Colour Panels (6 Pieces)	27
Electronic Components	66
Front Brake 3-Wheel - Metro	45
Front Brake 3-Wheel - Pro	47
Front Brake 4-Wheel - Metro/ Pro	49
Front Fork 3-Wheel - Metro	55
Front Fork 3-Wheel - Pro	57
Front Frame 3-Wheel - Metro	5
Front Frame 3-Wheel - Pro	9
Front Frame 4-Wheel - Metro	13
Front Frame 4-Wheel - Pro	17
Front Shroud 3-Wheel - Metro	7
Front Shroud 3-Wheel - Pro	11
Front Shroud 4-Wheel - Metro	15
Front Shroud 4-Wheel - Pro	19
Front Suspension 4-Wheel - Metro	59


TABLE of CONTENTS

Description	Page Number
Front Suspension 4-Wheel - Pro	61
Light System	29
Other Accessories	82
Powerbox and Limit Switch Trigger	68
Propulsion Metro 12 km/h	41
Propulsion Pro 15 km/h	43
Rear Bumper Shroud	25
Rear Shroud	23
Rear Swing	21
Rollator Bracket	76
Safety	70
Seat Bracket	32
Seat Post	31
Seating	72
Seats	35
Slide Assembly	33
Steering Complete	53
Steering Head	51
Storage	78
Walking Support	74
Wheels - 5" Rim	36

TABLE of CONTENTS

Description	Page Number
Wheels - 6" Rim	38

Front Frame 3-Wheel - Metro




SCP_093A

Front Frame 3-Wheel - Metro

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1638663	Frame, cpl. - Metro 3-Wheel	--	1

Front Shroud 3-Wheel - Metro




SCP_094A

Front Shroud 3-Wheel - Metro

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1638851	Front Shroud - Metro 3-Wheel	--	1

Front Frame 3-Wheel - Pro




SCP_095A

Front Frame 3-Wheel - Pro

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1638664	Frame, cpl. - Pro 3-Wheel	--	1

Front Shroud 3-Wheel - Pro




SCP_096A

Front Shroud 3-Wheel - Pro

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1638852	Front Shroud - Pro 3-Wheel	--	1

Front Frame 4-Wheel - Metro




SCP_097A

Front Frame 4-Wheel - Metro

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1638665	Frame, cpl. - Metro 4-Wheel	--	1

Front Shroud 4-Wheel - Metro




SCP_098A

Front Shroud 4-Wheel - Metro

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1638853	Front Shroud Metro 4-Wheel	--	1

Front Frame 4-Wheel - Pro




SCP_099A

Front Frame 4-Wheel - Pro

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1638666	Frame, cpl. - Pro 4-Wheel	--	1

Front Shroud 4-Wheel - Pro




SCP_100A

Front Shroud 4-Wheel - Pro

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1638854	Front Shroud - Pro 4-Wheel	--	1

Rear Swing




SCP_101A

Rear Swing

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1638668	Rear Swing, cpl.	--	1
2		1638669	Shock absorber, cpl.	--	1


Rear Shroud



Rear Shroud

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1638971	Rear Shroud, cpl. - Metro (w/o model labels)	--	1
1		1638972	Rear Shroud, cpl. - Pro (w/o model labels)	--	1
2		1638844	Fender, cpl. - right	--	1
3		1638845	Fender, cpl. - left	--	1
	A	1637867	Label Kit - Metro	--	1
	A	1637868	Label Kit - Pro	--	1
NOTE: A - No Shown.					

Rear Bumper Shroud




SCP_103A


Rear Bumper Shroud

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1638670	Cover, bumper rear (w/o lights)	--	1

Colour Panels (6 Pieces)




SCP_104A



Colour Panels (6 Pieces)

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1608821	Colour panels - Marble White (Classic Line)	--	1
2		1608822	Colour panels - Jade Green (Classic Line)	--	1
3		1608823	Colour panels - Sapphire Blue (Classic Line)	--	1
4		1608831	Colour panels - Diamond Silver (Classic Line)	--	1
5		1608832	Colour panels - Jasper Red (Classic Line)	--	1
6		1637415	Colour panels - Luxor Beige (Premium Line)	--	1
7		1637416	Colour panels - Crystal Copper (Premium Line)	--	1
8		1637417	Colour panels - Carmine Red (Premium Line)	--	1
9		1637418	Colour panels - Electric Blue (Premium Line)	--	1
10		1637419	Colour panels - Ionic Chrome (Premium Line)	--	1

Light System




SCP_105A

Light System

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1527473	Head Light, Complete	--	1
2		1638865	Turn signal 3-Wheel, cpl. - left	--	1
3		1638866	Turn signal 3-Wheel, cpl. - right	--	1
4		1518034	Rear Light incl. Turn Signal - LH	--	1
5		1518035	Brake Light, cpl.	--	1
6		1518033	Rear Light incl. Turn Signal - RH	--	1
7		1522030	Bulb, Turn Signal (24V/10 W)	--	4
8		1522034	Bulb, Rear Light (24V/5W)	--	1
9		1517733	Turn Signal - Left	--	1
10		1517734	Turn Signal - Right	--	1
11		1527474	Bulb, Head Light (12V/15W)	--	2


Seat Post



SCP_106A

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1521691	Seat Post Bracket, Standard	--	1


Seat Bracket



SCP_107A

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1636855	Seat Bracket Assembly, cpl.	--	1



Slide Assembly



SCP_108A

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1521690	Slide Assembly	--	1

Armrests




SCP_109A

Armrests


Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1636857	Handrail - RH	--	1
2		1636858	Handrail - LH	--	1

Seats



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1586116	Seat, Standard (20") w/ Double White Stitching	--	1
1		1586115	Seat, HD (26"), Double White Stitching	--	1
2		1529380	Headrest Standard,Seat	--	1
2		1602195	Headrest HD Seat	--	1

Wheels - 5" Rim




SCP_111A

Wheels - 5" Rim

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1586108	Front Wheel 11", pneumatic, 5" rim, black tire - 3-Wheel Metro	--	1
1		1608680	Front Wheel 11", puncture proof, 5" rim, black tire - 3-Wheel	--	1
2		1586111	Metro Front Wheel 11", pneumatic, 5" rim, black tire - 4-Wheel	--	1
2		1608681	Metro Front Wheel 11", puncture proof, 5" rim, black tire - 4-Wheel	--	1
2		1608593	Pro Front Wheel 12", pneumatic, 5" rim, black tire - 4-Wheel	--	1
2		1608683	Pro Front Wheel 12", puncture proof, 5" rim, black tire - 4-Wheel	--	1
3		1586112	Pro Rear Wheel 11", pneumatic, 5" rim, black tire Metro	--	1
3		1608697	Rear Wheel 11", puncture proof, 5" rim Metro	--	1
3		1637428	Rear Wheel 12", pneumatic, 5" rim, RH Pro	--	1
3		1637430	Rear Wheel 12", puncture proof, 5" rim, RH Pro	--	1
4		1608674	Rear Wheel 12", pneumatic, 5" rim, LH Pro	--	1
4		1608773	Rear Wheel 12", puncture proof, 5" rim, LH Pro	--	1
5		1638780	5" Rim Kit, Front - 3-Wheel Metro	--	1
5		1638862	5" Rim Kit, Front - 3-Wheel Pro	--	1
6		1638863	5" Rim Kit, Front - 4-Wheel Metro & Pro	--	1
7		1638838	5" Rim Kit, Rear Metro LH/RH & Pro RH	--	1
8		1638864	5" Rim Kit, Rear Pro LH	--	1
9		1586109	Tire with Tube 11" (3.50/4.10-5) - Black Metro	--	1
9		1586110	PU Tire (puncture proof) 11" (3.50/4.10-5) - black Metro	--	1
9		1640319	Tire (pneumatic) 12" (4.10/3.50-5) - black Pro	--	1
9		1638976	PU Tire (puncture proof) 12" (3.50-5) - black Pro	--	1
10		1527606	Tube 11" (3.50/4.10-5) Metro	--	1
10		1640320	Tube 12" (4.10/3.50-5) Pro	--	1

Wheels - 6" Rim




SCP_112A

Wheels - 6" Rim

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1608776	Front Wheel 11", pneumatic, Low Profile, 6" rim, Orion Metro 3-Wheel	--	1
1		1608794	Front Wheel 11", pneumatic, Low Profile, 6" rim, Orion Metro 4-Wheel	--	1
1		1608799	Rear Wheel 11", pneumatic, Low Profile, 6" rim - Orion Metro	--	1
1		1608796	Front Wheel 12", pneumatic, Low Profile, 6" rim - Orion Pro 3 & 4-Wheel	--	1
1		1637432	Rear Wheel 12", pneumatic, Low Profile, 6" rim, RH - Orion Pro	--	1
2		1608812	Rear Wheel 12", pneumatic, Low Profile, 6" rim, LH - Orion Pro	--	1
3		1638977	Tire 11x3.50-6	--	1
3		1638974	Tire 12x4.00-6	--	1
4		1638978	Tube 11x3.50-6	--	1
4		1638975	Tube 12x4.00-6	--	1
5		1638835	6" Rim Kit, Front - Metro 3-Wheel	--	1
5		1638840	6" Rim Kit, Front - Metro 4-Wheel	--	1
6		1638836	6" Rim Kit, Rear - Metro LH/RH & Pro RH	--	1
7		1638861	6" Rim Kit, Rear - Pro LH	--	1


Anti-tipper Wheel



SCP_113A

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1522021	Wheel, Anti-Tipper	--	1

Propulsion Metro 12 km/h




SCP_114A

Propulsion Metro 12 km/h

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1638833	Differential 12 km/h	--	1
2		1638791	Motor 12 km/h - Metro	--	1
3		1638969	Clutch & Brake	--	1
4		1505672	Disclutching Lever 113.2 mm	--	1
5	A	1571584	Guide	--	1
6		1527845	Carbon Brush (4 pcs)	--	1
NOTE: A - For 12 km/h Metro					

Propulsion Pro 15 km/h




SCP_115A

Propulsion Pro 15 km/h

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1638834	Differential 15 km/h - Pro	--	1
2		1638831	Motor 15 km/h - Pro	--	1
3		1636891	Clutch & Brake 15 km/h (Pro)	--	1
4		1505672	Disclutching Lever 113.2 mm	--	1
5		1571584	Guide	--	1
6		1636874	Speed Sensor, cpl.	--	1
7		1636877	Carbon brush 15 km/h (Pro)	--	1

Front Brake 3-Wheel - Metro




SCP_116A

Front Brake 3-Wheel - Metro

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1638758	Brake assembly - Metro 3-Wheel	--	1
2		1638841	Cable assembly, Front brake	--	1
3		1518418	Brake Shoes, cpl. (95mm)	--	1

Front Brake 3-Wheel - Pro




SCP_117A

Front Brake 3-Wheel - Pro

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		SP1522025	Brake Assembly 3 Wheeler	--	1
2		1518419	Cable, Brake, Three Wheel	--	1
3		1518418	Brake Shoes, cpl. (95mm)	--	1

Front Brake 4-Wheel - Metro/ Pro




SCP_118A

Front Brake 4-Wheel - Metro/ Pro

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1521589	Brake Assembly, Four Wheel	--	1
2		1518422	Cable, Brake, Four Wheel	--	1
3		1518418	Brake Shoes, cpl. (95mm)	--	1


Steering Head



Steering Head

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1638667	Steering Bar Assembly, Pegasus Metro	--	1
2		1636816	Bearing Kit	--	1



Steering Complete



Steering Complete

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1608184	Steering	--	1
2		1636767	Grip, Handle Bar, cpl.	--	1
3		1608201	Cable, steering adj.	--	1
4		1636812	Joint, steering adjustment, cpl.	--	1
5		1608205	Front steering cover	--	1
6		1636804	Lever, Brake - RH	--	1
7		1636813	Lever, steering adjustment	--	1
8		1636814	Rear steering cover	--	1
9		1608202	Steering adjustment Assembly, cpl.	--	1

Front Fork 3-Wheel - Metro





SCP_121A

Front Fork 3-Wheel - Metro

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1638757	Front fork assembly, cpl. - Metro	--	1
2		1638779	Front Axle, cpl. - Metro 3-Wheel	--	1

Front Fork 3-Wheel - Pro




SCP_122A

Front Fork 3-Wheel - Pro

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1638781	Front fork assembly, cpl. - Pro 3-Wheel	--	1
2		1638759	Shock absorber, cpl. - Pro 3-Wheel	--	1
3		1638778	Front Axle, cpl. - Pro 3-Wheel	--	1


Front Suspension 4-Wheel - Metro



Front Suspension 4-Wheel - Metro

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1518078	Rod Assembly, Rear	--	1
2		1518237	Axle, Front - Left	--	1
3		1518238	Axle, Front - Right	--	1
4		1518077	Rod Assembly, Front	--	1
5		1505487	Kit, Bearing, Scooter Frame	--	1
6		1518196	Spring, Front Wheel	--	1
7		1518076	Swing holder, front	--	1
8		1563328	Suspension arm fixation	--	1


Front Suspension 4-Wheel - Pro



Front Suspension 4-Wheel - Pro

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1638787	Front Axle Kit LH - Pro 4-Wheel	--	1
2		1638782	Front Axle Kit RH - Pro 4-Wheel	--	1
3		1638788	Shock absorber, cpl. - Pro 4-Wheel	--	1
4		1636872	Rod Assembly, cpl.	--	1
5		1636859	Bearing Kit	--	1

Batteries and Cables




Batteries and Cables

Batteries and Cables


Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
2		1636819	Fuse battery cable (70A) PRO	--	1
3		1636823	Fuse 70A	--	1
4		1608209	Wire, Battery Output PRO	--	1
5		SP1521644	Off Board Charging Cable	--	1
6		1521758	Battery Strap cpl.	--	1
For Item #1 12V/50Ah or 12V/75Ah Battery see Scooter Price List. Scooters use (2) Batteries.					

Battery Charger



Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1123249	Charger, Battery, 8 Amp. 115vac (Off Board)	--	1
2		1121086	Single Mode Round Plug w/ AC Power Line Cord Cord, AC Power, 115VAC, 3 Prong On-Board (6' 6-1/2" L)	--	1

Electronic Components





SCP_127A

Electronic Components

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1506153	Buzzer	--	1
2		1523534	Cover, Charger Port	--	1
3		1523543	Charger Port	--	1
4		1523544	Bracket, Charger-Inlett	--	1
5		1505451	Ignition, Key Cylinder (Standard)	--	1
6		1638848	P.C. Board LED	--	1
7		1638849	Meter & Lamp, P.C. Board LED	--	1
8		1505462	Variable Resistor, Complete (30.0KO)	--	1
9		1638847	Display cover & rubber keys (LED Display)	--	1
10		1608207	Knob, Speed control	--	1
11		1608204	Temperature sensor	--	1
12		1608185	Steering Upper Assembly LCD	--	1
13		1608200	Wig Wag, cpl.	--	1
14		1640514	Wire Harness - Metro	--	1
14		1608212	Wire harness - Pro	--	1
15		1638963	Steering Upper Assembly LED	--	1
16		1637215	Display cover & rubber keys (LCD Display)	--	1
17		1608208	P.C. Board LCD	--	1

Powerbox and Limit Switch Trigger




SCP_128A

Powerbox and Limit Switch Trigger

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1556281	Speed Limit Switch Box Comet Pro 15 km/h	--	1
2		1556281	Speed Limit Switch Box Comet Pro 15 km/h	--	1
3		1608230	Trigger, limit switch	--	1
		1581766	Adjustment Gauge, Speed Reduction Switch	--	1
NOTE: A - Not Shown. When mounting the Speed Reduction Switch, the Adjustment Gauge part #SP1581766 is needed.					

Safety




SCP_129A



Safety

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		SP1518241	Rear-View Mirror Left/Right, Black	--	1
2		1637911	High-end rear mirror set	--	1
3		1518274	Lap Belt with Auto Lock Comet & Orion	--	1
4		SP1521592	Safety Vest	--	1

Seating



- | | |
|---|--------------------|
| ① | 508mm / max. 100kg |
| ② | 533mm / max. 100kg |
| ③ | 558mm / max. 100kg |
| ④ | 508mm / max. 160kg |
| ⑤ | 533mm / max. 160kg |
| ⑥ | 558mm / max. 160kg |





SCP_130A

Seating

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1640422	Seat Suspension 201.5mm 201lbs	--	1
2		1640423	Seat Suspension 201.5mm 285lbs	--	1
3		1640424	Seat Suspension 226.5mm 201lbs	--	1
4		1640425	Seat Suspension 226.5mm 285lbs	--	1
5		1640426	Seat Suspension 251mm 201lbs	--	1
6		1640427	Seat Suspension 251mm 285lbs	--	1
7		1521684	Seat Cover Rain & Dirt protection (Not HD Seat!)	--	1

Walking Support




SCP_131A

Walking Support

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1528274	Crutchholder, Complete	--	1
2		1505278	Walking Frame Holder, Complete	--	1








Rollator Bracket



Rollator Bracket

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1553851	Rollator Holder	--	1
2		1560465	Cane peel, cpl.	--	1
3		1553890	Rollator hooks	--	1
4		1553891	Belt, cpl.	--	1
5		1552594	Extension tubes (Comet only)	--	1
Item #6 is not available					

Storage






SCP_133A

Storage

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1505924	Bracket, Front Basket, cpl.	--	1
2		1637915	Bigger front basket	--	1
3		1637910	Front Box Assembly - black	--	1
4		1637913	Lockabel rear storage box	--	1
5		1505727	Rear Basket, cpl.	--	1
6		1505731	Bag, Seat Rear	--	1
7		1505729	Oxygen Bottle Holder (complete with 3 screws)	--	1

Bumpers









SCP_134A

Bumpers

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1528275	Front Bumper, 3-Wheeler	--	1
2		1638965	Front steel bumper - 4-Wheel	--	1
3		1638966	Rear Steel Bumper	--	1

Other Accessories



SCP_135A

Other Accessories

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1		1543160	Windscreen	--	1
2		1518240	Seat Cover, Comet - Black Only	--	1
3		1637914	Footrest, adjustable	--	1
4		1521594	Beverage Holder	--	1
5		1637917	USB charger	--	1
6		1637912	Mobile phone holder	--	1

Invacare Canada LP

www.invacare.ca



570 Matheson Blvd E;8
Mississauga, Ontario,
L4Z 4G4 Canada
800-668-5324

Technical Support
8006685324 x2655
Customer Service
800-668-5324

Literature
800-668-5324
or www.invacare.ca

© 2022 Invacare Corporation. All rights reserved. Republication, duplication or modification in whole or in part is prohibited without prior written permission from Invacare. Trademarks are identified by ™ and ®. All trademarks are owned by or licensed to Invacare Corporation or its subsidiaries unless otherwise noted. (Revised 01/25/22)

Form No. 18-075C