


# One Piece Adjustable Angle Footplate



## One Piece Adjustable Angle Footplate

Item Number	Note	Part Number	Description	Quantity Per	
				Package	Assembly
1	A	1046409-NLA	NLA - Tube, Support, Footplate, Black	1	2
1	B	1056296-NLA	NLA - Support Tube, Footplate	1	2
2		1046421-NLA	NLA - Screw, Socket Flat Head (1/4-28 x 1")	1	4
3		1013971	Washer, Nylon (1/4 x 7/16 x 5/64")	1	8
4		1025765	Locknut (1/4-28)	1	4
5	A	1046420-NLA	NLA - Hardware, Footplate Mounting	1	1
5	B	1058667-NLA	NLA - Hardware, Footplate Mounting	1	1
6		See Chart	#6 One Piece Adjustable Angle Footplate	1	1

NOTE: A - Tapered Front  
B - "V" Front

### #6 One Piece Adjustable Angle Footplate

Front Frame Type	14" Wide	15" Wide	16" Wide	17" Wide	18" Wide	19" Wide	20" Wide
Tapered	1046461-NLA	1046462	1046463-NLA	1046464-NLA	1046465-NLA	1046466	1046466
"V" Front	1046462	1046463-NLA	1046464-NLA	1046465-NLA	1046466	1046466	1046466